

MEDISANA®

Made for Life

PRODUCT CATALOGUE 2017



Mobile Health Health Control Sport Wellness Therapy
Healthy Home Personal Care Kitchen Scales

Live healthier, live better.

MEDISANA – Made for Life.

Mark, 28 & Tom, 18
Training together for the next triathlon.

Mia, 10
Has just overcome a cold.

Alina, 37
Enjoying her pregnancy.

Paula, 41
Feels perfectly fit.

Ralph, 72
Wants to live 100 years.

Health has been our main objective for more than 30 years. Healthcare is our core competence. In recent years, health awareness has evolved into a mega trend that is increasingly changing society. On the list of the most important values for people, health ranks in at first place. This new mind-set decisively influences people's behaviour and lifestyle. Confidence and self-responsibility are the driving forces for a new health culture.

MEDISANA regards itself as the leading specialist in the field of healthcare, and feels it has a special responsibility. We are constantly aiming to improve our products and also optimize the quality, and we are continuously developing products and services that address the newly emerged trend towards taking an active responsibility for the own body and health. Our products and services are Made for Life. Made for people in all life situations.



Trust is created through quality.

The healthcare market is a trust-based market. For us this means, MEDISANA products must meet a broad range of requirements. Therefore, we place a key emphasis on product quality, ease of use and safety. We meet the criteria of the Medical Devices Directive, MDD and the IVD Directive (in vitro diagnostics), but also the criteria of various consumer organisations. This is emphasized by our many certificates and awards. The most important criteria, however, is the trust our customers have in our products. That is our benchmark.

Contents

Mobile Health	8
Medisana App-World	10
VitaDock+ App 2.0	12
Health Control	14
Sleep Monitor with Bluetooth 4.0 Sleepace	16
Upper Arm Blood Pressure Monitor BU 530 connect	18
Upper Arm Blood Pressure MTP	20
Bluetooth Thermometer TM 735	21
Body Analysis Scale TargetScale 3	22
Personal Scale PS 430	24
Product overview	26
Sport	42
Foam Roller PowerRoll Soft	45
Foam Roller PowerRoll Pro	46
Foam Roller VarioBall	47
Activity Tracker ViFit touch	48
Product overview	50
Wellness	54
Shiatsu Massage Seat Cover MCN Pro	56
Neck Massage Vibration NM 865	58
Mini-Hand Massager HM 858	59
Heat Pad HP 625 comfort	60
Product overview	62
Therapy	76
Daylight Lamp LT 460	78
Daylight Lamp LT 470	79
Daylight Lamp LT 480	80
Ultrasonic Inhaler IN 525	81
Wake-Up Alarm Clock WL 444	82
Wake-Up Alarm Clock WL 460	83
Product overview	84
Healthy Home	90
Air Humidifier AH 660	92
Air Humidifier AH 662	93
Aroma Diffuser AD 620	94
Aroma Diffuser AD 630	95
Aroma Diffuser AD 640	96
Room Air Purifier Medisana Air	97
Product overview	98
Personal Care	102
Manicure and Pedicure Device MP 840	104
Cellulite Massager AC 855	105
Facial Cleansing Brush FB 885	106
Cosmetic Mirror CM 835	108
Product overview	110
Kitchen Scales	116
Digital Kitchen Scale with App KS 250	118
Product overview	120
Accessories	124

Marc, 32
Always has his family doctor at his side.

Patrick, 31
Monitors his vitals himself.

Frida, 27
Uses mobile diagnostics.

More quality for life.

Mobile health management – Made for Life.

For nearly 90% of those over the age of 14 in Germany, health is the core quality of life aspect*. More and more people taking their health into their own hands. Checks, prevention, fitness, wellness – modern man actively structures his own health. The key factor is the adequate information about the own vital signs. MEDISANA supports the self-conducted health management with integrated solutions – mobile and online.

* Source: Philips Health Study 2015. How health is becoming the driving force of a new health culture.

Getting to know your own body.

The MEDISANA App-World enables qualified health monitoring and promotes personal health. All apps have the goal of promoting self-responsibility for the own health.

Knowledge is helpful. Sufficient knowledge of the vital signs, health values, sleep activity, the nutritional values of the food with which you want to cook, or the intensity of massage programmes. The MEDISANA App World gives you more freedom to consciously adapt your ambitions to your own body.




VitaDock+ App 2.0

The free VitaDock+ App for all vital signs, which can then be clearly displayed in the form of daily, weekly, monthly and yearly views.




Kitchen scales App

A balanced diet is our common goal. The kitchen scales App provides reliable and comprehensive information on the key nutritional values.



Bodytoner App

What type of massage is good? And with what intensity can you attain the best result? With the help of the Body Toner app, the form and intensity of the massage can be individually controlled.



Sleepspace App

How well do you sleep? Sleep depth, duration, heart and respiratory rates and movements during sleep – all conveniently recorded.



VitaTherm App

The continuous monitoring of body temperature offers special comfort in the case of illness. An alarm will be triggered immediately if the temperature increases above a set value.



Thermo App

The body temperature in a continuous real-time monitoring. All values at a glance. All risks are minimised.



Exciting. All vital signs at a glance.

Activity, blood pressure, weight, blood sugar and blood oxygen. With using the new the VitaDock+ App it is possible to control all vital signs even better online and mobile.

The VitaDock+ App not only records all vital data. The user can add comments to any value, and have measurements and readings visualized. This makes health management

even more personalized. Another innovation: The simple online registration via Facebook.



Download the new VitaDock+ App 2.0:

Android

iTunes

NEW VitaDock+ App 2.0

- Newly developed user interface
- List of measurement results
- Comment functions
- Registration via Facebook
- Simplified connection of Bluetooth devices



Sandra, 27
Always up-to-date.

Max, 38
Always on top of his vitals.



Trendy. For the lifestyle and health.

Innovative health management – Made for Life.

Peter, 44
Has his diabetes under control.

Simon, 39
Regularly checks his blood pressure.

Isabelle, 37
Always has her weight in check.

Barbara, 40
Knows her body well.

My lifestyle, my everyday life, my health. People will in the future be making their lifestyles an indicator in determining health. My measurement values and needs determine how I deal with my body and my health. The willingness to take responsibility for oneself is growing as is the obvious need to keep an eye on one's own health. MEDISANA offers innovative products to support this willingness to check and evaluate vital signs. For easy and simple handling of your own health.



Daniel, 35
Starts fit into the day.

Mia, 32
Soundly asleep every night.

Sleep better.
Live healthier.


We sleep about 3000 hours per year. During this time, the cells regenerate, the metabolism is regulated, and the immune system strengthened. Good sleep keeps us fit. Why do we pay so little attention to it?

The Sleep Monitor Sleepace can monitor our sleep on our behalf. It analyses your sleeping habits and monitors a broad range of values. In the connected Sleepace App, measurement results are stored and evaluated.



Automatic data transfer
via Bluetooth Smart to Sleepace App

Download
the Sleepace App:



Sleep better with Sleepace
sleep monitoring.



NEW Sleep Monitor
Bluetooth 4.0 Sleepace

- Accurate and comprehensive analysis of the sleeping habits
- Monitors depth and duration of sleep, heart rate, respiration rate and movements when asleep



Annette, 42
Everything is totally fine.

Upper Arm Blood Pressure Monitor BU 530 connect

- Accurate blood pressure measurement
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Display items: systolic, diastolic, pulse, date, time
- Classification of the measurement values according to WHO*
- Certified medical device**



* World Health Organization ** Certified according to MDD (Medical Device Directive)

A tedious chore will become a pleasant task

Blood pressure is vital to the human circulatory system. Too low or too high – both can be hazardous to health. Therefore the monitoring of blood pressure is highly recommended, not only for diagnosed patients.

What could be a tedious control measure, will become a pleasant task with the upper arm blood pressure monitor BU 530. Precise measurement, intuitive operation and thanks to large digits easy to read.

The measurement data can be stored and transferred to the VitaDock+ App via Bluetooth. Another highlight: the automatic calculation of the average value from the last three measurements.

Excellent readability through extra large digits.



Dr. med. Johannes Wimmer helps users to understand the measurements in VitaDock+ App better and thus live healthier.



Android



iTunes



Download
the new
VitaDock+
App 2.0:



Precision provides safety.

Measuring blood pressure at the doctor’s office is a routine examination. But not at home. With the upper arm blood pressure monitor MTP, the measurement will become a precise, safe and comfortable examination.

NEW, now with WHO* and arrhythmia detection



TOP-SELLER
1.000.000
units sold

According to the World Health Organization (WHO), one in three adults in Germany over 35 suffers from hypertension – often without knowing it, because high blood pressure often doesn’t cause any noticeable health problems over a long period of time. The MTP provides more than a million people with accurate information about their systolic and diastolic values. Simple, convenient and secure.

The MTP also gives risk patients a reliable safety.



Easy to use one-button operation

NEW Upper Arm Blood Pressure Monitor MTP

- NEW: Measurements are classified according to the WHO*
- NEW: Arrhythmia detection
- Accurate blood pressure measurement
- Extra large, easy-to-read display
- 99 memory spaces for each of 2 users
- Display items: Systole, diastole, pulse, date, time
- Certified medical device**

ETM TESTMAGAZIN - URTEIL	
MEDISANA Blutdruckmessgerät MTP	
GUT	89,3 %
Im Test: 9 Oberarm-Blutdruckmessgeräte	
Testurteile: 4x sehr gut, 5x gut	
Heft 02/2015	



Benjamin, 6 weeks
Well protected.

NEW Bluetooth Thermometer TM 735 connect

- Continuous real-time monitoring
- Alarm is triggered when exceeding the set temperature zone
- Clear measurement display
- Certified medical device**

Knowledge gives peace of mind.

Continuous real-time monitoring of the body temperature.

It is very important to monitor the body temperature in real-time during a fever – but it has been a lot more cumbersome until now. The Bluetooth thermometer TM 735 makes the monitoring very efficient and comfortable. Fixed directly onto the skin below the armpit using dermatologically compatible pads, the TM 735 will send all measurement values directly to the smartphone.

Another highlight: the memory function for the medication intake.

Always maintain an overview - and know that you can respond immediately.

Download the Thermo App:

Android

iTunes

* World Health Organization ** Certified according to MDD (Medical Device Directive)

The goal – my desired weight.

This will ease the way to reach your desired weight. The body analysis scale TargetScale 3 records all body values and visualizes the difference to the target weight via light rings.

TargetScale 3 supports the attainment of your desired weight with extensive body value measurements: Weight, Body Mass Index, kcal per day, percentage of body fat and body water, as well as muscle and bone mass. Up to 240 measurements per user can be stored. By connecting to the free app VitaDock+ and the online platform, the values can be analysed and, for example, shared with physicians.

Josephine, 26
Wants precise values.

Advanced functions via the Bluetooth connection to the free VitaDock+ App for iOS and Android



Beautifully designed, innovative functions.



Android



iTunes



Download the new VitaDock+ App 2.0:



One scale for everyone. Automatic recognition of up to eight users.



NEW Body Analysis Scale TargetScale 3

- NEW: Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app
- NEW: Storage of 240 measurements
- NEW: Storage and graphic display of measurement results over time
- Precise determination of weight, body mass index, kcal per day, Amount of body fat and body water, as well as muscle and bone mass
- Lighted rings depict the difference between current and your target weight



Maria, 35
Loves her slender legs.

Frank, 37
Loves Mary's entire shape.

NEW Personal Scale PS 430

- Anti-slip surface for safe foothold
- Energy saving auto-off in six seconds
- 4 high precision strain gauge sensors for accurate results



ANTI-SLIP
SURFACE

Stand securely. Safely weigh yourself.

Quickly weigh yourself on the scale after the morning shower. For many a ritual – that should be safe. The personal scale PS 430 provides a secure footing with the anti-slip surface.

Four high-precision strain gauge sensors cannot be fooled. The exact measurement results on the easily readable LCD display provide precise information on

the actual weight. With its 6 mm safety glass which has an anti-slip surface, this scale makes the use in the bathroom safe and convenient.



Beautifully designed, elegant and safe. With its ultra-thin design, the PS 430 will integrate perfectly into any bathroom.



Upper Arm Blood Pressure Monitor
BU 530 CONNECT



- Accurate blood pressure measurement on the upper arm
 - Arrhythmia indicator
 - 120 memory spaces for each of 2 users
 - Very large digits for easy reading
 - Display items: systolic, diastolic, pulse, date, time
 - Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*
-
- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
 - Calculates average value of data of the last 3 results
 - Measurement for guests without storage of reading possible
 - Battery operation (power adapter optional available)
 - Very large cuff (arm circumferences 22-36 cm)
 - Incl. storage bag
 - Incl. batteries (4 x 1.5V AAA)
 - Certified medical device***

Size:	approx. 15 x 10 x 6 cm
Weight:	approx. 0.4 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 0-299 mmHg Pulse: 40-199 beats/min.
Batteries:	4 x 1.5V (AAA)
Article No.:	51174
EAN:	40 15588 51174 5
Warranty:	3 Years



Upper Arm Blood Pressure Monitor
BU 550 CONNECT



- Accurate blood pressure measurement on the upper arm
 - Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
 - 500 memory spaces for each of 2 users
 - Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*
-
- Arrhythmia indicator
 - Very large digits for easy reading
 - Display items: systolic, diastolic, pulse, date, time
 - Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
 - Calculates average value of data of the last 3 results
 - Measurement for guests without storage of reading possible
 - Battery operation (power adapter optional available)
 - Ultra slim design
 - Very large cuff (arm circumferences 22-42 cm)
 - Incl. storage bag
 - Incl. batteries (4 x 1.5V AAA)
 - Certified medical device***

Size:	approx. 9 x 14 x 3.3 cm
Weight:	approx. 0.35 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 40-230 mmHg Pulse: 40-199 beats/min.
Batteries:	4 x 1.5V (AAA)
Article No.:	51290
EAN:	40 15588 51290 2
Warranty:	3 Years

2in1
Upper Arm Blood Pressure Monitor
BU 575 CONNECT



- Accurate blood pressure measurement on the upper arm
 - Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
 - 180 memory spaces for each of 2 users
 - Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*
-
- Test reminder thanks to integrated alarm clock
 - Rechargeable lithium-ion battery operation (power adapter optional available)
 - Incl. USB charging cable
 - Further technical data like BU 550 connect
- Alarm clock mode**
- Display of date, time, alarm
 - 2 alarm settings adjustable
 - “Snooze” function
 - Alarm with increasing volume

Size:	approx. 6 x 11 x 3 cm
Weight:	approx. 0.35 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 40-230 mmHg Pulse: 40-199 beats/min.
Batteries:	Lithium ion accumulator
Article No.:	51296
EAN:	40 15588 51296 4
Warranty:	3 Years

Upper Arm Blood Pressure Monitor
MTX CONNECT



- Higher accuracy of measurements for severely fluctuating blood pressure thanks to patented 3MAM technology (weighted average of 3 consecutive measurements)
 - 99 memory spaces for each of 2 users
 - Very large digits for easy reading
 - Bluetooth® data transfer to to the VitaDock App for iOS and Android and to VitaDock Online
-
- Arrhythmia indicator
 - Display items: systolic, diastolic, pulse, date, time
 - Easy-to-use „EasyForm“ cuff (arm circumferences 22-32 cm)
 - Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
 - Calculates average value of all stored data
 - Automatic time and date setting via VitaDock app
 - Battery operation (power adapter optional available)
 - Incl. storage bag
 - Certified medical device***

Size:	approx. 12,5 x 16,3 x 6,5 cm
Weight:	approx. 0.45 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 30-280 mmHg Pulse: 40-200 beats/min.
Batteries:	4 x 1.5 V (AA)
Article No.:	51085
EAN:	40 15588 51085 4
Warranty:	3 Years

Upper Arm Blood Pressure Monitor
MTX



- Accurate blood pressure measurement on the upper arm
 - 99 memory spaces for each of 2 users
 - Very large digits for easy reading
 - USB port for data transfer
-
- Arrhythmia indicator
 - Display items: systolic, diastolic, pulse, date, time
 - Easy-to-use „EasyForm“ cuff (arm circumferences 22-32 cm)
 - Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
 - Calculates average value of all stored data
 - Automatic time and date setting via integrated radio clock
 - Battery operation (power adapter optional available)
 - Incl. USB cable
 - Incl. storage bag
 - Certified medical device***

Size:	approx. 12,5 x 16,3 x 6,5 cm
Weight:	approx. 0.45 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 30-280 mmHg Pulse: 40-200 beats/min.
Batteries:	4 x 1.5 V (AA)
Article No.:	51083
EAN:	40 15588 51083 0
Warranty:	3 Years

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3 **World Health Organisation *** Certification according to MDD (Medical Device Directive)

Upper Arm Blood Pressure Monitor

MTP




- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 99 memory spaces for each of 2 users
- Extra large, easy-to-read display
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of all stored data
- 1-button operation
- Battery operation (power adapter optional available)
- Large cuff (arm circumferences 22-32 cm)
- Incl. storage bag
- Incl. blood pressure pass
- Incl. batteries (4 x 1.5V AA)
- Certified medical device***

Size:	approx. 12 x 16,3 x 6,6 cm
Weight:	approx. 0.5 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 30-280 mmHg Pulse: 40-200 beats/min.
Batteries:	4 x 1.5 V (AA)
Article No.:	51090
EAN:	40 15588 51090 8
Warranty:	3 Years

Upper Arm Blood Pressure Monitor

BU 505



- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 120 memory spaces for each of 4 users
- Very large digits for easy reading (up to 31,3 mm height)
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of all stored data
- Battery level indicator
- Very large cuff (arm circumferences 22-36 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AA)
- Certified medical device***

Size:	approx. 9,8 x 13,7 x 4,5 cm
Weight:	approx. 0.35 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 30-250 mmHg Pulse: 40-180 beats/min.
Batteries:	4 x 1.5 V (AAA)
Article No.:	51185
EAN:	40 15588 51185 1
Warranty:	3 Years

Upper Arm Blood Pressure Monitor

BU 510




- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 90 memory spaces for each of 2 users
- Very large digits for easy reading (up to 22 mm height)
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of all stored data
- Battery operation (power adapter optional available)
- Very large cuff (arm circumferences 22-36 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AA)
- Certified medical device***

Size:	approx. 13 x 11 x 6 cm
Weight:	approx. 0.3 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 40-260 mmHg Pulse: 40-180 beats/min.
Batteries:	4 x 1.5 V (AA)
Article No.:	51160
EAN:	40 15588 51160 8
Warranty:	3 Years

Upper Arm Blood Pressure Monitor

BU 535




- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Very large cuff (arm circumferences 22-36 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AAA)
- Certified medical device***

Size:	approx. 15 x 9,5 x 4,1 cm
Weight:	approx. 0.23 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 40-260 mmHg pulse: 40-180 beats/min.
Batteries:	4 x 1.5 V (AAA)
Article No.:	51176
EAN:	40 15588 51176 9
Warranty:	3 Years

Upper Arm Blood Pressure Monitor

BU 535 VOICE





- Accurate blood pressure measurement on the upper arm
- Voice output in selectable languages: DE, GB, NL, ES, FR, IT
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Very large cuff (arm circumferences 22-36 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AAA)
- Certified medical device***

Size:	approx. 15 x 9,5 x 4,1 cm
Weight:	approx. 0.23 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 40-260 mmHg pulse: 40-180 beats/min.
Batteries:	4 x 1.5 V (AAA)
Article No.:	51178
EAN:	40 15588 51178 3
Warranty:	3 Years

Upper Arm Blood Pressure Monitor

MTS




- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 60 memory spaces for each of 2 users
- Very large digits for easy reading
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Automatic switch-off after 1 min.
- Battery operation (power adapter optional available)
- Large cuff (arm circumferences 22-30 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AA)
- Certified medical device***

Size:	approx. 19,5 x 22,5 x 26,5 cm
Weight:	approx. 0.3 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 0-300 mmHg pulse: 30-180 beats/min.
Batteries:	4 x 1.5 V (AA)
Article No.:	51152
EAN:	40 15588 51152 3
Warranty:	3 Years

Wirst Blood Pressure Monitor
BW 300 CONNECT




- Accurate blood pressure measurement on the wrist
- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
- 180 memory spaces for each of 2 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- Arrhythmia indicator
- Very large digits for easy reading
- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Measurement for guests without storage of reading possible
- Ultra slim design
- Cuff for wrist circumferences 13.5-21.5 cm
- With rechargeable lithium-ion battery and USB charging cable
- Certified medical device***

Size:	approx. 6.5 x 8 x 1.3 cm
Weight:	approx. 0.1 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 40-230 mmHg pulse: 40-199 beats/min. Lithium ion accumulator
Batteries:	
Article No.:	51294
EAN:	40 15588 51294 0
Warranty:	3 Years

Wirst Blood Pressure Monitor
BW 310




- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 60 memory spaces for each of 2 users
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Incl. storage box
- Incl. battery (2 x 1.5V AAA)
- Certified medical device***

Size:	approx. 7.2 x 7.5 x 3 cm
Weight:	approx. 0.15 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 30-250 mmHg pulse: 40-180 beats/min. 2 x 1.5 V (AAA)
Batteries:	
Article No.:	51070
EAN:	40 15588 51070 0
Warranty:	3 Years

Wirst Blood Pressure Monitor
BW 333




- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Extra large, easy-to-read display
- Display items: systolic, diastolic, pulse, date, time
- Ultra slim design – only 1,5 cm

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Automatic switch-off
- Battery level indicator
- Cuff for wrist circumferences 13.5-23 cm
- Incl. battery (2 x 1.5V AAA)
- Certified medical device***

Size:	approx. 8.4 x 7.7 x 1.5 cm
Weight:	approx. 0.12 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 30-250 mmHg pulse: 40-180 beats/min. 2 x 1.5 V (AAA)
Batteries:	
Article No.:	51075
EAN:	40 15588 51075 5
Warranty:	3 Years


Wirst Blood Pressure Monitor
BW 335



- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading
- Display items: systolic, diastolic, pulse, date, time
- Ultra slim design

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Cuff for wrist circumferences 14-19.5 cm
- Incl. storage box
- Incl. batteries (2 x 1.5V AAA)
- Certified medical device***

Size:	approx. 8.9 x 6 x 3.1 cm
Weight:	approx. 0.07 kg
Measuring method:	Oscillometric
Measuring range:	Blood pressure: 0-299 mmHg pulse: 40-180 beats/min. 2 x 1.5 V (AAA)
Batteries:	
Article No.:	51077
EAN:	40 15588 51077 9
Warranty:	3 Years



Sleep Monitor
SLEEPSPACE





- Monitors depth and duration of sleep, heart rate, respiration rate and movements when asleep
- Accurate and comprehensive analysis of sleeping habits







- Ultra-thin piezoelectric band (1.5 mm) for a comfortable fit between the mattress and the sheets
- Patented activation system
- Rechargeable lithium-ion battery
- Incl. USB charging cable
- Sleepspace App functions**
- Automatic data transfer via Bluetooth Smart to the Sleepspace App
- Saves and displays all measurement results clearly over time: day, week or month views

Size sensor:	approx. 8,8 x 6,5 x 1,5 cm
Size band:	approx. 79,5 x 6,5 x 0,15 cm
Weight:	approx. 160 g
Battery:	Rechargeable 3,7V lithium-ion battery
Article No.:	79490
EAN:	40 15588 79490 2
Warranty:	3 Years

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3 **World Health Organisation *** Certification according to MDD (Medical Device Directive)

Body Analysis Scale
BS 444 CONNECT





- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- High-quality stainless steel electrodes
- Convenient "step-on" switch
- Large easy-to-read LCD display
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)
- Maximal load capacity: 180 kg
- Incl. batteries (3 x 1.5 V AAA)

Size: approx. 30 x 30 x 2.3 cm

Weight: approx. 1.6 kg

Display size: 7.3 x 4.8 cm

Weight display: 5 bis 180 kg

Graduation: 100 g, 0.2 lb, 1/4 lb

Batteries: 3 x 1.5 V (AAA)

Article No.: 40444

EAN: 40 15588 40444 3

Warranty: 3 Years



INFO Body Analysis Scale
BS 444 connect







It is so easy to comprehensively measure and analyse your own vital signs.

Is the BMI value correct? What is the situation with your body fat? How many kilocalories does the body need? The Body analysis scales BS 444 connect give reliable and precise information. The scales measure five different vital signs, determine the BMI value and allow convenient evaluation and analysis of the measurement values via Bluetooth transmission to the new VitaDock+ App.



Body Analysis Scale
BS 430 CONNECT





- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- High-quality stainless steel electrodes
- Convenient "step-on" switch
- Large easy-to-read LCD display with backlight
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)
- Maximal load capacity: 180 kg
- Incl. batteries (4 x 1.5 V AAA)

Size: approx. 32.5 x 32.5 x 2.5 cm

Weight: approx. 2.15 kg

Display size: 7 x 4.9 cm

Weight display: 5 bis 180 kg

Graduation: 100 g, 0.2 lb, 1/4 lb



Batteries: 4 x 1.5 V (AAA)








Article No.: 40422

EAN: 40 15588 40422 1

Warranty: 3 Years

Body Analysis Scale
BS 440 CONNECT





- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- High-quality ITO electrodes
- Convenient "step-on" switch
- Large easy-to-read LCD display with backlight
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)
- Maximal load capacity: 180 kg
- Incl. batteries (4 x 1.5 V AAA)

Size: approx. 32 x 32 x 2.7 cm

Weight: approx. 2 kg

Display size: 7 x 4.9 cm

Weight display: 5 bis 180 kg

Graduation: 100 g, 0.2 lb, 1/4 lb



Batteries: 4 x 1.5 V (AAA)








Article No.: 40423

EAN: 40 15588 40423 8

Warranty: 3 Years

Body Analysis Scale
TARGETSCALE 3





- Precise determination of weight (up to 180 kg), Body Mass Index (BMI), kcal per day, calculates percentage of body fat and water as well as muscle and bone mass
- Multiuser: Automatic detection of up to 8 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- Lighted rings depict the difference between current and your target weight
- Large easy-to-read LCD display with white luminous digits
- Athlete mode for sportspersons
- Stores up to 240 measurements in the scale memory
- Maximal load capacity: 180 kg
- Incl. batteries (4 x 1.5 V AA)

Size: approx. ø 36 x 4.75 cm

Weight: approx. 2 kg

Display size: 7 x 5 cm

Weight display: 5 bis 180 kg

Graduation: 100 g, 0.2 lb, 1/4 lb






Batteries: 4 x 1.5 V (AA)

Article No.: 40413








EAN: 40 15588 40413 9

Warranty: 3 Years

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3

Personal Scale	Body Analysis Scale
PS 430	ISA
	
	
	
<ul style="list-style-type: none">• Anti-slip surface for safe foothold• Energy saving auto-off in six seconds• Easily readable LCD display	<ul style="list-style-type: none">• Measurement of weight, body fat, body water and muscle mass• Determination of the BMI value• Large easy-to-read LCD display with backlight
<ul style="list-style-type: none">• Convenient "step-on" switch• Treads of the standing area made of 6mm safety glass with rounded corners• Extra slim design• Toggles between kg, lb and st• Overload indication• Low battery display• Four high-precision strain gauge sensors for accurate results• Maximal load capacity: 180 kg• Incl. battery (1 x CR2032)	<ul style="list-style-type: none">• 10 memory slots• Automatic switch-off after 8 seconds• Convenient "step-on" switch• Tread area of 8 mm safety glass• Modern blue LED edge lighting• High-quality stainless steel electrodes• Non-slip tread area• Toggles between kg, lb and st• Four high-precision strain gauge sensors for accurate readings• Maximal load capacity: 180 kg• Incl. batteries (3 x 1.5 V AA)
<div><div>Size:</div><div>Weight:</div><div>Display size:</div><div>Weight display:</div><div>Graduation:</div><div>Battery:</div></div> <div><div>approx. 30 x 30 x 2.1 cm</div><div>approx. 1.5 kg</div><div>3.0 x 7.4 cm</div><div>bis 180 kg, 396 lb, 28 st 4 lb</div><div>100 g, 0.2 lb, 1/4 lb</div><div>1 x CR2032</div></div> <div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>40458</div><div>40 15588 40458 0</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Display size:</div><div>Weight display:</div><div>Graduation:</div><div>Batteries:</div></div> <div><div>approx. 31.9 x 31.9 x 2.15 cm</div><div>approx. 2.6 kg</div><div>10.5 x 4.4 cm</div><div>bis 180 kg, 396 lb, 28 st 4 lb</div><div>100 g, 0.2 lb, 1/4 lb</div><div>3 x 1.5 V (AA)</div></div> <div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>40480</div><div>40 15588 40480 1</div><div>3 Years</div></div>



XL Personal Scale	Body Analysis Scale	Personal Scale
PS 460	PSM	PSD
		
		
		
<ul style="list-style-type: none">• Particularly large tread area• Elegant design of high quality stainless steel• Large easy-to-read LCD display with red luminous digits• Maximal load capacity: 200 kg	<ul style="list-style-type: none">• Measurement of weight, body fat, body water and muscle mass• Integrated calorie analysis• Large easy-to-read LCD display with white luminous digits	<ul style="list-style-type: none">• Trendy retro design• Battery-free, mechanical operation• Large full-view analogue dial
<ul style="list-style-type: none">• Automatic switch-off after 8 seconds• Convenient "step-on" switch• Toggles between kg, lb and st• Four high-precision strain gauge sensors for accurate readings• Incl. batteries (6 x 1.5 V AAA)	<ul style="list-style-type: none">• 10 memory slots• Automatic switch-off after 8 seconds• Convenient "step-on" switch• Tread area of 8 mm safety glass• High-quality stainless steel electrodes• Toggles between kg, lb and st• Four high-precision strain gauge sensors for accurate readings• Maximal load capacity: 180 kg• Incl. batteries (3 x 1.5 V AAA)	<ul style="list-style-type: none">• Non-slip tread area• Solid metal casing
<div><div>Size:</div><div>Weight:</div><div>Display size:</div><div>Weight display:</div><div>Graduation:</div><div>Batteries:</div></div> <div><div>approx. 35 x 34 x 2 cm</div><div>approx. 2.6 kg</div><div>8.8 x 4.3 cm</div><div>bis 200 kg, 440 lb, 31 st 6.9 lb</div><div>100 g, 0.2 lb, 1/4 lb</div><div>6 x 1.5 V (AAA)</div></div> <div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>40433</div><div>40 15588 40433 7</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Display size:</div><div>Weight display:</div><div>Graduation:</div><div>Batteries:</div></div> <div><div>approx. 32 x 32 x 2.6 cm</div><div>approx. 2.5 kg</div><div>10.5 x 4.4 cm</div><div>bis 180 kg, 396 lb, 28 st 4 lb</div><div>100 g, 0.2 lb, 1/4 lb</div><div>3 x 1.5 V (AAA)</div></div> <div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>40446</div><div>40 15588 40446 7</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Weight display:</div><div>Graduation:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 28.5 x 42.5 x 8.5 cm</div><div>approx. 3.6 kg</div><div>bis 150 kg</div><div>500 g</div><div>40461</div><div>40 15588 40461 0</div><div>3 Years</div></div>

Bluetooth Thermometer
TM 720 CONNECT



Q1/2017

- Continuous real-time monitoring of body temperature
- The display shows the current temperature (Function also without app)
- Visual fever alarm by a red LED
- Universal user group: for small children to adults
- Bluetooth® Smart data transfer to the VitaTherm app for iOS and Android

- Pleasant wearing comfort
 - Armband with Velcro closer 20 cm long (Including extension strip 15 cm)
 - Application: On the arm, in the vicinity of the armpit
 - Including extension band
 - Incl. battery (1 x CR2032)
 - Certified medical device*
- Functions of the app**
- Neatly arranged display of measurement values in chronological order
 - Alarm is triggered when exceeding the set temperature zone
 - Reminding function for intake of medication
 - Alarm when Bluetooth® connection is disconnected (radius max. 10 m)
 - Unit of measurement can be selected (°C or °F)

Size:	approx. 3.6 x 3 x 1 cm Wrist band: 20 x 2.4 cm
Weight:	approx. 25 g
Display:	approx. 1.8 cm
Measuring range:	32-43°C
Accuracy:	± 0.1°C (35-42°C)
Battery:	1 x CR2032
Article No.:	76160
EAN:	40 15588 76160 7
Warranty:	3 Years

Bluetooth Thermometer
TM 735 CONNECT



Q1/2017

- Continuous real-time monitoring of body temperature
- Universal user group: for small children to adults
- Bluetooth® Smart data transfer to the ThermoApp for iOS and Android

- Where to use: on the skin, in the armpit
 - If using on children below 2 years of age, consult a doctor
 - Incl. battery (1x CR2025)
 - Certified medical device*
- Functions of the app**
- Neatly arranged display of measurement values in chronological order
 - Alarm is triggered when exceeding the set temperature zone
 - Reminding function for intake of medication
 - Alarm when Bluetooth® connection is disconnected (radius max. 10 m)
 - Unit of measurement can be selected (°C or °F)

Size:	approx. 6 x 3 x 0.6 cm
Weight:	approx. 16 g
Measuring range:	25-45°C
Accuracy:	± 0.05°C (35-38.5°C) und ± 0.1°C (< 35°C und > 38.5°C)
Battery:	1 x CR2025
Article No.:	76165
EAN:	40 15588 76165 2
Warranty:	3 Years



Incl. 10 x 3M pads

Bluetooth Thermometer
TM 740 CONNECT



Q1/2017

- Continuous real-time monitoring of body temperature
- Coloured LEDs on the thermometer indicate the temperature range – green: normal temperature, Yellow: increased temperature, Red: Fever (from 38 °C), Blue: Bluetooth connection to the smart-phone has been set up
- Bluetooth® Smart data transfer to the VitaTherm app for iOS and Android

- Very flat thermometer for pleasant wearing comfort
 - Where to use: on the skin, in the vicinity of the armpit
 - Universal user group: for small children to adults
 - Incl. battery (1 x CR2032)
 - Certified medical device*
- Functions of the app**
- Neatly arranged display of measurement values in chronological order
 - Alarm is triggered when exceeding the set temperature zone
 - Reminding function for intake of medication
 - Alarm when Bluetooth® connection is disconnected (radius max. 10 m)
 - Unit of measurement can be selected (°C or °F)

Size:	approx. 4.2 x 4.2 x 0.9 cm
Weight:	approx. 16 g
Measuring range:	32-42°C
Accuracy:	± 0.1°C (35-42°C)
Battery:	1 x CR2032
Article No.:	76163
EAN:	40 15588 76163 8
Warranty:	3 Years



Incl. 10 x 3M pads



INFO Bluetooth Thermometer
TM 740 connect

It is reassuring. Measuring and continuously monitoring the body temperature in real time.

The multifunction thermometer TM 740 makes it easy for you. It is tolerated well dermatologically when applied to the skin. The temperature is displayed in three coloured areas. The measurement values are continuously sent to your Smartphone via Bluetooth. Monitoring with a certified medical product is especially reassuring with small children.

Multifunctional Thermometer
TM 770










Q1/2017

- Precise measurement of the body temperature on the forehead in seconds
- Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100° C)
- Visual fever alarm with colour change – green: temperature normal, red: raised temperature

- Acoustic signal at the beginning and end of temperature reading
- Automatic storage of the last 32 measurement results
- Automatic power off after 30 seconds
- Unit selectable: °C or °F
- Incl. battery (2 x 1,5 AA)
- Certified medical device**

Size:	approx. 9.3 x 15.3 x 4.1 cm
Weight:	approx. 160 g
Measuring range:	Körper: 32-42.9°C Objekt: 0-100°C
Accuracy:	Stirn: ± 0.3°C (32-42.9°C) Objekt: ± 0.1°C (0-100°C)
Battery:	2 x 1.5V (AA)
Article No.:	76150
EAN:	40 15588 76150 8
Warranty:	3 Years

* Certification according to MDD (Medical Device Directive)


38	<div>Infrared Multi-Functional Thermometer</div> <div>TM 750 CONNECT</div> <div></div> <div></div> <div>Q1/2017</div> <div><ul style="list-style-type: none">• Precise measurement of the body temperature in the ear and on the forehead in seconds• Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100° C)• Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*</div> <div><ul style="list-style-type: none">• Automatic storage of the last 30 measurement results• Visual fever alarm with colour change – green: temperature normal, red: raised temperature• Acoustic signal at the beginning and end of temperature reading• Display of date and time• Incl. storage pouch• Incl. batteries (2 x 1.5V AAA)• Certified medical device**</div> <div><div>Size: approx. 13,4 x 4,2 x 4,7 cm</div><div>Weight: approx. 65 g</div><div>Measuring range: Körper: 34-43°C Objekt: 0-100°C</div><div>Accuracy: Stirn: ± 0,3°C (34-43°C) Ohr: ± 0,3°C (35,5-42°C) Objekt: ≤30 ±1.5°C / >30 ±5%</div><div>Batteries: 2 x 1.5V (AAA)</div><div>Article No.: 76145</div><div>EAN: 40 15588 76145 4</div><div>Warranty: 3 Years</div></div>	<div>Infrared Multi-Functional Thermometer</div> <div>TM 750</div> <div></div> <div></div> <div><ul style="list-style-type: none">• Precise measurement of the body temperature in the ear and on the forehead in seconds• Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100° C)• Automatic storage of the last 30 measurement results</div> <div><ul style="list-style-type: none">• Visual fever alarm with colour change – green: temperature normal, red: raised temperature• Acoustic signal at the beginning and end of temperature reading• Display of date and time• Incl. storage pouch• Incl. batteries (2 x 1.5V (AAA)• Certified medical device**</div> <div><div>Size: approx. 13,4 x 4,2 x 4,7 cm</div><div>Weight: approx. 65 g</div><div>Measuring range: Körper: 34-43°C Objekt: 0-100°C</div><div>Accuracy: Stirn: ± 0,3°C (34-43°C) Ohr: ± 0,3°C (35,5-42°C) Objekt: ≤30 ±1.5°C / >30 ±5%</div><div>Batteries: 2 x 1.5V (AAA)</div><div>Article No.: 76140</div><div>EAN: 40 15588 76140 9</div><div>Warranty: 3 Years</div></div>
----	--	---



<div>Infrared Thermometer</div> <div>FTN</div> <div></div> <div></div> <div><ul style="list-style-type: none">• Precise non-contact measurement of the body temperature• Additional function: measurement of ambient temperature, liquids and surfaces• Measurement in a few seconds</div> <div><ul style="list-style-type: none">• Automatic storage of the last 30 measurement results• Visual fever alarm with colour change in the display• Incl. storage box• Incl. battery (2 x 1.5V AAA)• Certified medical device**</div> <div><div>Size: approx. 14,7 x 3,8 x 2,1 cm</div><div>Weight: approx. 48 g</div><div>Measuring range: Körper: 34-42.2°C Objekt: 0-100°C</div><div>Skalierung: 0.1°C</div><div>Battery: 2 x 1.5V (AAA)</div><div>Article No.: 76120</div><div>EAN: 40 15588 76120 1</div><div>Warranty: 3 Years</div></div>	<div>Thermometer</div> <div>TM 700</div> <div></div> <div></div> <div><ul style="list-style-type: none">• Precise measurement of the body temperature: oral, axillary, rectal• Flexible tip• Measurement in 10 seconds• Water resistant</div> <div><ul style="list-style-type: none">• Acoustic fever alarm• Acoustic signal at the beginning and end of temperature reading• Automatic storage of the last measurement result• Automatic power off after 10 minutes• Digital LCD display• Easy battery replacement• Incl. storage box• Incl. battery (1x 1.5V AG3/LR41)• Certified medical device**</div> <div><div>Size: approx. 13,9 x 2,2 x 1,3 cm</div><div>Weight: approx. 11 g</div><div>Measuring range: 32-42.9°C</div><div>Skalierung: 0.1°C</div><div>Battery: 1x 1.5V (AG3/LR41)</div><div>Article No.: 77040</div><div>EAN: 40 15588 77040 1</div><div>Warranty: 3 Years</div></div>	<div>Thermometer</div> <div>FTC</div> <div></div> <div></div> <div><ul style="list-style-type: none">• Precise measurement of the body temperature: oral, axillary, rectal• Acoustic signal at the beginning and end of temperature reading• Water resistant</div> <div><ul style="list-style-type: none">• Acoustic fever alarm• Automatic storage of the last measurement result• Automatic power off after 8 minutes• Digital LCD display• Easy battery replacement• Incl. storage box• Incl. battery (1x 1.5V AG3/LR41)• Certified medical device**</div> <div><div>Size: approx. 12,8 x 1 x 1,9 cm</div><div>Weight: approx. 10 g</div><div>Measuring range: 32-42°C</div><div>Skalierung: 0.1°C</div><div>Battery: 1x 1.5V (AG3/LR41)</div><div>Article No.: 77030</div><div>EAN: 40 15588 77030 2</div><div>Warranty: 3 Years</div></div>
---	--	---

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3

** Certification according to MDD (Medical Device Directive)



Paul, 34
Monitors his performance.

Tom, 27
Strictly monitoring his calorie consumption.

Alina, 31
Focuses on her goals.

Without a goal,
there won't be any
measurable progress.

Sport products – Made for Life.

Active sport is more than just a hobby. Futurologists refer to sport as the main source of self-awareness. The desire to move, the feeling of wanting to do something for the own health, has fundamentally changed the understanding of the work-life-balance. A new determination to stay fit and healthy through sports, has succeeded randomly practiced leisure sports. MEDISANA utilizes innovative fitness and medical devices to promote this new ambition.

The gentle form of regeneration.

The pleasant way to stretch before and after exercising, and to recover.

Everyone who loves their sport knows: stretch properly before the sporting activity, and massage away the tension afterwards. The PowerRoll will promote your desire to stretch before and massage after the sports activity. With four gentle intensity levels it is

possible to reach even deeper muscles. The warm-up is accelerated, and tensions effectively dissolved. The rechargeable, vibrating foam roller is the choice of many professional athletes.



Svenja Thoes, 24 Triathlon athlete
Knows what effectively helps against tensions.



The comprehensive path
to a better fitness:
Activation. Regeneration.
Phase training. Massage.

NEW Foam Roller PowerRoll Soft

- With depth vibration – also reaches the deeper muscles
- 4 Intensity levels, adapted to female needs
- 2 rechargeable Lithium-Ion batteries
- Includes training DVD and exercise plan



With deep-reaching
vibration

Jeanne, 28
Faster recovery
after intense sports



Faster ready for more.

NEW Foam Roller
PowerRoll Pro

- With depth vibration – also reaches the deeper muscles
- 4 intensity levels
- 2 rechargeable Lithium-Ion batteries
- Includes training DVD & training plan



Frank, 37
In best form.



Fitness recovery to the best possible degree.

The foam roller with deep-reaching vibration accelerates the warm-up, and quickly relieves tension after sports activities sustainably. The special highlight: The combination of pressure massage through the own body weight with a depth-effective vibration massage.



POWER VIBRATION
BY PRESSING THE BUTTON



The 1st team of HaSpo Bayreuth knows the effective secret: Effective warm up, faster fitness recovery.

Get rid of many tension.

One for many cases. The VarioRoll for the activation and regeneration.

Every body is different. Luckily, the VarioRoll can be adjusted through a clearance range from 20 to 24 cm between the roller balls, for a perfect adaption to the different muscle and body areas. Tensions, trigger points and adhesions are resolved faster and easier.



NEW
VarioRoll

- Particular suitable for punctual self-massage – for example, back, neck etc.
- For the targeted application of fascia, trigger points and muscles
- Gentle to the spine thanks to the adjustment function
- Incl. training plan




Celine, 19
Sees success.

Counts every move. Because every move counts.


Life is movement. The Activity and Sleep Tracker ViFit touch will help you become more aware of your daily activities - and will help you to better manage, analyse and optimise them.

ViFit touch shows what kind of movement you are capable of. Five parameters are measured in the activity mode, and two in sleep mode. The highly sensitive motion sensor accurately responds only to the steps actually travelled. Another highlight: ViFit touch displays incoming calls and messages.



Simple operation via touch display.
Good overview through luminous digits.

More control and visibility through data transfer via Bluetooth Smart to the VitaDock+ App.








Leonie, 18
Maintains full control during sports.










Activity Tracker ViFit touch




- Activity and Sleep Tracker in one
- Displays incoming calls and messages
- Records the number of steps, distance covered, burned calories, activity time in hours and daily target in %
- Records the movement activity during sleep and the sleep duration
- Movement memory






Foam Roller POWERROLL PRO	Foam Roller POWERROLL SOFT
	
<div><div></div><div></div></div>	<div><div></div><div></div></div>
<div><div></div><div>Q4/2016</div></div>	
<div><div><ul style="list-style-type: none">• PowerRoll with depth vibration - also reaches the deeper muscles• For activation and regeneration as well as for strength and stability exercises and fasciae training• 4 intensity levels• Simple One-touch operation</div></div>	<div><div><ul style="list-style-type: none">• PowerRoll with depth vibration - also reaches the deeper muscles• For activation and regeneration as well as for strength and stability exercises and fasciae training• 4 intensity levels• Simple One-touch operation</div></div>
<div><div><ul style="list-style-type: none">• LED indication of the selected intensity level• Surface structure supports the massage action• Surface is resistant to moisture• 2 integrated rechargeable Lithium-Ion batteries (charger included)• Max. operating time with fully charged battery: ca. 180 min. (depending on the selected intensity level)• Auto-Stop after 10 minutes of operation</div></div>	<div><div><ul style="list-style-type: none">• LED indication of the selected intensity level• Surface structure supports the massage action• Surface is resistant to moisture• 2 integrated rechargeable Lithium-Ion batteries (charger included)• Max. operating time with fully charged battery: ca. 260 min. (depending on the selected intensity level)• Auto-Stop after 10 minutes of operation</div></div>
<div><div></div><div>Includes training DVD and training plan</div></div>	<div><div></div><div>Includes training DVD and training plan</div></div>
<div><div><div><div>Size:</div><div>Weight:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div><div><div>approx. 31 cm, ø 15 cm</div><div>approx. 0.75 kg</div><div>Lithium-ion battery</div><div>79466</div><div>40 15588 79466 7</div><div>3 Years</div></div></div></div>	<div><div><div><div>Size:</div><div>Weight:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div><div><div>approx. 31 cm, ø 15 cm</div><div>approx. 0.75 kg</div><div>Lithium-ion battery</div><div>79518</div><div>40 15588 79518 3</div><div>3 Years</div></div></div></div>











Foam Roller SOLIDROLL	Foam Roller VARIOROLL	Fitness Mat FITNESSMAT
		
<div><div></div><div></div></div>	<div><div></div><div></div></div>	
<div><div></div><div>Q4/2016</div></div>		
<div><div><ul style="list-style-type: none">• For activation and regeneration as well as for strength and stability exercises• Ideal for functional fasciae training• Surface structure supports the massage action• Filled core to withstand heavy loads</div></div>	<div><div><ul style="list-style-type: none">• Particular suitable for punctual self-massage – for example, of the back, neck etc.• For the targeted application of fascia, trigger points and muscles• Especially gentle to the spine thanks to a pull-out mechanism• For extended applications such as the arms, legs etc.</div></div>	<div><div><ul style="list-style-type: none">• Fitness mat for MEDISANA massage roller• Impact-cushioning (shock-absorbent), soft surface• Slightly corrugated surface</div></div>
<div><div><ul style="list-style-type: none">• Helps prevent muscle tension• Resistant to moisture• Includes a training plan• Designed in Germany</div></div>	<div><div><ul style="list-style-type: none">• For activation and regeneration• Resistant to moisture• Designed in Germany</div></div>	<div><div><ul style="list-style-type: none">• Including a strap for easy transport• Surface resistant to moisture - hygienic to wipe</div></div>
<div><div></div><div>Includes training plan</div></div>	<div><div></div><div>Includes training plan</div></div>	
<div><div><div><div>Size:</div><div>Weight:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div><div><div>approx. 31 cm, ø 15 cm</div><div>approx. 0.14 kg</div><div>79467</div><div>40 15588 79467 4</div><div>3 Years</div></div></div></div>	<div><div><div><div>Size:</div><div>Weight:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div><div><div>approx. 20 x 10 x 10 cm, ø 15 cm</div><div>approx. 0.13 kg</div><div>79516</div><div>40 15588 79516 9</div><div>3 Years</div></div></div></div>	<div><div><div><div>Size:</div><div>Weight:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div><div><div>approx. 60 x 190 x 1.5 cm</div><div>approx. 0.85 kg</div><div>79468</div><div>40 15588 79468 1</div><div>3 Years</div></div></div></div>





Activity Tracker VIFIT TOUCH	
	
	
	
<ul style="list-style-type: none">• Activity and Sleep Tracker in one• High-quality OLED display• Displays incoming calls or messages• Individual daily goals and movement reminder• Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*	
<ul style="list-style-type: none">• Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in %• Sleep mode: Tracks activity during sleep and hours slept• Highly sensitive motion sensor for precise activity measurement: reacts to steps (not vibrancy)• Splashproof• Vibration alarm clock• Battery life: 5 - 7 days• Rechargeable lithium-ion battery• Displays date, time and battery• Incl. USB cable	
<div><div>Size:</div><div>Weight:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 23 x 2 x 0.9 cm</div><div>approx. 22 g</div><div>Lithium-ion battery</div><div>79486</div><div>40 15588 79486 5</div><div>3 Years</div></div>	

Activity Tracker VIFIT TOUCH	
	
	
	
<ul style="list-style-type: none">• Activity and Sleep Tracker in one• High-quality OLED display• Displays incoming calls or messages• Individual daily goals and movement reminder• Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*	
<ul style="list-style-type: none">• Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in %• Sleep mode: Tracks activity during sleep and hours slept• Highly sensitive motion sensor for precise activity measurement: reacts to steps (not vibrancy)• Splashproof• Vibration alarm clock• Battery life: 5 - 7 days• Rechargeable lithium-ion battery• Displays date, time and battery• Incl. USB cable	
<div><div>Size:</div><div>Weight:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 23 x 2 x 0.9 cm</div><div>approx. 22 g</div><div>Lithium-ion battery</div><div>79487</div><div>40 15588 79487 2</div><div>3 Years</div></div>	

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3

Activity Tracker VIFIT CONNECT MX3	
	
	
	
<ul style="list-style-type: none">• Activity and Sleep Tracker in one• High-quality OLED display• Easy 1-button operation• Stores the results of the last 7 days• Flight mode (Bluetooth deactivated)• Splashproof• Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*	
<ul style="list-style-type: none">• Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in %• Sleep mode: Tracks movement activity during sleep and duration of sleep• Motion sensor for activity measurement: counts actual steps• Vibration alarm clock• Scratch and abrasion resistant housing• Battery life: up to 7 days• Rechargeable lithium-ion battery• Display of time, date and battery indicator	
	
Incl. 3 wristband in 2 lengths (black, aqua, coral)	
<div><div>Size:</div><div>Weight:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 2 x 1 x 21.5 cm (short wristband)</div><div>approx. 2 x 1 x 25.2 cm (long wristband)</div><div>approx. 22 g (Incl. wristband)</div><div>Lithium-ion battery</div><div>79790</div><div>40 15588 79790 3</div><div>3 Years</div></div>	

Activity Tracker VIFIT CONNECT	
	
	
	
<ul style="list-style-type: none">• Activity and Sleep Tracker in one• High-quality OLED display• Easy 1-button operation• Stores the results of the last 15 days• Modern and attractive design: can be worn discreetly in your trouser pocket or shirt pocket• Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*	
<ul style="list-style-type: none">• Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in %• Sleep mode: Tracks movement activity during sleep and duration of sleep• Highly sensitive motion sensor for precise activity measurement: reacts to steps (not vibrancy)• Scratch and abrasion resistant housing• Rechargeable lithium-ion battery• Battery life: up to 5-7 days• Time and battery indicator• Incl. USB cable	
	
Incl. wristband	
<div><div>Size:</div><div>Weight:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 4.6 x 2,3 x 1 cm</div><div>approx. 30 g (Incl. Armband)</div><div>Lithium-ion battery</div><div>79415</div><div>40 15588 79415 5</div><div>3 Years</div></div>	

Activity Tracker VIFIT	
	
	
	
<ul style="list-style-type: none">• Activity and Sleep Tracker in one• High-quality OLED display• Easy 1-button operation• Stores the results of the last 15 days• Modern and attractive design: can be worn discreetly in your trouser pocket or shirt pocket	
<ul style="list-style-type: none">• Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in %• Sleep mode: Tracks movement activity during sleep and duration of sleep• Highly sensitive motion sensor for precise activity measurement: reacts to steps (not vibrancy)• Scratch and abrasion resistant housing• Rechargeable lithium-ion battery• Battery life: up to 5-7 days• Time and battery indicator• Incl. USB cable	
	
Incl. wristband for the night	
<div><div>Size:</div><div>Weight:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 5.8 cm x 2.1 cm x 1.5 cm</div><div>approx. 12 g</div><div>Lithium-ion battery</div><div>79410</div><div>40 15588 79410 0</div><div>3 Years</div></div>	

Mary, 67
Likes it warm and cozy.

Werner, 72
Feels like 35.

Enjoy life,
totally relaxed.

Sophisticated wellness products – Made for Life.

Today, health is understood in a very diverse manner by different people. Health is more than just being free of illness. More and more people perceive personal health satisfaction as a high value - a physical and mental state that people desire and appreciate. People are willing invest an effort to promote their health, but they also want to enjoy life. Wellness is one of these aspects, and MEDISANA supports this basic feeling of health satisfaction with innovative products that specifically promote the individual well-being.

Releases tension from the body and provides relaxation.

The good feeling when tensions and strains are released. The Shiatsu massage seat cover MCN pro gives you exactly that feeling – with a precise massage.

What would you like today? The entire back, the neck, or optionally the upper or lower back? With red light and heat? With vibration in the seat? The intensive Shiatsu massage provides rapid and lasting well-being.



Maintaining agility and feeling completely well.

Martin, 42
Deeply relaxed like never before.



NEW
Shiatsu Massage Seat Cover MCN Pro


- 3 selectable massage zones: total, upper or lower back
- Additional red light and heat function
- Removable washable neck cover
- Individually adjustable neck massage
- Vibration massage in the seat

ADJUSTABLE
**NECK
MASSAGE**



The small neck tension relief for in-between and on-the-go.


Relaxation after a stressful day. Release tension to clear the mind again.



NEW Neck Massager Vibration NM 865

- 2-stage vibration massage
- Additional red light and heat function
- AC and battery operation, also ideally suited for travel
- Wireless operation

Susan, 28
Relaxed like never before.



The Neck Massager NM 865 gently releases tension using vibration, red light and heat. The targeted massage acts upon the neck and shoulders. Ideal for the tension release after a long workday at the desk. An additional highlight: The simple wireless operation with an integrated control panel.

A lot of office work? NM 865 relieves tension.



You have the well-being at your fingertips.

NEW Hand Massager HM 858


- Vigorous vibration massage
- Infinitely variable regulation of the massage intensity
- Additional red light and heat function
- Various massage attachments: Acupressure and spot massage



For specific massage:
Acupressure and point massage

HM 858 makes everyone a massage professional. Easy operation for a precise massage.

Certain types of tension are persistent because they can't be reached properly. The hand massager HM 858 supports individual massage. The massage intensity can be regulated variably, and the massage effect is enhanced even further through various massage attachments.



Peter, 63
Always relaxed.



Treat yourself and your partner.

Feel warmth with a new intensity.

Discover warmth more efficiently: The new 4D heat flow technology transports heat faster, better and to where you need it.

The innovative 4D heat flow technology ensures an effective heat generation: Up to 8% higher temperature, 15% faster heat emission and a more than 200% better heat transfer. The matrix sensor system

and the flexible filament provide a uniform heat distribution. A visual and audible temperature control make the HP 625 extremely safe.




Applied heat wherever is a need.

Better heat transmission + 200 %

Fast warm-up time + 15 %

Higher temperature + 8 %



Paul, 57
Loves Margit's warmth.

Margit, 59
Loves cold winter days.

Heat Pad HP 625 comfort

- Innovative 4D heat flow technology – for an effective heat development
- Flexible heating wire and matrix sensor system for an even heat distribution
- 4 temperature levels



4D
Heat Flow
Technology



Shiatsu Massage
Seat Cover with Gel Cushions

MC 830



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Massage heads made of high-quality Technogel® for a softer and healthier massage – made in Germany
- Spares spine and hip bones
- 3 selectable massage zones: total, upper or lower back
- Separate vibration massage in the seat with 3 selectable intensity levels

- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Handy remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 105 x 43 x 8,5 cm
Weight: approx. 3.4 kg
Power supply: 230 V/50 Hz
Power output: 18 Watts
Article No.: 88944
EAN: 40 15588 88944 8
Warranty: 3 Years

Shiatsu Acupressure
Massage Seat Cover

MC 825



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Novel acupressure massage in the waist and thigh area
- Separate, height-adjustable neck massage
- 3 selectable massage zones: total, upper or lower back

- 3 intensity levels of acupressur massage selectable
- Intensive shiatsu massage in the back area
- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 110 x 42.5 x 18 cm
Weight: approx. 5.9 kg
Power supply: 220-240V/50-60 Hz
Power output: 40 Watts
Article No.: 88939
EAN: 40 15588 88939 4
Warranty: 3 Years

Shiatsu Massage Seat Cover

MC 822



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- 3 types of massage in the back area: Shiatsu, rolling and spot massage
- Separate, height-adjustable neck massage
- 3 selectable massage zones: total, upper or lower back
- Separate vibration massage in the seat with 3 selectable intensity levels

- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 118 x 40 x 9 cm
Weight: approx. 6.5 kg
Power supply: 220-240V/50-60 Hz
Power output: 40 Watts
Article No.: 88922
EAN: 40 15588 88922 6
Warranty: 3 Years

Shiatsu Massage Seat Cover

MCN



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- 3 selectable massage zones: total, upper or lower back
- Separate neck massage
- Separate vibration massage in the seat with 3 selectable intensity levels

- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 118 x 40 x 9 cm
Weight: approx. 6.5 kg
Power supply: 220-240V/50-60 Hz
Power output: 48 Watts
Article No.: 88930
EAN: 40 15588 88930 1
Warranty: 3 Years

Shiatsu Massage Seat Cover

MCN PRO



NEW



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Separate adjustable neck massage
- 3 selectable massage zones: total, upper or lower back
- Separate vibration massage in the seat with 3 selectable intensity levels

- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 118 x 40 x 9 cm
Weight: approx. 6.5 kg
Power supply: 220-240V/50-60 Hz
Power output: 48 Watts
Article No.: 88970
EAN: 40 15588 88970 7
Warranty: 3 Years

Shiatsu Massage Seat Cover

RBI



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- 3 selectable massage zones: total, upper or lower back
- Separate vibration massage in the seat with 3 selectable intensity levels

- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Handy remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair


Size: approx. 105 x 45 x 9 cm
Weight: approx. 4.5 kg
Power supply: 230V/50 Hz
Power output: 18 Watts
Article No.: 88914
EAN: 40 15588 88914 1
Warranty: 3 Years

Massage Seat Cover Vibration MC 810	Massage Seat Cover Vibration MCH	Presentation Chair
		
<div>HD</div> <div></div> <div>HEAT</div> <div>9 PROGRAMS</div> <div>3</div> <div></div>	<div>HD</div> <div></div> <div>HEAT</div> <div>5 PROGRAMS</div> <div>3</div> <div></div>	
<div><div>• 4 massage zones in the back and seat area</div><div>• 3 selectable massage intensities</div><div>• 9 selectable massage programmes</div><div>• Intensive heat function with or without massage possible</div></div>	<div><div>• 4 massage zones in the back and seat area</div><div>• 3 selectable massage intensities</div><div>• 5 selectable massage programmes</div><div>• Additional heat function</div></div>	<div><div>• Ergonomic swivel chair for presentation of massage seat covers</div><div>• Optimum geometry for seat back massage cushion</div></div>
<div><div>• 6 vibration motors</div><div>• Auto-Stop after 15 minutes of operation</div><div>• Handy remote control</div><div>• Easy to clean material</div><div>• Very handy and easy to store</div><div>• Also ideal for office use or in the car</div><div>• Incl. car power adapter</div></div>	<div><div>• 5 vibration motors</div><div>• Auto-Stop after 15 minutes of operation</div><div>• Handy remote control</div><div>• Easy to clean material</div><div>• Very handy and easy to store</div><div>• Also ideal for office use or in the car</div><div>• Incl. car power adapter</div></div>	<div><div>• With attachment for remote control</div><div>• Suitable for the massage seat covers by Medisana</div><div>• Easy to assemble</div></div>
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 107 x 46 x 4 cm</div><div>approx. 1.4 kg</div><div>100-240 V / 50-60 Hz, 12 V DC 2500mA</div><div>30 Watts</div><div>88937</div><div>40 15588 88937 0</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 107 x 47 x 3.5 cm</div><div>approx. 1.1 kg</div><div>100-240 V / 50-60 Hz, 12V DC 800mA</div><div>9,6 Watts</div><div>88935</div><div>40 15588 88935 6</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Max. load capacity:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 73 x 66 x 111 cm</div><div>approx. 13.3 kg</div><div>max. 120 kg</div><div>88913</div><div>40 15588 88913 4</div><div>3 Years</div></div>

INFO Shiatsu Massage Seat Cover
MC 822

Forget tension. The MC 822 shiatsu massage cushion relaxes you with three different types of massage.

Incorrect sitting, incorrect movements - there are many reasons for stubborn tensions. You can now counteract tensions in a targeted way by using the MC 822 shiatsu massage cushion. Three types of massage in the back area, neck massage that is adjustable for height and separate vibration massage in the sitting area.

Massage Mat MM 825

<div>HD</div> <div></div> <div>HEAT</div> <div>5 PROGRAMS</div> <div>2</div> <div></div>
<div><div>• 4 selectable massage zones: upper and lower back, upper and lower leg</div><div>• 2 intensity levels for each massage zone</div><div>• Comfortable vibration massage for the whole body</div><div>• Smooth heat function in the back area</div></div>
<div><div>• 5 selectable massage programmes</div><div>• Auto-Stop after 15 minutes of operation</div><div>• Handy remote control</div><div>• Extra soft and cosy fleece cover</div><div>• Cushion in the head area of the mat for lying comfortably</div></div>
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 180 x 62 x 10 cm</div><div>approx. 1.8 kg</div><div>100-240 V / 50-60 Hz, 12V DC 1200mA</div><div>14,4 Watts</div><div>88955</div><div>40 15588 88955 4</div><div>3 Years</div></div>

Shiatsu Neck Massager

NM 860



HD

RED LIGHT + HEAT

INTEGRATED

COVER

• Intensive shiatsu massage in the neck area; also applicable on shoulder, back, abdomen, upper and lower leg

• 4 rotating massage heads

• Massage intensity individually adjustable through pulling the holding straps

• 2 levels of rotation speed selectable

• Additional light and heat function in the upper massage heads

• Simple to use and cordless thanks to the integrated operating panel

• Removable, washable cover over the massage heads

Size:

Weight:

Power supply:

Power output:

Article No.:

EAN:

Warranty:

approx. 58 x 34 x 20 cm

approx. 1.55 kg

100-240 V / 50-60 Hz,
12V DC 1.5A

18 Watts

88942

40 15588 88942 4

3 Years

Neck Massager

NM 865



2

RED LIGHT + HEAT

INTEGRATED

• 2-levelled vibration massage for the neck and shoulders

• Optional red light and warmth function (when connected to the mains)

• Mains and battery powered – perfect for travelling too

• Simple to use and cordless thanks to the integrated operating panel

• High-definition massage: targeted, intensive massage for tangible sense of well-being

Size:

Weight:

Power supply:

Power output:

Article No.:

EAN:

Warranty:

approx. 35.5 x 24 x 12 cm

approx. 0.4 kg

100-240 V / 50-60 Hz,
6V DC 600mA

3.6 Watts

88945

40 15588 88945 5

3 Years

Neck Massager

NM 870



Cosy

INTEGRATED

• Neck Massager made of comfortable molded foam – fits perfectly to the contours of the head and neck

• Pleasant vibration massage for the neck and shoulders

• Ideal for travelling and at home

• With integrated "One Touch"- control panel

• Extra soft and fluffy fleece cover

• With integrated storage bag – Space-saving storage or on the go

• Incl. 2 batteries (1.5 V AA)

Size:

Weight:

Batteries:

Article No.:

EAN:

Warranty:

approx. 28 x 26 x 9 cm

approx. 0.4 kg

2 x 1.5 V (AA)

88946

40 15588 88946 2

3 Years

Neck Massager with Tapping Function

MNT



5

3 PROGRAMS

INTEGRATED

• Vigorous tapping massage

• 5 selectable massage intensities

• 3 selectable massage programs

• Ergonomic shape

• Suitable for neck, shoulders, legs and hips

• 2 massage heads

• Position of the individual tapping elements are easy movable

• Handy remote control

Size:

Weight:

Power supply:

Power output:

Article No.:

EAN:

Warranty:

approx. 61 x 43 x 6 cm

approx. 2.6 kg

230 V, 50 Hz

75 Watts

88940

40 15588 88940 0

3 Years

Neck Massager

MNV



2

RED LIGHT + HEAT

INTEGRATED

• 2-levelled vibration massage for the neck and shoulders

• Optional red light and warmth function (when connected to the mains)

• Mains and battery powered – perfect for travelling too

• Simple to use and cordless thanks to the integrated operating panel

Size:

Weight:

Power supply:

Power output:

Article No.:

EAN:

Warranty:

approx. 34 x 24 x 13.5 cm

approx. 0.8 kg

100-240 V / 50-60 Hz,
6V DC 600mA

3.6 Watts

88941

40 15588 88941 7

3 Years

Product Overview | Wellness

67

Shiatsu Massage Cushion
MC 840



Q4/2016

- Soothing Shiatsu massage for neck, shoulder, back and legs
- 4 rotating massage heads
- Additional red light and heat function
- Particularly flat design for flexible, convenient operation

- Handy cable remote control
- Convenient velcro for fastening at the back of the chair
- Easy to clean material

INFO Handmassagegerät
HM 858

Release tensions at the exact location. The HM 858 hand massage device supports you in this.

What better way to deal with tensions than by tackling them at the source? The HM 858 hand massage device has various add-on massage sets for effective acupressure or point-type massage. You can increase the massage effect by switching on the red light and warmth. The intensity of the massage is continuously variable.

- Size:** approx. 31 x 23 x 11 cm
Weight: approx. 1.3 kg
Power supply: 220-240 V / 50-60 Hz
12V DC 1000mA
Power output: 12 Watts
Article No.: 88949
EAN: 40 15588 88949 3
Warranty: 3 Years

Hand Massager
HM 855



- Comfortable vibration massage
- Additional red light and heat function

- Easy to use thanks to its ergonomic shape
- Light and handy design for convenient self-massage

- Size:** approx. 28 x 6 x 7 cm
Weight: approx. 0.3 kg
Power supply: 220-240 V / 50-60 Hz
Power output: 15 Watts
Article No.: 88297
EAN: 40 15588 88297 5
Warranty: 3 Years

Hand Massager
HM 858



Q1/2017

- Vigorous vibration massage
- Continuously adjustment of massage intensity
- Additional red light and heat function
- Various massage attachments: acupressure and spot massage

- Easy to use thanks to its ergonomic shape
- Non-slip grip
- Extra-long cable



Incl. 2 x 2 massage attachments

- Size:** approx. 32,5 x 17 x 14,3 cm
Weight: approx. 1.8 kg
Power supply: 220-240 V / 50-60 Hz
Power output: 28 Watts
Article No.: 88292
EAN: 40 15588 88292 0
Warranty: 3 Years

Mini-Hand Massager
HM 840



- Comfortable vibration massage
- With removable and washable massage pad made of natural loofah
- Waterproof for use in the shower or when bathing

- Intensive peeling effect when using with loofah pad
- Light and handy design for convenient self-massage
- Incl. 3 batteries (1.5V AAA)



Incl. loofah pad

- Size:** approx. 8 x 8 x 7 cm
Weight: approx. 0.2 kg
Batteries: 3 x 1.5V (AAA)
Article No.: 88521
EAN: 40 15588 88521 1
Warranty: 3 Years

FS 881



- Relaxing spa and vibration massage
- Soothing foot spa
- 3 in 1: water heating, bubble- and vibrating massage
- Intensive reflexology

- Includes massage attachment
- Non-slip rubber feet
- Automatic heat function
- Non-slip rubber feet
- Comfortable Size (suitable up to shoe size 44)

Size: approx. 39 × 33 × 13 cm
Weight: approx. 1.22 kg
Power supply: 220-240V/50Hz
Power output: 60 Watts

Article No.: 88380
EAN: 40 15588 88380 4
Warranty: 3 Years

Foot Spa

FS 883



NEW

Q4/2016

- Relaxing spa and vibration massage
- Pleasant red light area
- 3 in 1: water heating, bubble- and vibrating massage
- Function for keeping the water warm
- 3 exchangeable pedicure attachments (Cornea remover, brush, massage attachment)

- Massage rollers for foot reflexology massage
- Massage-aiding footbed
- Dry massage function
- Non-slip rubber feet
- Removable splash protection
- Comfortable Size (suitable up to shoe size 46)

Size: approx. 41 × 35 × 18 cm
Weight: approx. 1.5 kg
Power supply: 220-240V/50-60Hz
Power output: 80 Watts

Article No.:	88391
EAN:	40 15588 88391 0
Warranty:	3 Years

Foot Spa

FS 885 COMFORT




- Relaxing spa and vibration massage
- 5 in 1: Bubble- and vibrating massage, red light and heat function and massage rollers in the middle of the foot bubble bath
- Water can be heated to the preferred temperature (20°C to 48°C)
- Timer function for 20/30/40/50/60 min.

- Massage rollers in the middle of the foot bubble bath to massage the inner side of the foot
- Digital display for timer setting and water temperature
- Red light and heat function
- Ergonomically arranged air bubble nozzles
- Removable massage rollers in the area of the soles
- Good heat insulation through double-walled casing
- Intensive reflexology
- Comfortable Size (suitable up to shoe size 46)

Size: approx. 38.6 x 24.3 x 45.9 cm
Weight: approx. 2.9 kg
Power supply: 220-240V/50Hz
Power output: 390 Watts

Article No.: 88378
EAN: 40 15588 88378 1
Warranty: 3 Years



Bath Spa

BBS



- Relaxing air bubble bath to stimulate blood circulation and to relieve muscles
- Massage/bubble function with 3 intensity levels
- Intensive heat function selectable
- Extra quiet

- Automatic heat function
- Intensity levels adjustable on the device or via remote control
- Timer function: 10/20/30 minutes
- Non-slip soft mat for lying comfort
- Suitable for any bath tub
- Integrated fragrance dispenser for scented oils
- Double safety valve for air hose and device
- Incl. 3 aroma dispensers
- Incl. remote control

Size:	Soft mat: approx. 120 x 36 cm Air hose: approx. 2.4 m Motor: approx. 31 x 21 x 31 cm
Weight:	approx. 4.3 kg (gesamt)
Power supply:	220-240V/50 Hz
Power output:	Remote control: 3 V CR2032 570 Watts

Article No.:	88386
EAN:	40 15588 88386 6
Warranty:	3 Years

Bath Spa

MBH






- Relaxing air bubble bath to stimulate blood circulation and to relieve muscles
- Massage/bubble function with 3 intensity levels
- Auto-Stop after 20 minutes of operation

- Automatic heat function
- Intensity levels adjustable on the device or via remote control
- Non-slip soft mat for lying comfort
- Suitable for any bath tub
- Integrated fragrance dispenser for scented oils
- Double safety valve for air hose and device
- Incl. 3 aroma dispensers
- Incl. remote control


Size:	Soft mat: approx. 120 x 36 cm Air hose: approx. 2,4 m Motor: approx. 25 x 16 x 17 cm approx. 2,7 kg (gesamt)
Weight:	
Power supply:	220-240V/50Hz
Power output:	Remote control: 3 V CR2032 450 Watts

Article No.:	88375
EAN:	40 15588 88375 0
Warranty:	3 Years

Heat Pad HP 605	Heat Pad for Neck and Back HP 610
	
<div><div><div>Turbo Heat</div><div>90 min AUTO OFF</div><div>WASHABLE</div></div></div>	<div><div><div>Turbo Heat</div><div>90 min AUTO OFF</div><div>WASHABLE</div></div></div>
	
<div><div>• Ultra-soft breathable cushion</div><div>• Safety system with overheating protection</div><div>• Quick heating with 100 watts</div><div>• 4 temperature settings</div></div>	<div><div>• Ultra-soft breathable cushion</div><div>• Safety system with overheating protection</div><div>• Quick heating with 100 watts</div><div>• 4 temperature settings</div></div>
<div><div>• Automatic switch-off after 90 minutes</div><div>• Ergonomically designed controller with LED control display on both sides for right and left-handed persons</div><div>• Removable cover of cotton</div><div>• Cover machine-washable at 30°C</div><div>• Cushion washable thanks to detachable controller</div></div>	<div><div>• Automatic switch-off after 90 minutes</div><div>• Ergonomically designed controller with LED control display on both sides for right and left-handed persons</div><div>• Removable cover of cotton</div><div>• Cover machine-washable at 30°C</div><div>• Cushion washable thanks to detachable controller</div><div>• Perfect fit for the back and neck</div><div>• Convenient velcro strap for individual fastening of the cushion</div></div>
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 40 x 33 cm</div><div>approx. 0.5 kg</div><div>220-240 V / 50 Hz</div><div>100 Watts</div><div>61147</div><div>40 15588 61147 6</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 40 x 60 cm</div><div>approx. 0.5 kg</div><div>220-240 V / 50 Hz</div><div>100 Watts</div><div>61167</div><div>40 15588 61167 4</div><div>3 Years</div></div>

Heat Pad

HP 615



APS
ACTIVE PROTECTION SYSTEM

Turbo Heat

90 min
AUTO OFF


XL SIZE

MICRO FIBRE

WASHABLE

- Ultra-soft breathable cushion
- “APS” technology – active protection system: visual and acoustic temperature control
- Flexible heating wires and a matrix sensor system for even heat distribution
- 4 temperature settings

- Turbo heat with powerful 100 watts
- Automatic switch-off after 90 minutes
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Removable cover of cosy and soft upper material
- Cover machine-washable at 40°C
- Cushion washable thanks to detachable controller



Incl. exchangeable microfibre cover

Size:

Weight:

Power supply:

Power output:

Article No.:

EAN:

Warranty:

approx. 41,5 x 53 cm

approx. 0.5 kg


230V / 50 Hz

100 Watts



61144



40 15588 61144 5

3 Years

Heat Pad for Neck and Shoulder HP 622

<div><div><div>Turbo Heat</div><div>90 min AUTO OFF</div><div>Cosy</div><div>WASHABLE</div></div></div>

<div><div>• Flexible heating wires and a matrix sensor system for even heat distribution</div><div>• 6 temperature settings</div><div>• Ergonomic shape for pleasant warmth at shoulder and neck</div></div>
<div><div>• Turbo heat with powerful 100 watts</div><div>• Automatic shut-off after 90 minutes</div><div>• Handy controller</div><div>• Machine washable at 40°C</div><div>• Oeko-Text® Standard 100</div></div>
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>ca. 56 x 52 cm</div><div>ca. 0.5 kg</div><div>220-230 V / 50 Hz</div><div>100 Watts</div><div>61155</div><div>40 15588 61155 1</div><div>3 Years</div></div>

Heat Pad HP 625 COMFORT

<div><div><div>4D TECHNOLOGY</div><div>APS ACTIVE PROTECTION SYSTEM</div><div>Turbo Heat</div><div>90 min AUTO OFF</div><div>Cosy</div><div>WASHABLE</div></div></div>

<div><div>• Ultra-soft breathable cushion</div><div>• Innovative 4D heat flow technology – for an effective heat development</div><div>• Flexible heating wire and matrix sensor system for an even heat distribution</div><div>• 4 temperature levels</div></div>
<div><div>• “APS” technology – active protection system: visual and acoustic temperature control</div><div>• Turbo heat with powerful 100 watts</div><div>• Automatic switch-off after 90 minutes</div><div>• Ergonomically designed controller with LED control display on both sides for right and left-handed persons</div><div>• Cover machine-washable at 40°C</div><div>• Removable cover of cosy and soft upper material</div><div>• Oeko-Text® Standard 100</div></div>
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 46 x 35 cm</div><div>approx. 0.5 kg</div><div>230V / 50 Hz</div><div>100 Watts</div><div>61140</div><div>40 15588 61140 7</div><div>3 Years</div></div>

Heat Pad for Shoulder and Back HP 630

<div><div><div>4D TECHNOLOGY</div><div>APS ACTIVE PROTECTION SYSTEM</div><div>Turbo Heat</div><div>90 min AUTO OFF</div><div>Cosy</div><div>WASHABLE</div></div></div>

<div><div>• Ultra-soft breathable cushion</div><div>• Innovative 4D heat flow technology – for an effective heat development</div><div>• Flexible heating wire and matrix sensor system for an even heat distribution</div><div>• 4 temperature levels</div></div>
<div><div>• “APS” technology – active protection system: visual and acoustic temperature control</div><div>• Turbo heat with powerful 100 watts</div><div>• Automatic switch-off after 90 minutes</div><div>• Ergonomically designed controller with LED control display on both sides for right and left-handed persons</div><div>• Machine-washable at 40°C</div><div>• Ergonomic shape for pleasant warmth at shoulder and neck</div><div>• Cosy and soft upper material</div><div>• Oeko-Text® Standard 100</div></div>
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 55 x 65 cm</div><div>approx. 0.76 kg</div><div>230V / 50 Hz</div><div>100 Watts</div><div>61157</div><div>40 15588 61157 5</div><div>3 Years</div></div>

Heat Belt
with Rechargeable Batteries

HS 680



Heating Blanket

HDW



- Heat belt with rechargeable batteries for use indoors and outdoors
- Soothing heat for the back and abdomen
- 4 temperature settings
- Maximum freedom of movement – can also be worn under clothing

- Easy to attach thanks to Velcro system
- Additional elastic belt for stabilization
- Overheating protection
- Automatic switch-off after 90 minutes
- Ideal at home and when travelling
- Long-lasting heat thanks to the powerful lithium-ion battery – min. 3 hours at max. temperature
- Machine washable at 30°C

Size: approx. 125 x 20.5 cm
Weight: approx. 0.5 kg
Kapazität: 7.2V DC 2200 mAh
Article No.: 61160
EAN: 40 15588 61160 5
Warranty: 3 Years

Size: approx. 180 x 130 cm
Weight: approx. 1.7 kg
Power supply: 220-240V / 50 Hz
Power output: 120 Watts
Article No.: 60227
EAN: 40 15588 60227 6
Warranty: 3 Years

Foot Warmer

FWS



- Extra-soft breathable core
- Safety system with overheating protection
- Quick heating with 100 watts
- 3 temperature settings
- Large foot space (suitable up to shoe size 46)

- Automatic switch-off after 90 minutes
- Ergonomically designed controller with LED control display
- Removable fluffy soft lining
- Lining machine washable at 30°C
- Oeko-Tex® Standard 100

Size: approx. 30 x 30 x 23 cm
Weight: approx. 0.6 kg
Power supply: 220-240V / 50 Hz
Power output: 100 Watts
Article No.: 60257
EAN: 40 15588 60257 3
Warranty: 3 Years

Heated Underblanket

HU 650



- Convenient heated underblanket with stretch function for easy attachment to the mattress
- Safety system with overheating protection
- 4 temperature settings
- Quick heating with 100 watts

- Two warming zones for body and feet
- Automatic switch-off after 180 minutes
- Ergonomically designed controller with LED control display
- Washable thanks to detachable controller
- Upper material of high-quality cotton
- Suitable for mattress sizes 200 x 90-100 cm

Size: approx. 200 x 90-100 cm
Weight: approx. 1,4 kg
Power supply: 230V / 50 Hz
Power output: 100 Watts
Article No.: 61210
EAN: 40 15588 61210 7
Warranty: 3 Years

Heated Underblanket

HU 662



- 6 temperature settings
- Safety system with overheating protection
- 6 temperature settings
- Quick heating with 100 watts

- Two warming zones for body and feet
- Automatic switch-off after 180 minutes
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Washable thanks to detachable controller
- Pleasantly soft upper material
- Suitable for mattress sizes 200 x 90-100 cm
- Oeko-Tex® Standard 100

Size: approx. 150 x 80 cm
Weight: approx. 1 kg
Power supply: 220-230V / 50 Hz
Power output: 100 Watts
Article No.: 61220
EAN: 40 15588 61220 6
Warranty: 3 Years

Heated Underblanket


HU 665



- Safety system with overheating protection
- 3 temperature settings
- Washable thanks to detachable controller

- Pleasantly soft upper material
- Automatic switch-off after 180 minutes
- Handy controller
- Suitable for mattress sizes 200 x 90-100 cm
- Oeko-Tex® Standard 100

Size: approx. 150 x 80 cm
Weight: approx. 1 kg
Power supply: 230V / 50 Hz
Power output: 60 Watts
Article No.: 60217
EAN: 40 15588 60217 7
Warranty: 3 Years



Beate, 34
No longer feels tired
during dark days.

Paula, 6
Has a sensitive lung.

People who actively do something for their health, feel better. Today, people place their focus on the individual quality of life, and the feeling that they are actively in control of their own health. This also applies for the therapeutic area. The aim is not to replace the physician. But many people are realizing that they can do a lot for their own prosperity and well-being. MEDISANA supports the desire of people to prudently treat their own health, with therapeutic products which provide a gentle and specific assistance.

Good. To help
yourself better.

Effective wellness products – Made for Life.

Dominic, 32
Always in a good mood.

Lea, 26
Feels well.



NEW

Daylight lamp LT 460

- Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- Uniform light distribution with LED technology
- 2 light intensities can be selected
- Certified medical device*



Light up for a good mood.

People shouldn't be depressed by gloomy weather or long winter days.

Energising daylight applied in a specific manner lifts the mood, and acts against the so-called „winter blues“. The two daylight lamps LT 460 and LT 470 are able to effectively replenish the positive mood.

NEW

Daylight Lamp LT 470

- Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- Light intensity up to 10,000 Lux
- Certified medical device*



A good mood is independent
of the weather.

* Certified according to MDD (Medical Device Directive)

Victoria, 37
Starts the day in a good mood.



NEW Daylight Lamp LT 480

- Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- Flicker-free
- Certified medical device*


Certified Medical Device


UV free


LUX 10.000


3 Years Warranty

Despite the wind and weather – the day will turn out good.

The day can take its course. I have ensured my good mood.

Daylight lamps simulate the light spectrum of daylight or sunlight. Weather-sensitive people often suffer during grey, cloudy weather. The daylight lamp LED 480 specifically counteracts this condition by simulating the sunlight spectrum. The body reacts to this and releases „activity hormones“.



The Certified medical device works with flicker-free light.

Improves breathing instantly.


Regular inhalations help people with sensitive lungs and bronchi, especially during the cold and wet season.

With the inhalator IN 525, daily inhalation will no longer be an unpleasant duty, but rather a comforting support for health issues regarding the upper and lower respiratory tract. The inhalator is particularly effective through to high respirable aerosol content. Additionally, it is very quiet and turns itself off when the active ingredient tank is empty.



Ultrasonic Inhalator IN 525

- Targeted treatment of diseases of the upper and lower respiratory tract
- Use for a short time in the supine position is also possible
- Automatic switch-off with an empty solution tank
- Certified medical device*



The IN 525 is an indispensable aid for people who suffer from allergies.



It's far easier to laugh with a clear lung.

* Certified according to MDD (Medical Device Directive)

Awakened by the sun.



Pete, 41
Nowadays likes to get up.

Sophie, 33
Likes it romantic.

NEW

Wake-Up Light WL 444

- Natural awakening for a soft start in the day
- Integrated FM Radio
- Long-lasting, energy-saving LED lighting





Full of energy early in the morning.



Relax and recharge your batteries with different colours.

The right light for a good sleep.

Fall sleep and wake up with the proper light?
The Wake-Up Light WL 460 makes it possible.

The proper light has a strong impact on the well-being. The Wake-Up Light WL 460 will gently wake you from the deep sleep phase with a simulated sunrise, 30 minutes prior to the wake up time. And it facilitates the natural falling asleep with red light. Eight different natural sounds can be played.



NEW

Wake-Up Light WL 460

- Natural awakening for a soft start in the day
- Natural falling asleep through red light
- Long-life, energy-saving LED lighting
- Integrated Radio
- 2 alarms can be adjusted



Natural falling asleep with red light
Because red light promotes the release of sleep hormones.

Daylight Lamp LT 460	Daylight Lamp LT 470	Daylight Lamp LT 480
 <div>NEW</div>	 <div>NEW</div>	 <div>NEW</div>
<div><div><div><div><div></div><div>LUX 10.000</div></div><div><div></div><div></div></div><div><div></div><div></div></div><div><div></div><div>LED</div></div><div><div></div><div>MDD 93/42</div></div></div></div></div>	<div><div><div><div><div></div><div>LUX 10.000</div></div><div><div></div><div></div></div><div><div></div><div>UV FREE</div></div><div><div></div><div></div></div><div><div></div><div>MDD 93/42</div></div></div></div></div>	<div><div><div><div><div></div><div>LUX 10.000</div></div><div><div></div><div></div></div><div><div></div><div>UV FREE</div></div><div><div></div><div></div></div><div><div></div><div>MDD 93/42</div></div></div></div></div>
Q4/2016	Q4/2016	Q4/2016
<div><div><div>• Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression</div><div>• Uniform light distribution with LED technology</div><div>• 2 light intensities can be selected</div></div></div>	<div><div><div>• Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression</div><div>• UV Filter</div><div>• Flicker-free</div></div></div>	<div><div><div>• Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression</div><div>• UV Filter</div><div>• Flicker-free</div></div></div>
<div><div><div>• Light intensity: up to 10,000 Lux</div><div>• Comfortable one-button operation</div><div>• Can be used horizontally and vertically</div><div>• Ideal for the home and on the road</div><div>• LED technology</div><div>• Including storage bag</div><div>• Certified medical device*</div></div></div>	<div><div><div>• Light intensity: up to 10,000 Lux</div><div>• Comfortable one-button operation</div><div>• Particularly silent</div><div>• Folding base (wall mounting possible)</div><div>• Energy-saving 45W lamp</div><div>• Certified medical device*</div></div></div>	<div><div><div>• Light intensity: up to 10,000 Lux</div><div>• Comfortable one-button operation</div><div>• Particularly silent</div><div>• Folding base (wall mounting possible)</div><div>• Two 36 W lamps</div><div>• Certified medical device*</div></div></div>
<div><div><div>Size: approx. 26 x 15 x 1.7 cm</div><div>Weight: approx. 0.5 kg</div><div>Leuchtfäche: 20.5 x 12 cm</div><div>Power supply: Eingang: 100-240V~, 50/60 Hz, 500mA</div><div>Ausgang: 12 V - 1 A</div><div>Power output: 12 Watts</div><div>Cable length: approx. 1.8 m</div><div>Article No.: 45220</div><div>EAN: 40 15588 45220 8</div><div>Warranty: 3 Years</div></div></div>	<div><div><div>Size: approx. 27.3 x 20.7 x 7.4 cm</div><div>Weight: approx. 0.8 kg</div><div>Leuchtfäche: 20 x 21 cm</div><div>Power supply: 220-240V~, 50 Hz</div><div>Power output: 45 Watts</div><div>Cable length: approx. 1.3 m</div><div>Article No.: 45222</div><div>EAN: 40 15588 45222 2</div><div>Warranty: 3 Years</div></div></div>	<div><div><div>Size: approx. 48 x 34 x 10 cm</div><div>Weight: approx. 2 kg</div><div>Leuchtfäche: 31 x 45 cm</div><div>Power supply: 220-240V~, 50 Hz</div><div>Power output: 72 Watts</div><div>Cable length: approx. 1.5 m</div><div>Article No.: 45224</div><div>EAN: 40 15588 45224 6</div><div>Warranty: 3 Years</div></div></div>



Wake Up Light WL 444	Wake Up Light WL 460
 <div>NEW</div>	 <div>NEW</div>
<div><div><div><div><div></div><div>SUNRISE</div></div><div><div></div><div>8x</div></div><div><div></div><div>WELLNESS LIGHT</div></div><div><div></div><div>LED</div></div><div><div></div><div></div></div></div></div></div>	<div><div><div><div><div></div><div>SUNRISE</div></div><div><div></div><div>8x</div></div><div><div></div><div>WELLNESS LIGHT</div></div><div><div></div><div>LED</div></div><div><div></div><div></div></div></div></div></div>
Q4/2016	Q1/2017
<div><div><div>• Natural awakening for a soft start in the day</div><div>• Alarm function through an integrated FM radio or 8 natural sounds</div><div>• Long-lasting, energy-saving LED lighting</div></div></div>	<div><div><div>• Natural awakening for a soft start in the day</div><div>• Natural falling asleep through red light* (Timer can be set: 30, 60 or 90 min)</div><div>• Integrated radio</div><div>• Long-lasting, energy-saving LED lighting</div></div></div>
<div><div><div>• Increasing light intensity: 30 minutes before alarm</div><div>• 2 alarms can be adjusted</div><div>• Touch button: snooze function</div><div>• Backup battery</div><div>• 8 radio-memory slots</div><div>• Light intensity can be adjusted</div><div>• Colour-changing wellness light with 7 colours</div><div>• Blue LCD display with backlight</div><div>• Incl. antenna</div></div></div>	<div><div><div>• Increasing light intensity: 30 minutes before alarm</div><div>• 2 alarms can be adjusted</div><div>• Snooze function and zizz button (extra timer for alarm after 15 to 90 minutes)</div><div>• 8 natural sounds with increasing volume: tweeting birds, croaking frogs, running brook, rushing waves, crackling fireplace, falling rain, chirping crickets, Asian wooden gong</div><div>• Colour-changing wellness light with 7 colours</div><div>• Large, easy-to-read display with dimmer function (adjustable light + dark)</div><div>• MP3 player connection is possible</div><div>• Incl. antenna</div></div></div>
<div><div><div>Size: approx. 17.5 x 17.5 x 14.5 cm</div><div>Weight: approx. 0.37 kg</div><div>Power supply: Netzadapter</div><div>Eingang: 100-240V~, 50/60Hz, 350mA</div><div>Ausgang: 12V = 1000mA</div><div>Article No.: 45110</div><div>EAN: 40 15588 45110 2</div><div>Warranty: 3 Years</div></div></div>	<div><div><div>Size: approx. 19 x 14.5 x 27.5 cm</div><div>Weight: approx. 0.9 kg</div><div>Power supply: 100-240 V/ 50-60 Hz, 12V DC 1A</div><div>Power output: 12 Watts</div><div>Cable length: approx. 2 m</div><div>Article No.: 45115</div><div>EAN: 40 15588 45115 7</div><div>Warranty: 3 Years</div></div></div>



* Certification according to MDD (Medical Device Directive) ** Red light stimulates the discharge of sleep hormones and makes it easier to fall asleep

Inhalator
IN 500 COMPACT



- Nebulisation with compressed air technology
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories: Inhalation with mask, mouthpiece or nose adapter – ideal for children

- Highly efficient due to very precise nebulisation
- Particularly effective thanks to high proportion of respirable aerosols
- Extra long hose – 2 m
- Certified medical device*



Incl. mouthpiece, masks and nose adapters for children and adults, connecting piece, 5 replacement filters

Size:	approx. 19.5 x 13.5 x 19.2 cm
Weight:	approx. 1.3 kg
Power supply:	230V~/50 Hz
Power output:	60 Watts
Filling volume:	max. 10 ml
Nebulizing output:	≤ 0.2 ml/min
Particle size:	>60% ≤ 5 µm
Hose length:	approx. 2 m
Article No.:	54520
EAN:	40 15588 54520 7
Warranty:	3 Years

Inhalator
IN 550 PRO



- Nebulisation with compressed air technology
- Including 3 attachments to produce different particle sizes (3,5/4,0/5,0 µm)
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories: Inhalation with mask, mouthpiece or nose adapter – ideal for children

- Highly efficient due to very precise nebulisation
- Particularly effective thanks to high proportion of respirable aerosols
- Extra long hose – 2 m
- Integrated accessory tray
- Certified medical device*



Incl. mouthpiece, masks and nose adapters for children and adults, connecting piece, 3 attachments to produce different particle sizes, 5 replacement filters, 2 replacement fuses

Size:	approx. 33 x 20 x 13.8 cm
Weight:	approx. 2 kg
Power supply:	230V~/50 Hz
Power output:	60 Watts
Filling volume:	max. 10 ml
Nebulizing output:	≤ 0.2 ml/min
Particle size:	4,5~6,0 µm (Mouth and throat) 3,5~4,5 µm (For respiratory tract and lungs front) 3,0~4,0 µm (For the deeper areas of the lung)
Hose length:	approx. 2 m
Article No.:	54530
EAN:	40 15588 54530 6
Warranty:	3 Years

Ultrasonic Inhalator
USC



- High efficiency through innovative micro-membrane nebulisation (mesh technology) with ultrafine droplets
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories

- Particularly effective thanks to high proportion of respirable aerosols
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulisation
- Incl. storage case
- Incl. 4 batteries (1.5 V AA)
- Certified medical device*



Incl. mouthpiece, masks for children and adults and connecting piece

Size:	approx. 14.2 x 7.5 x 7.5 cm
Weight:	approx. 0.3 kg (without batt.)
Batteries:	4 x 1.5 V (AA)
Filling volume:	max. 8 ml
Nebulizing output:	≤ 0.4 ml/min
Particle size:	Average <4,8 µm
Article No.:	54105
EAN:	40 15588 54105 6
Warranty:	3 Years

Ultrasonic Inhalator
IN 525



NEW



- High efficiency through innovative micro-membrane nebulisation (mesh technology) with ultrafine droplets
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories
- Use for a short time in the supine position is also possible

- Particularly effective thanks to high proportion of respirable aerosols
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulisation
- Low power consumption (no more than 0.9 Watts)
- Automatic switch-off with an empty solution tank
- Incl. storage case
- Incl. 2 batteries (1.5 V AA)
- Certified medical device*



Incl. mouthpiece and masks for children and adults

Size:	approx. 7.2 x 4.0 x 10.8 cm
Weight:	approx. 0.1 kg (without batt.)
Batteries:	2 x 1.5 V (AA)
Power output:	0,9 Watts
Filling volume:	max. 8 ml
Nebulizing output:	≤ 0.2 ml/min
Particle size:	Average <5 µm
Article No.:	54115
EAN:	40 15588 54115 5
Warranty:	3 Years

Q1/2017

* Certification according to MDD (Medical Device Directive)

Electronic Body Toner
BT 850



- Stimulates blood flow in the area of use
- Relaxes tense muscles
- Helps to build up muscles
- Compatible with: Bluetooth Smart 4.0 from iPhone 4s and Android 4.3 onwards



- 5 different application programmes for the waist, shoulder and neck, back, limbs
- 7 intensity levels to choose from
- Different user accounts can be saved with individual settings (name, age, height, weight, gender)
- With timer function
- Including 4 self-adhesive electrodes, storage box and replacement battery
- Only works in conjunction with the app
- Incl. storage box
- Incl. 2 batteries (CR2032)



Incl. 4 self-adhesive electrodes

Size:	Box: 12 x 6.5 x 2.5 cm Electrodes: approx. ø 40mm
Weight:	approx. 25 g (without storage box)
Batteries:	1 x 3 V (CR2032)
Cable length:	approx. 40 cm
Article No.:	88330
EAN:	40 15588 88330 9
Warranty:	3 Years

Electric Lice Comb	Abdominal Muscle Stimulator	Abdominal Muscle Stimulator
LCS	BOB	AM 880
		
		
<ul style="list-style-type: none">• Simple combing frees hair from head lice and nits• No side effects, and no use of chemical agents• Also suitable for dogs and cats	<ul style="list-style-type: none">• Flexible abdominal belt with inbuilt electrodes• 10 pre-installed programmes• Intensity can be set over 25 levels• Can be worn under clothing thanks to the flat design	<div>Q1/2017</div> <ul style="list-style-type: none">• Flexible Neoprene waist belt with incorporated electrodes• Choice of 6 abdominal muscle training programmes of 30 minutes each• Intensity can be set over 40 levels• Can be worn under clothing thanks to the flat design
<ul style="list-style-type: none">• The lice are destroyed by an electrical current upon contact with the teeth of the comb• Incl. batterie (1.5 V AA)	<ul style="list-style-type: none">• Automatic switch-off at the end of the programme• Simple operation with a large LCD display• Device is integrated into the belt• 1 channel with 2 electrodes• Waist belt with Velcro closer• Incl. 3 batteries (1.5 V AAA)• Certified medical device*	<ul style="list-style-type: none">• Automatic switch-off at the end of the programme• Simple operation with a large LCD display• Device is integrated into the belt• 1 channel with 4 electrodes• Waist belt with Velcro closer• Water contact electrodes, no need for contact gel• Incl. 3 batteries (1.5 V AAA)
Size: approx. 17.3 x 6.3 x 2.3 cm Weight: approx. 75 g Power supply: 1 x 1.5 V (AA) Article No.: 41015 EAN: 40 15588 41015 4 Warranty: 3 Years	Size: approx. 110 x 14.4 x 2.2 cm Weight: approx. 0.23 kg Batteries: 3 x 1.5 V (AAA) Channels: 1 Curve form: Bi-phase square wave pulse Intensity: 0–90mA Frequency: 3–100 Hz Pulse width: 150–200 µs Article No.: 88320 EAN: 40 15588 88320 0 Warranty: 3 Years	Size: approx. 110 x 16 x 3 cm Weight: approx. 0.31 kg Batteries: 3 x 1.5 V (AAA) Channels: 1 Curve form: Bi-phase square wave pulse Intensity: 0–70mA Frequency: 30–80 Hz Pulse width: 200–300 µs Article No.: 88322 EAN: 40 15588 88322 4 Warranty: 3 Years

Infrared Lamp	Infrared Lamp	Infrared Lamp
IR 885	IRH	IRL
		
		
<ul style="list-style-type: none">• Intense infrared light to relax muscles or treat colds• Stimulates blood circulation• Relieves pain• Durable lamp (halogen)	<ul style="list-style-type: none">• Intense infrared light to relax muscles or treat colds• Stimulates blood circulation• Relieves pain• Durable lamp (Philips)	 <ul style="list-style-type: none">• Intense infrared light to relax muscles or treat colds• Stimulates blood circulation• Relieves pain• Durable lamp (Philips)
<ul style="list-style-type: none">• Adjustable tilt angle (35°)• On/Off button• Timer (15 to 1 minute)• Fan caster for cooling• Cable rewind• Carrying handle• Large-surface application• Incl. lamp (halogen)• Certified medical device*	<ul style="list-style-type: none">• Adjustable tilt angle• Made in EU• Certified medical device*	<ul style="list-style-type: none">• Effective heat therapy through pressed glass bulb• Adjustable tilt angle• Made in EU• Certified medical device*
Size: approx. 19.8 x 18.6 x 29.5 cm Weight: approx. 1,5 kg Power supply: 230 V/50 Hz Power output: 300 Watts Article No.: 88257 EAN: 40 15588 88257 9 Warranty: 3 Years	Size: approx. 16.5 x 13.5 x 20.5 cm Weight: approx. 0.36 kg Power supply: 230-240 V/50-60 Hz Power output: 100 Watts Article No.: 88258 EAN: 40 15588 88258 6 Warranty: 3 Years	Size: approx. 19 x 14 x 22 cm Weight: approx. 0.65 kg Power supply: 230 V/50 Hz Power output: 150 Watts Article No.: 88254 EAN: 40 15588 88254 8 Warranty: 3 Years

* Certification according to MDD (Medical Device Directive)

Well-being.
With us, you can experience
it in your own home.

Soothing ambient air – Made for Life.

Mia, 9
Enjoys an asthma free day

Sabine, 38
Enjoys the freshness.

Elias, 44
Feels fit for everyday life.

Linus, 7
Always sleeps with the light on.

The subject of health is increasingly gaining importance. Even in your own home. Whether it concerns ecological building materials, energy from renewable sources or a healthy living room climate - people want to live healthy, and are therefore design their homes according to health criteria more than ever before. With indoor air innovation products, MEDISANA can ensure that everyone can create their own personal well-being and healthy feeling at home.

This is how room air should be.
Not too dry. Not too damp.

Christian, 46 & Annette, 42
Enjoy their time at home.

Air Humidifier AH 660

- Ultrasound technology to improve the room air quality
- Infinitely adjustable nebulisation intensity
- Highly effective through a micro fine atomisation
- Low-noise and energy-saving



Good, healthy air has a humidity of 40 - 60%.

The ideal humidity in rooms is 40 – 60%.
The AH 660 improves the room air.

People inhale about 20,000 litres of air every day. It is good if this air is as clean and soothing as possible. Using the finest ultrasound technology, the air humidifier AH 660 can ensure a pleasant feel-good climate.

There's something in the air.

AROMA
DISPENSER
FOR FRAGRANCE OILS

Nathalie, 25
Just feels well.



The AH 662 will perform convincingly in any environment.

The fragrance makes the difference.

Stimulating, calming or activating scents. There are many fragrance oils with very specific scents. The air humidifier AH 662 conveys the moisture into the room unobtrusively through a micro fine atomisation.

Air Humidifier
AH 662

- Air humidifier with ultrasound technology to improve the room air quality
- High efficiency through micro-fine atomisation
- Aroma dispenser for fragrance oils



Tim, 29 & Susan, 28
Enjoy the freshness.

Frank, 37
Has a clear head again.

Aroma Diffuser AD 620

- A fine mist of fragrance that spreads throughout your entire room
- High efficiency through micro-fine atomisation
- Wellness light with a colour change in 5 colours

AROMA

Wellness light

3 Years Warranty



Inhaling with a sigh of relief.

The elegantly designed Aroma Diffuser AD 620 provides a relaxing fragrance experience.

The Aroma Diffuser AD 620 has a gentle, soothing and relaxing effect. The aroma of scented oils is carried into the room through a micro-atomisation. The built-in water tank holds 300 ml to provide many hours of fragrant air.



With the interplay of 5 wellness colours, the AD 620 will ensure a good mood in the morning

Good air. Good atmosphere.

You will hardly notice it, but you will feel it. Relaxation and tranquillity will emerge. Enjoy the time.

Wellness Light demonstrably has a strong calming effect, and can relieve winter depression. Combined with fine aroma oils, the relaxing effect is intensified even more. You'll feel better and more concentrated.

The proper light and air for better concentration.



Aroma Diffuser AD 630

- Using ultrasound technology, creates a fine mist of fragrance that spreads throughout your entire room
- Wellness light with a colour change in 6 colours



The proper wellness feeling.

Well-being at the push of a button.

How much good air would you like today? The Aroma Diffuser AD 640 gives you complete freedom.

Whether it's a relaxing half hour break, or up to 3 hours of feel-good atmosphere. Thanks to the intuitive touch buttons, you will have full control. Up to 15 ml are atomized per hour. That is the proper amount for a pleasant indoor climate.

Christian, 46
Loves good air and a good atmosphere.

Annette, 42
Loves her home.



Aroma Diffuser AD 640

- Aroma diffuser with touch control buttons
- LED time display
- Timer for programming operation
- Wellness light with colour change in 3 colours



What is more pleasant than good air for relaxation?



Light has a stimulating or soothing effect. The smooth change of three colours subconsciously creates a good feeling.

Visibly cleaner air.

Hugo, 3
Loves fresh air.

Finn, 6 months.
Extremely active.



Room Air Purifier Medisana Air

- Cleans with an antibacterial particulate filter (F9>PM2.5)
- Frees the air of particulates, dust, pollens, odours and bacteria
- Particularly effective and quiet with a low energy consumption
- Allergy friendly
- LED filter change indicator












Free the air of dust, pollen, dust, odours and bacteria.









Beautifully designed and elegant, but a power pack in performance. Medisana Air cleans the air with an antibacterial dust filter, and an additional activated carbon filtration. If you desire, a 60 m² room once in 24 hours, or a 20 m² room three times in 24 hours. The choice is yours.



Easy handling: The colours indicate whether the filter has to be changed.

Air Humidifier	Air Humidifier
AH 660	AH 662
	
	
	
<ul style="list-style-type: none">• Air humidifier with ultrasound technology to improve the breathing air quality• High efficiency through micro-fine atomisation• Steplessly adjustable misting levels	<ul style="list-style-type: none">• Air humidifier with ultrasound technology to improve the breathing air quality• High efficiency through micro-fine atomisation• Integrated fragrance dispenser for fragrant oils
<ul style="list-style-type: none">• Removable tank• Automatic shut-off when the tank is empty• Low-noise and energy-saving: 30 watts• Incl. water filter• Water tank 4,500 ml• For rooms up to approx. 30 m²	<ul style="list-style-type: none">• Automatic shut-off when the tank is empty• Low-noise and energy-saving• Water level indicator (full=blue, empty=red)• Water tank 900 ml• Incl. cleaning brush
	
Incl. water filter	
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Nebulizing output:</div><div>Operating time:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 24 x 17 x 36 cm</div><div>approx. 1.3 kg</div><div>220-240 V/50-60 Hz</div><div>30 Watts</div><div>max. 300 ml/h</div><div>approx. 15 Std.</div><div>60067</div><div>40 15588 60067 8</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Nebulizing output:</div><div>Operating time:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 15 x 15 x 22 cm</div><div>approx. 0.45 kg</div><div>100-240V~, 50/60Hz, 1.6A</div><div>30 Watts</div><div>max. 60 ml/h</div><div>approx. 15 Std.</div><div>60077</div><div>40 15588 60077 7</div><div>3 Years</div></div>

Air Humidifier
UHW



<ul style="list-style-type: none">• Air humidifier with ultrasound technology to improve the breathing air quality• High efficiency through micro-fine atomisation• Steplessly adjustable misting levels
<ul style="list-style-type: none">• Automatic shut-off when the tank is empty• Low-noise and energy-saving• Water tank 4,200 ml• For rooms up to approx. 30 m²
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Nebulizing output:</div><div>Operating time:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 26,5 x 33 x 20 cm</div><div>approx. 1.46 kg</div><div>220-240 V/50-60 Hz</div><div>30 Watts</div><div>max. 300 ml/h</div><div>variable, depending on setting</div><div>60065</div><div>40 15588 60065 4</div><div>3 Years</div></div>

Air Humidifier
MEDIBREEZE



<ul style="list-style-type: none">• Air humidifier with ultrasound technology to improve the breathing air quality• High efficiency through micro-fine atomisation• 2 intensitiy levels to choose from
<ul style="list-style-type: none">• Switchable lights• Automatic shut-off when the tank is empty• Low-noise and energy-saving• Water tank 680 ml
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Nebulizing output:</div><div>Operating time:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 16 x 11 x 20 cm</div><div>approx. 1 kg</div><div>100-240 V/50-60 Hz/0.5 A, 13.5V DC 1A</div><div>14 Watts</div><div>max. 80 ml/h</div><div>Stufe 1: approx. 9 Std. Stufe 2: approx. 6 Std.</div><div>60003</div><div>40 15588 60003 6</div><div>3 Years</div></div>

Air Humidifier with Timer Function
ULTRABREEZE



<ul style="list-style-type: none">• Air humidifier with ultrasound technology to improve the breathing air quality• High efficiency through micro-fine atomisation• Timer for programming the misting start and duration• Pre-heating to kill germs and bacteria
<ul style="list-style-type: none">• Switchable lights• 2 intensitiy levels to choose from• Automatic shut-off when the tank is empty• Low-noise and energy-saving• Integrated clock with alarm function• Water tank 1,300 ml
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Nebulizing output:</div><div>Operating time:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 19.6 x 13.5 x 25.6 cm</div><div>approx. 1.1 kg</div><div>100 - 240V/50 Hz/1.5 A, 24 V DC 1.5A</div><div>36 Watts</div><div>max. 65 ml/h</div><div>Stufe 1: approx. 20 Std. Stufe 2: approx. 30 Std.</div><div>60050</div><div>40 15588 60050 0</div><div>3 Years</div></div>

Hygrometer
HG 100


<ul style="list-style-type: none">• Displays humidity and room temperature for a healthy climate indoors• Saves the highest and lowest values• For hanging on the wall (magnetic mount) or standing upright
<ul style="list-style-type: none">• Temperature measurement range: -10°C ~ +50°C (+14° ~ +122 F°)• Humidity measurement range: 20% ~ 90%• Comfort zone indicator (smiley display)• Time display• Alarm function• Incl. battery (1.5 V AAA)
<div><div>Size:</div><div>Weight:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 8.2 x 7 x 2 cm</div><div>approx. 0.07 kg</div><div>1 x 1.5 V (AAA)</div><div>60079</div><div>40 15588 60079 1</div><div>3 Years</div></div>

Room Air Purifier
MEDISANA AIR®



- Cleans room air using an antibacterial particulate filter (F9>PM2.5) from 3M® and additional activated carbon filtering; allergy-friendly
- Frees the air of particulates, dust, pollens, odours and bacteria
- Particularly effective and quiet with low power consumption

- Simple and safe to use thanks to the LED filter change indicator
- Switches automatically between the day and night modes (even quieter)
- Operated using the included power supply unit and USB cable
- For rooms up to approx. 60 m²
- Incl. USB cable
- Digital motor with an extra long operating life – 5 year warranty*



Incl. 3M® filter

Size:	approx. 14.8 x 14.8 x 30.7 cm
Weight:	approx. 0.6 kg
Power supply:	DC 5V/300mA
Power output:	1,5 Watts
Purifying output:	60m²/24h
Cable length:	approx. 1 m
Article No.:	60300
EAN:	40 15588 60300 6
Warranty:	5 Years*

* 3 year manufacturer's warranty plus 2 year warranty extension when you register online



Aroma Diffusor
AD 610



- Aroma diffuser for fragrant oils
- Using ultrasound technology, creates a fine mist of fragrance that spreads throughout your entire room
- Switchable lights
- Operation via mains cable or USB cable

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Incl. measure cup
- Incl. USB cable

Size:	approx. 12 x 9 x 7 cm
Weight:	approx. 0.12 kg
Power supply:	100-240V/50-60Hz, DC 5V/800mA; 4W
Water tank:	80 ml
Nebulizing output:	max. 15 ml/h
Operating time:	approx. 5-8 Std.
Article No.:	60080
EAN:	40 15588 60080 7
Warranty:	3 Years

Aroma Diffusor
AD 620



- Aroma diffuser for fragrant oils
- Using ultrasound technology, creates a fine mist of fragrance that spreads throughout your entire room
- Wellness light with a colour change in 5 colours
- Timer for programming operation: 2 and 4 hours

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Incl. measure cup

Size:	approx. 13.2 x 13.2 x 14.5 cm
Weight:	approx. 0.25 kg
Power supply:	100-240V/50-60Hz, DC 24V/500mA; 12W
Wassertank:	300 ml
Nebulizing output:	max. 30 ml/h
Operating time:	approx. 9-27 Std.
Article No.:	60082
EAN:	40 15588 60082 1
Warranty:	3 Years

Aroma Diffusor
AD 630



- Aroma diffuser for fragrant oils
- Using ultrasound technology, creates a fine mist of fragrance that spreads throughout your entire room
- Wellness light with a colour change in 6 colours
- Timer for programming operation: 1/2/3 hours

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Incl. measure cup

Size:	approx. 15 x 15 x 14 cm
Weight:	approx. 0.26 kg
Power supply:	AC 100-240 V/DC 12V/1A, 9W
Wassertank:	80 ml
Nebulizing output:	max. 18 ml/h
Operating time:	approx. 3 Std.
Article No.:	60084
EAN:	40 15588 60084 5
Warranty:	3 Years

Aroma Diffusor
AD 640



- Aroma diffuser for fragrant oils
- Using ultrasound technology, creates a fine mist of fragrance that spreads throughout your entire room
- Wellness light with a colour change in 3 colours
- Timer for programming operation: 30/60/120/180 minutes

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Touch control buttons
- LED time display
- Incl. measure cup

Size:	approx. 11.6 x 11.2 x 16.2 cm
Weight:	approx. 0.43 kg
Power supply:	100-240V/50-60Hz, DC 24V/12W
Wassertank:	100 ml
Nebulizing output:	max. 15 ml/h
Operating time:	approx. 8 Std.
Article No.:	60086
EAN:	40 15588 60086 9
Warranty:	3 Years

Marie, 32
Likes groomed legs.

Anna, 31
Has beautiful feet.

Laura, 27
Makes a good
impression everywhere.

Timeless.
Beauty and style.
Smart care products – Made for Life.

The trend towards individualisation is increasing in society. The desire to do something good for the body, to keep it fit, to eat healthy food - these are expressions of a new self-consciousness. Of course, this includes the growing trend for more body care - for women and men alike. It's part of our lifestyle and an expression of our joy of life. MEDISANA provides care products that meet everyone's care requirements.

Nail care – not just for the summer.

Small details with big impact on the attractiveness.

Finger and toe nails need a lot of attention to look good and stay healthy. With the manicure and pedicure device MP 840, you can pamper your nails and protect them from discolouration and cracking. The LED light supports the precise and correct application.



NEW

Manicure and Pedicure Device MP 840

- Care and treatment of nails, cuticles and slight calluses
- 7 high-quality accessory parts for various applications
- Additional nail-polishing attachment with 2 smoothing and polishing rollers
- With an LED light for a precise application



Celine, 19
Pays attention to her appearance.

Well groomed and healthy nails - the best foundation for beautiful nail polish.

Gentle on the skin. Tough on cellulite.

Lara, 26
Feels good about her appearance.

Sophie, 33
Likes showing her legs.

NEW Cellulite Massager AC 855

- Effective suction and rolling massage for stimulating the tissue for tangibly tauter skin
- Infrared heating function stimulates the circulation
- Additional vibration enhances the massage effect



Noticeably
firmer
skin

A massage with the AC 855, a thoroughly stimulating experience.

The cellulite massage device does not only stimulate the skin unilaterally, but revives it through suction, vibration and the supportive infrared. This makes the massage much more effective. The blood circulation is stimulated and the connective tissue is firmed.



A good feeling. Firmer skin,
better circulation in the connective tissue.



Hannah, 27
Purity in perfection.

NEW Facial Cleansing Brush FB 885

- Facial cleansing brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products
- Rechargeable lithium-ion battery



Clean from the very start.

Our face is exposed to many stimuli and environmental challenges throughout the entire day. It is therefore beneficial if the facial skin is cleansed and refreshed regularly.



Sun, wind, rain, dust. Our facial skin is exposed to many forms of stress every day. With the facial cleaning brush FB 885, the face can be thoroughly cleaned within 60 seconds. The gentle vibration removes dead skin cells, stimulates the circulation and prepares the skin perfectly for care products.



Incl. 4 cleansing attachments
(Sensitive skin, regular skin, peeling, cosmetics sponge)



3 cleaning attachments for normal and sensitive skin

Mirror, mirror in the hand...

The light is often a problem during the application of make-up. Too bright, too dark, too one-sided. The CM 835 can perfectly illuminate the face.

Have fun applying make-up alone or with friends. The results are convincing. The cosmetic mirror CM 835 perfectly illuminates the face, thanks to the ring illumination comprised of 12 white LED - and on two mirror surfaces. The face is reflected normal, or magnified by a factor of five.



Secure pedestal.
Now the make-up fun can begin.

Eva, 52
Always an eye-catcher.

NEW Cosmetic Mirror CM 835

- Illuminated cosmetic mirror
- 2 swivelling mirrors:
regular and 5x magnification
- Ring illumination from 12 white LEDs



Hair Removal System
SILHOUETTE IPL 800



- Effective technology for permanent hair reduction
- Suitable for the face, arms, legs, underarms, bikini zone, back, chest
- Up to 50% hair reduction already after 4 weeks
- Incl. 2 cartridges with 50,000 light pulses each

- Disables the hair root deep under the skin and prevents further hair growth
- Easy and safe to use at home
- 5 light intensity levels
- Safe and effective for the treatment of skin types I - IV and black to to dark blond hair
- Integrated UV filter
- Flash-glide function for comfortable use
- Sensor starts treatment when skin contact is made
- Automatic shut-off after 10 min.
- Mains adapter included



Replacement cartridge with 50.000 light pulses incl.

Size:	approx. 14 x 8.1 x 4.8 cm
Weight:	approx. 0.2 kg
Lighting area:	2.7 cm ²
Intensity:	4J /cm ²
Wavelength:	475-1,200 nm
Power supply:	100-240V, 50/60Hz 0.8A-DC 12V, 3A
Article No.:	88580
EAN:	40 15588 88580 8
Warranty:	3 Years

INFO Hair Removal System
Silhouette IPL 800

Gentle and long term. This is how modern hair removal functions with the IPL 800 for use at home.

Many people simply cannot stand the conventional methods for hair removal, such as shaving, epilation or waxing. On the other hand, hair removal with light pulses is also used by doctors. The light pulses deactivate the root of the hair deep under the skin and thus prevent the hair from growing back.

Facial Cleansing Brush
FB 880



- Facial cleansing brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products

- Cleanses the skin in 60 seconds (Timer)
- 2 speed settings: normal or intensive cleansing
- Switches off automatically after 1 minute
- Incl. 3 cleansing attachments (Sensitive skin, regular skin, peeling)
- batteries last up to 70 min.
- Splash-proof
- Incl. protective cover, storage/carrying bag and batteries (1.5 V AA)



Incl. 3 cleansing attachments (Sensitive skin, regular skin, peeling)

Size:	approx. 15.7 x 5.7 x 7.9 cm
Weight:	approx. 0.22 kg
Batteries:	2 x 1.5 V (AA)
Article No.:	88560
EAN:	40 15588 88560 0
Warranty:	3 Years

Facial Cleansing Brush
FB 885



NEW



- Facial cleansing brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products
- Rechargeable lithium-ion battery

- Cleanses the skin in 60 seconds (Timer)
- 4 speed settings
- Switches off automatically after 1 minute
- Incl. 4 cleansing attachments (Sensitive skin, regular skin, peeling, cosmetics sponge)
- Splash-proof
- Incl. protective cover, storage/carrying bag and charge base



Incl. 4 cleansing attachments (Sensitive skin, regular skin, peeling, cosmetics sponge)

Size:	approx. 18.6 x 5.0 x 2.7 cm
Weight:	approx. 0.27 kg
Batteries:	Lithium-ion battery
Power output:	500 mAh
Cable length:	1.30 m
Article No.:	88565
EAN:	40 15588 88565 5
Warranty:	3 Years

Facial Sauna
FSS



Q4/2016

- Opens the pores, cleans the skin and relaxes facial lines with humid heat
- With aroma attachment for beneficial aroma treatment
- Easy cleaning and handling

- For cosmetic use
- Diffuser for even steam distribution
- Additional nose piece
- High-quality heating element
- Indicator light
- Measuring cup for easy filling



Incl. nose piece and measuring cup

Size:	approx. 30 x 20 x 19.5 cm
Weight:	approx. 0.56 kg
Power supply:	100-240 V/50-60 Hz
Power output:	max. 100 Watts
Article No.:	88245
EAN:	40 15588 88245 6
Warranty:	3 Years

Manicure and Pedicure Device
MP 840



Q4/2016

- Care and treatment of nails, cuticles and slight calluses
- With LED light for accurate application
- 3 speed levels for the individual requirements of care-giving


- Clockwise and counter-clockwise rotation
- High quality and durable attachments
- Stops automatically when the pressure is too high – also suitable for diabetics
- Additional nail polishing attachment with 2 smoothing and polishing rollers
- Natural, durable nail gloss (without additional care products)
- Simple to use, easy replacement of the attachments
- With storage case



Incl. 7 attachments: Sapphire cone (coarse), Sapphire disc (fine), Cylindrical burr, Flame burr, Felt cone, Needle burr as well as an additional attachment for polishing nails

Size:	approx. 3.8 x 16.1 cm
Weight:	approx. 0.27 kg
Power supply:	230 V, 50 Hz, DC18V/500mA
Rotation speed:	3 Stufen: 5,000/6,000/7,000
Power output:	9 Watts
Cable length:	3 m
Article No.:	85155
EAN:	40 15588 85155 1
Warranty:	3 Years

Manicure and Pedicure Device
MEDISTYLE L



- Care and treatment of nails, cuticles and slight calluses
- With LED light for accurate application
- Continuous speed adjustment

- Clockwise and counter-clockwise rotation
- High quality and durable sapphire attachments
- Stops automatically when the pressure is too high – also suitable for diabetics
- With storage case



Incl. 10 attachments:
Manicure: Flame burr, Sapphire disc (fine), Sapphire disc (coarse), Felt cone, Cuticle pusher. **Pedicure:** Sapphire cone (coarse), Sapphire cone (long), Sapphire grinder (cylindrical), Cylindrical burr, Needle burr

Size:	approx. 15 x 4 x 4 cm
Weight:	approx. 0.13 kg
Power supply:	100-240 V, 50/60 Hz, 18V DC 500 mA
Rotation speed:	1,600-5,800 U/min
Article No.:	85130
EAN:	40 15588 85130 8
Warranty:	3 Years

Manicure and Pedicure Device
MPS



- Care and treatment of nails, cuticles and slight calluses
- Continuous speed adjustment

- Clockwise and counter-clockwise rotation
- High quality and durable sapphire attachments
- Stops automatically when the pressure is too high – also suitable for diabetics
- With storage case



Incl. 4 attachments: Sapphire cone (coarse), Sapphire disc (fine), Flame burr, Felt cone

Size:	approx. 15 x 4 x 4 cm
Weight:	approx. 0.13 kg
Power supply:	100-240 V, 50/60 Hz, 18V DC 500 mA
Rotation speed:	3,800-5,500 U/min
Article No.:	85126
EAN:	40 15588 85126 1
Warranty:	3 Years

Electric Nail Polisher
NP 860



Q1/2017

- For glossy, natural-looking nails
- Professional nail care at home
- Durable nail gloss (without additional care products)

- Simple to use on fingernails and toenails
- Small and handy – ideal even when travelling
- 2 speed levels for individual needs of nail care
- Ergonomic design
- Smoothing and polishing – in two steps to nails with the “Wow” effect
- Incl. 2 batteries (1.5 V AA)



Incl. 2 smoothing and polishing rollers

Size:	approx. 15.5 x 3.6 x 3.3 cm
Weight:	approx. 0.15 kg (incl. batteries)
Batteries:	2 x 1.5 V (AA)
Rotation speed:	2 Stufen: 2,000/2,500 U/min
Article No.:	85160
EAN:	40 15588 85160 5
Warranty:	3 Years




Electric Callus Remover
CR 860







- Gently and effectively remove excess calluses
- For smoother hands, feet and elbows
- Apply to dry skin

- 2 speed levels
- Ergonomic, non-slip handle
- Incl. cleaning brush, protective cover, storage bag and batteries



Incl. 2 attachments (fine and coarse)

Size:	approx. 15.7 x 6.7 x 4.6 cm
Weight:	approx. 0.175 kg (incl. batteries)
Batteries:	2 x 1.5 V (AA)
Rotation speed:	2 Stufen: 1,600/1,800 U/min
Article No.:	88570
EAN:	40 15588 88570 9
Warranty:	3 Years

Cellulite Massager AC 850	Cellulite Massager AC 855
	 <div>NEW</div>
	
	Q4/2016
<ul style="list-style-type: none">• Effective massage to stimulate the connective tissue for noticeably firmer skin• Can help improve blood circulation• 2 selectable massage intensities• Light and handy design for convenient self-massage	<ul style="list-style-type: none">• Effective suction and rolling massage for stimulating the tissue for tangibly tauter skin• Infrared heating function stimulates the circulation• Additional vibration enhances the massage effect• Promotes lymphatic drainage and strengthens the tissue
<ul style="list-style-type: none">• 3 × 2 rotating massage rollers• Individually adjustable hand strap• Single-hand operation• Particularly powerful thanks to mains operation• Easy cleaning of the removable massage rollers	<ul style="list-style-type: none">• Supports the treatment of cellulite• The massaging suction intensity is continuously variable• Light and handy design for a comfortable self-massage• Single-hand operation• Particularly powerful thanks to mains operation• Removable handle
<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Cable length:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 13 x 11 x 8.5 cm</div><div>approx. 0.35 kg</div><div>110-240V/50-60Hz; 12V DC 1200 mA</div><div>14.4 Watts</div><div>2 m</div><div>88540</div><div>40 15588 88540 2</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Power supply:</div><div>Power output:</div><div>Cable length:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 12.6 x 10.6 x 8.2 cm</div><div>approx. 0.26 kg</div><div>110-240V/50-60Hz; 9V DC 600 mA</div><div>5.4 Watts</div><div>1.80 m</div><div>88545</div><div>40 15588 88545 7</div><div>3 Years</div></div>



Cosmetic Mirror CM 835	2in1 Cosmetic Mirror CM 840	3in1 Cosmetic Mirror CM 845
 <div>NEW</div>		
		
	Q4/2016	
<ul style="list-style-type: none">• Illuminated cosmetic mirror• 2 swivelling mirrors: Normal and 5x magnification• Ring illumination from 12 white LEDs• Diameter of mirror approx. 12 cm	<ul style="list-style-type: none">• Illuminated cosmetics mirror• 2 swivelling mirrors: Normal and 5x magnification• Ring illumination from 17 white LEDs• Diameter of mirror approx. 13 cm	<ul style="list-style-type: none">• Illuminated cosmetics mirror incl. wall bracket and stand• 2 swivelling mirrors: Normal and 5x magnification• Ring illumination fromwith 17 white LEDs• Diameter of mirror approx. 13 cm
<ul style="list-style-type: none">• High-quality chrome finish• Stable base• Incl. 3 batteries (1.5 V AAA)	<ul style="list-style-type: none">• High-quality chrome finish• Stable base• Incl. 4 batteries (1.5 V AAA)	<ul style="list-style-type: none">• High-quality chrome finish• Includes wall bracket, stand, batteries and screws for mounting onto the wall• Incl. 4 batteries (1.5 V AAA)
		
		Includes wall bracket and stand
<div><div>Size:</div><div>Weight:</div><div>Spiegelfläche:</div><div>Batteries:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 16.9 x 8 x 19 cm</div><div>approx. 0.375 kg</div><div>approx. ø 12 cm</div><div>3 x 1.5 V (AAA)</div><div>88554</div><div>40 15588 88554 9</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Spiegelfläche:</div><div>Batteries:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 20.5 x 13.5 x 31.5 cm</div><div>approx. 0.95 kg</div><div>approx. ø 13 cm</div><div>4 x 1.5 V (AAA)</div><div>88550</div><div>40 15588 88550 1</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Spiegelfläche:</div><div>Batteries:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 17.5 x 24/47 x 36.5 cm</div><div>approx. 1.55 kg</div><div>approx. ø 13 cm</div><div>4 x 1.5 V (AAA)</div><div>88552</div><div>40 15588 88552 5</div><div>3 Years</div></div>

Esther, 31
Likes trying out new ingredients.

Lukas, 34
Very accurate when cooking.

Lily, 3
Not very accurate.

Always. The best ingredients:
Accuracy, enjoyment, health.

Professional kitchen helpers – Made for Life.

In the future, everything will continue to revolve around enjoying food in the kitchen. Food will be more valuable to us. It should be healthy and nutritious for us. We aim to enjoy it together or alone at home. Cooking will become an adventure, dining together will become a highlight of the day. People rely on good recipes and good kitchen appliances. MEDISANA supports the desire for healthy good food with kitchen utensils, such as scales, which greatly ease the kitchen experience.

NEW Digital Kitchen Scale with App KS 250

- Determination and display of calories, protein, fat, cholesterol, carbohydrates, fibres and sodium in an app
- Sensor touch keypads
- Large LCD display
- High quality stainless steel weighing surface



**SMARTER
WEIGHING**
WITH
EXTENSIVE NUTRITIONAL
VALUE APP




How balanced the diet is, will be output the nutritional value App.



Strictly speaking, for the best recipe.

Enjoy your meal – thanks to good nutritional values. The kitchen scale with App will specify what and how much.

The high-quality kitchen scale is not just a visual highlight. When connected to the accompanying App you can receive reliable information on the key nutritional values, such as calories, protein, fat, cholesterol, carbohydrates, fibres and sodium. How? Select a food, and then weigh it. At the same time, it will maintain a nutritional diary of all past values.



Lisa, 6
Loves cake.

Emma, 5
Already a skilled cook.



Healthy food, good mood. The KS 250 provides the certainty that the diet will be balanced.

5 kg capacity, large LCD display, anti-slip feet, stainless steel weighing surface. What will be conjured up on the table today?



Digital Kitchen Scale

KS 200



- Energy-saving automatic shut-off function
- Overload indicator
- High-quality steel weighting surface
- Anti-slip feet

- Practical add-weight function (tare)
- On/off switch
- Low battery display
- Bears weights up to 3 kg
- Incl. battery (CR2032 3V)

Size: approx. 16 x 11,7 x 2,7 cm
Weight: approx. 0,24 kg
Weighting surface: approx. 11 x 11 cm
LCD Display: approx. 1,6 x 4,6 cm
Weight display: g, oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40469
EAN: 40 15588 40469 6
Warranty: 3 Years



Digital Glass Kitchen Scale

KS 210



- Touch sensor keys
- Energy-saving automatic shut-off function
- Flat construction
- Glass weighting surface – easy to clean

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- Bears weights up to 5 kg / 5000 ml
- Incl. battery (CR2032 3V)

Size: approx. 22.4 x 18.4 x 1.7 cm
Weight: approx. 0.6 kg
Weighting surface: approx. 22.4 x 18.4 cm
LCD Display: approx. 5.5 x 2.5 cm
Weight display: g, kg, ml, fl oz und lb oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40463
EAN: 40 15588 40463 4
Warranty: 3 Years

Digital Glass Kitchen Scale

KS 210



- Touch sensor keys
- Energy-saving automatic shut-off function
- Flat construction
- Glass weighting surface – easy to clean

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- Bears weights up to 5 kg / 5000 ml
- Incl. battery (CR2032 3V)

Size: approx. 22.4 x 18.4 x 1.7 cm
Weight: approx. 0.6 kg
Weighting surface: approx. 22.4 x 18.4 cm
LCD Display: approx. 5.5 x 2.5 cm
Weight display: g, kg, ml, fl oz und lb oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40471
EAN: 40 15588 40471 9
Warranty: 3 Years

Digital Glass Kitchen Scale

KS 210



- Touch sensor keys
- Energy-saving automatic shut-off function
- Flat construction
- Glass weighting surface – easy to clean

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- Bears weights up to 5 kg / 5000 ml
- Incl. battery (CR2032 3V)

Size: approx. 22.4 x 18.4 x 1.7 cm
Weight: approx. 0.6 kg
Weighting surface: approx. 22.4 x 18.4 cm
LCD Display: approx. 5.5 x 2.5 cm
Weight display: g, kg, ml, fl oz und lb oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40472
EAN: 40 15588 40472 6
Warranty: 3 Years

Digital Glass Kitchen Scale

KS 210



- Touch sensor keys
- Energy-saving automatic shut-off function
- Flat construction
- Glass weighting surface – easy to clean

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- Bears weights up to 5 kg / 5000 ml
- Incl. battery (CR2032 3V)



Size: approx. 22.4 x 18.4 x 1.7 cm
Weight: approx. 0.6 kg
Weighting surface: approx. 22.4 x 18.4 cm
LCD Display: approx. 5.5 x 2.5 cm
Weight display: g, kg, ml, fl oz und lb oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40473
EAN: 40 15588 40473 3
Warranty: 3 Years

Digital Glass Kitchen Scale	Digital Kitchen Scale	Digital Glass Kitchen Scale XL
KS 220	KS 230	KS 240
		
		
<ul style="list-style-type: none">• With time and temperature• Energy-saving automatic shut-off function• Timer Function• Glass weighting surface – easy to clean	<ul style="list-style-type: none">• Touch sensor keys• Energy-saving automatic shut-off function• Large, easy-to-read LCD display with white luminous digits• High-quality stainless steel weighting surface• Hanging fixture• Flat construction	<ul style="list-style-type: none">• Touch sensor keys• HOLD / Memory function• Energy-saving automatic shut-off function• Glass weighting surface – easy to clean• Hanging fixture
<ul style="list-style-type: none">• Practical add-weight function (tare)• Overload indicator• Anti-slip feet• On/off switch• Low battery display• Bears weights up to 5 kg• Incl. battery (CR2032 3V)	<ul style="list-style-type: none">• Practical add-weight function (tare)• Overload indicator• Anti-slip feet• On/off switch• Low battery display• Bears weights up to 5 kg / 5,000 ml• Flat construction• Incl. 2 batteries (CR2032 3V)	<ul style="list-style-type: none">• Practical add-weight function (tare)• Overload indicator• Anti-slip feet• On/off switch• Low battery display• Bears weights up to 20 kg / 20,000 ml• Flat construction• Incl. battery (CR2032 3V)
<div><div>Size:</div><div>Weight:</div><div>Weighting surface:</div><div>LCD Display:</div><div>Weight display:</div><div>Graduation:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 15 x 22.7 x 3.6 cm</div><div>approx. 0.54 kg</div><div>approx. ø 15 cm</div><div>approx. 2.4 x 5.3 cm</div><div>g, kg, oz und lb oz</div><div>1 g</div><div>1 x CR2032 3V</div><div>40467</div><div>40 15588 40467 2</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Weighting surface:</div><div>LCD Display:</div><div>Weight display:</div><div>Graduation:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 20.2 x 21.8 x 1.6 cm</div><div>approx. 0.34 kg</div><div>approx. 20.2 x 18.0 cm</div><div>approx. 2.3 x 5.3 cm</div><div>g, ml, oz, fl oz und lb oz</div><div>1 g</div><div>2 x CR2032 3V</div><div>40465</div><div>40 15588 40465 8</div><div>3 Years</div></div>	<div><div>Size:</div><div>Weight:</div><div>Weighting surface:</div><div>LCD Display:</div><div>Weight display:</div><div>Graduation:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 21 x 29.7 x 1.7 cm</div><div>approx. 0.82 kg</div><div>approx. 21 x 29.7 cm</div><div>approx. 6.5 x 3.5 cm</div><div>g, kg, ml, fl oz und lb oz</div><div>1 g</div><div>1 x CR2032 3V</div><div>40468</div><div>40 15588 40468 9</div><div>3 Years</div></div>

INFO Digital Glass Kitchen Scale XL
KS 240

Every kitchen should have them. Nicely shaped and able to take up to 20 kg.

Cooking and baking become even better. The KS 240 sensor-touch touch buttons quickly makes itself indispensable in the kitchen. Up to 20 kg can be weighed precisely. The overload display reacts at once to an overload. The glass surface is easy to clean. The weight displays at the top and on the side also help when weighing letters or parcels.

Digital Kitchen Scale with App
KS 250


<div>Q4/2016</div>
<ul style="list-style-type: none">• Determination and display of calories, protein, fat, cholesterol, carbohydrates, fibre and sodium in an app• Sensor touch keypads• Large, easy-to-read LCD display with white luminous numerals• High-quality stainless steel weighting surface
<ul style="list-style-type: none">• Practical add-weight function (tare)• Energy-saving automatic shut-off function• Overload indicator• Anti-slip feet• On/off switch• Low battery display• Bears weights up to 5 kg• Flat construction• Incl. 2 batteries (CR2032 3V)
<div><div>App functions</div><div><ul style="list-style-type: none">• Foods are simply selected in the app and then weighed• Nutrition diary: keeping a record of all the measured data by date and time• Multifunctional nutritional value determination• More than 500 foods saved in the app• Creating your own foods possible</div></div>
<div><div>Size:</div><div>Weight:</div><div>Weighting surface:</div><div>LCD Display:</div><div>Weight display:</div><div>Graduation:</div><div>Battery:</div><div>Article No.:</div><div>EAN:</div><div>Warranty:</div></div> <div><div>approx. 20.2 x 18.3 x 1.7 cm</div><div>approx. 0.34 kg</div><div>approx. 20 x 18 cm</div><div>approx. 2.3 x 5.3 cm</div><div>g, oz, lb oz, ml, fl oz</div><div>1 g</div><div>2 x CR2032 3V</div><div>40474</div><div>40 15588 40474 0</div><div>3 Years</div></div>

Upper arm cuffs
for blood pressure monitors



MTM MTP	
Large cuff (32-42 cm)	
Article No.:	51037
EAN:	40 15588 51037 3
Normal cuff (22-32 cm)	
Article No.:	51038
EAN:	40 15588 51038 0
Small cuff (17-22 cm)	
Article No.:	51034
EAN:	40 15588 51034 2

MTV MTC MTD MTS MTR CARDIOCOMPACT	
Large cuff (30-42 cm)	
Article No.:	51126
EAN:	40 15588 51126 4
Normal cuff (22-30 cm)	
Article No.:	51135
EAN:	40 15588 51135 6

MTX	
Große “Easy form” cuff (22-42 cm)	
Article No.:	51087
EAN:	40 15588 51087 8
Normale “Easy form” cuff (22-32 cm)	
Article No.:	51088
EAN:	40 15588 51088 5

BU 550 CONNECT BU 575 CONNECT	
Normal cuff (22-42 cm)	
Article No.:	51299
EAN:	40 15588 51299 5

BU 510 BU 530 CONNECT	
Large cuff (32-42 cm)	
Article No.:	51169
EAN:	40 15588 51169 1
Normal cuff (22-32 cm)	
Article No.:	51168
EAN:	40 15588 51168 4

CARDIODOCK CARDIODOCK 2	
Große “Easy form” cuff (30-42 cm)	
Article No.:	51284
EAN:	40 15588 51284 1
Normale “Easy form” cuff (22-32 cm)	
Article No.:	51283
EAN:	40 15588 51283 4

Mains adapters
for blood pressure monitors



MTM MTP MTX	
Article No.:	51036
EAN:	40 15588 51036 6

MTV MTC MTD MTS BU 510 BU 530	
Article No.:	51125
EAN:	40 15588 51125 7

CARDIOCOMPACT	
Article No.:	51095
EAN:	40 15588 51095 3

BU 550	
Article No.:	51013
EAN:	40 15588 51013 7

BU 575 CONNECT BW 300 CONNECT	
Article No.:	51014
EAN:	40 15588 51014 4

Protective sleeves
for thermometer



FTF	
50 pcs.	
Article No.:	76111
EAN:	40 15588 76111 9

Accessories
for blood glucose monitors



MEDITOUCH GLUCODOCK	
Control solution	
Article No.:	79029
EAN:	40 15588 79029 4

MEDITOUCH 2 MEDITOUCH 2 CONNECT	
Control solution	
Article No.:	79039
EAN:	40 15588 79039 3
Test stripes	
Article No.:	79038
EAN:	40 15588 73038 6

MEDITOUCH MEDITOUCH 2 MEDITOUCH 2 CONNECT GLUCODOCK	
Lancing device	
Article No.:	79002
EAN:	40 15588 79002 7
Lancets	
Article No.:	79028
EAN:	40 15588 79028 7

MEDITOUCH 2	
USB cable	
Article No.:	79049
EAN:	40 15588 79049 2

Wristbands
for activity tracker



VIFIT CONNECT	
Black	
Article No.:	79417
EAN:	40 15588 79417 9
Orange	
Article No.:	79418
EAN:	40 15588 79418 6
Blue	
Article No.:	79470
EAN:	40 15588 79470 4
Green	
Article No.:	79471
EAN:	40 15588 79471 1
Red	
Article No.:	79472
EAN:	40 15588 79472 8
Pink	
Article No.:	79473
EAN:	40 15588 79473 5



Accessories
for bath spas



BBS MBH	
Air hose	
Article No.:	30666
EAN:	40 15588 30666 2
BBS	
remote control	
Article No.:	88384
EAN:	40 15588 88384 2
Soft mat	
Article No.:	88383
EAN:	40 15588 88383 5

Loofah massage pad
for mini handmassager



HM 840	
Article No.:	88526
EAN:	40 15588 88526 6

Electrodes
for pain therapy devices



TDD TDP TED 4 Multiple electrodes (approx. 3 x 5 cm) Article No.: 87080 EAN: 40 15588 87080 4 2 Multiple electrodes (approx. 5 x 9 cm) Article No.: 87100 EAN: 40 15588 87100 9
TDD TDP Cable set Article No.: 88315 EAN: 40 15588 88315 6
BT 850 Adhesive pads (2 pcs.) Article No.: 88334 EAN: 40 15588 88334 7

Replacement lamp
for infrared lamps



IRH 100 Watts Article No.: 30101 EAN: 40 15588 30101 8
IRL 150 Watts Article No.: 88256 EAN: 40 15588 88256 2

Accessories
for Inhalators



IN 500 IN 550 Air filter and nebulizer kit with 3 nebulizer attachments Article No.: 54524 EAN: 40 15588 54524 5 Air filter and nebulizer kit Article No.: 54525 EAN: 40 15588 54525 2 Mouthpiece and nose adapter Article No.: 54526 EAN: 40 15588 54526 9 Air hose Article No.: 54527 EAN: 40 15588 54527 6 Mask for children Article No.: 54528 EAN: 40 15588 54528 3 Mask for adults Article No.: 54529 EAN: 40 15588 54529 0 Nebulizer attachments (3 pcs.) Article No.: 54531 EAN: 40 15588 54531 3

USC Mask set (2 masks + adapter) Article No.: 54101 EAN: 40 15588 54101 8 Mouthpiece Article No.: 54102 EAN: 40 15588 54102 5 Covering cap Article No.: 54103 EAN: 40 15588 54103 2 Nebulizer kit Article No.: 54104 EAN: 40 15588 54104 9

Aroma essences
for Aroma diffusers
and air humidifiers



Fragrance Lemon 10 ml Article No.: 60030 EAN: 40 15588 60030 2 Fragrance Eucalyptus 10 ml Article No.: 60031 EAN: 40 15588 60031 9 Fragrance Lavender 10 ml Article No.: 60032 EAN: 40 15588 60032 6 Fragrance Citronella 10 ml Article No.: 60033 EAN: 40 15588 60033 3 Fragrance Ylang-Ylang 10 ml Article No.: 60035 EAN: 40 15588 60035 7 Fragrance Peppermint 10 ml Article No.: 60036 EAN: 40 15588 60036 4 Fragrance Orange 10 ml Article No.: 60037 EAN: 40 15588 60037 1 Fragrance Pinewood 10 ml Article No.: 60038 EAN: 40 15588 60038 8 Fragrance Wildberries 10 ml Article No.: 60039 EAN: 40 15588 60039 5 Fragrance Oceanbreeze 10 ml Article No.: 60040 EAN: 40 15588 60040 1
--

Replacement filter
for air purifier



MEDISANA AIR® 1 pc. Article No.: 60390 EAN: 40 15588 60390 7 2 pcs. Article No.: 60391 EAN: 40 15588 60391 4

Replacement water filter
for air humidifier



AH 660 2 pcs. Article No.: 60069 EAN: 40 15588 60069 2
--

Replacement cartridge
für hair removal system



IPL 800 Replacement cartridge (50.000 light pulses) Article No.: 88581 EAN: 40 15588 88581 5
--

Replacement brushes
for facial cleansing brushes



FB 880 3 pcs. Article No.: 88561 EAN: 40 15588 88561 7 FB 885 4 pcs. Article No.: 88566 EAN: 40 15588 88566 2
--

Replacement rolls
for electric callus remover



CR 860 2 pcs. Article No.: 88572 EAN: 40 15588 88572 3
--

Attachments
for manicure and pedicure devices



MANILUX L MANILUX S MP 820 Set of attachments (10 pieces)* Article No.: 85450 EAN: 40 15588 85450 7 *contains: Coarse sapphire cone, coarse sapphire disc, flame burr, sapphire ball burr, felt cone, long sapphire cone, fine sapphire disc, needle burr, cylindrical burr, drainage grinder
--

MEDISTYLE L MEDISTYLE S MPS Felt cone Article No.: 85001 EAN: 40 15588 85001 1 Coarse sapphire cone Article No.: 85003 EAN: 40 15588 85003 5 Fine sapphire disc Article No.: 85004 EAN: 40 15588 85004 2 Sapphire ball burr Article No.: 85005 EAN: 40 15588 85005 9 Long needle burr Article No.: 85006 EAN: 40 15588 85006 6 Long sapphire cone Article No.: 85007 EAN: 40 15588 85007 3 Cuticle pusher Article No.: 85008 EAN: 40 15588 85008 0 Cylindrical burr Article No.: 85009 EAN: 40 15588 85009 7 Coarse sapphire disc Article No.: 85010 EAN: 40 15588 85010 3 Rotary flame burr Article No.: 85011 EAN: 40 15588 85011 0 Mandrel Article No.: 85015 EAN: 40 15588 85015 8
--

MEDISANA®

Jagenbergstraße 19
41468 Neuss

Phone: +49 2131/36 68-0
Fax: +49 2131/36 68-50 95

info@medisana.de
www.medisana.com

September 2016

Mobile Health Health Control Sport Wellness Therapy
Healthy Home Personal Care Kitchen Scales