

MEDISANA®



Made for life

PRODUCT CATALOGUE 2018

MOBILE HEALTH HEALTH CONTROL
SPORT WELLNESS THERAPY HEALTHY HOME
BODY CARE KITCHEN SCALES

MEDISANA®

Made for Life



Excellent products, excellent quality.

Health is the most important personal asset for most people. From this, they derive quality of life and satisfaction. That is why we have been consistently setting new standards in the fields of health, wellness and wellbeing for more than 30 years. We support our customers with our products in order to keep them healthy, fit and well-groomed.

With the highest degree of product quality, ease of use and safety. Thereby, it is important for us to always meet the criteria of the Medical Devices Directive (MDD), but also the criteria of various consumer organisations. However, the satisfaction of our customers with our products is always paramount for us.

Contents

Mobile health _____	6	Heat cape HP 626 _____	62
MEDISANA App-World _____	8	Cosy heating blanket XXL HB 675 _____	63
VitaDock+ App 2.0 _____	9	Product overview _____	64
Health control _____	10	Therapy _____	84
Upper arm blood pressure monitor BU 535 Voice _____	12	Daylight lamp LT 460 _____	86
Upper arm blood pressure monitor MTP Pro _____	14	Wake up light WL 444 _____	88
Wrist blood pressure monitor BW 335 _____	15	Abdominal muscle stimulator AM 880 _____	89
Glass personal scale PS 425 _____	16	Inhaler IN 500 compact _____	90
Motif personal scale PS 405 _____	18	Electric lice comb LC 870 _____	91
Design personal scale PS 435 _____	19	Product overview _____	92
Multifunctional thermometer TM 750 connect _____	20	Healthy home _____	98
Infrared thermometer TM 770 _____	21	Aroma diffusor AD 625 _____	100
Blood glucose meter MediTouch 2 connect _____	22	Air humidifier AH 661 _____	102
Product overview _____	24	Air humidifier AH 662 _____	103
Sport _____	38	Product overview _____	104
Foam ball PowerRoll-Ball _____	40	Body care _____	108
Foam roller PowerRoll / PowerRoll-XT _____	42	Hair removal system IPL 850 _____	110
Foam roller PowerRoll-Slim _____	44	Hair removal system IPL 840 _____	111
Foam roller VarioRoll _____	46	Manicure and pedicure device MP 815 _____	112
EMS Body Trainer _____	47	UV nail hardener ND 870 _____	113
Product overview _____	48	Cellulite massager AC 850 _____	114
Wellness _____	52	Cellulite massager AC 855 _____	115
Shiatsu massage seat cover MC 826 _____	54	Face cleansing brush FB 885 _____	116
Shiatsu massage seat cover MC 824 _____	55	Facial sauna FSS _____	117
Comfort massage seat cover MC 818 _____	56	Product overview _____	118
Shiatsu massage seat cover MC 815 _____	57	Kitchen scales _____	124
Shiatsu acupressure massage seat cover MC 825 _____	58	Digital kitchen scale with App KS 250 _____	126
Shiatsu massage seat cover MCN _____	59	Product overview _____	128
Shiatsu foot massager FM 885 _____	60	Accessories _____	132
Foot massager FM 888 _____	61		

Stefanie, 32
Always has her family doctor at her side.

Lara, 34
Uses the mobile diagnostics.

Patrick, 31
Checks himself.

Our health at a glance.

Mobile health management – Made for Life.

More and more people are leaving nothing to chance when it comes to health, fitness, prevention and wellness. They want to actively participate, control their vital signs and take care of their body. Everything begins with sufficient knowledge about one's own vital signs. MEDISANA keeps all of them up-to-date them with integrated solutions. Mobile and online. For the best possible personal health management.

The vital signs your own body.

There are many factors in the body which can be influenced in a targeted manner. The MEDISANA App World supports the self-conducted check and self-optimisation

Conscious living means to deliberately focus on the values and functions of one's own body. What are the values for blood pressure, blood oxygen and blood glucose? How is my activity developing? Answers are provided by the apps from the MEDISANA App-world. Reliable, comprehensive and easy to use.



 **VitaDock+ App**

Recording of all vital signs such as blood pressure, blood glucose, body temperature, weight and activity.



 **Bodytoner App**

Shape and intensity of the massage can be individually controlled.



 **VitaTherm App**

Continuous monitoring of body temperature.



 **Sleepace App**

How well do you sleep? The app for monitoring the sleep behavior.



 **Thermo App**

Continuous real-time monitoring of body temperature.



 **Nährwerte App**

The kitchen scale app informs about the most important nutritional values.

To have an optimal overview of the self.

The new VitaDock+ App has many new features, and makes it even easier for all those who want to keep an overview of their vital signs.

The VitaDock+ 2.0 is able to impress with a completely redesigned user interface. While the vital signs could already be supplemented with brief commentaries, the dashboard now also allows short commentaries on the last vital value in the weekly, monthly and yearly graphs – and this can be done with evaluations according to the recommendations of the WHO*. More results on the screen allow for better insights.





VitaDock+ App 2.0

- Simplified connection of Bluetooth devices
- Registration via Facebook
- List of measurement results
- Comment functions



Max, 31
Always has an overview of his values.

Download the new VitaDock+ App 2.0 here:



Karl, 72
Has his diabetes
under control.


Isabelle, 30
Always has her
weight in view.

Andreas, 35
Regularly checks
his blood pressure.

More safety for your own lifestyle.

Innovative health management – Made for Life.

Health is a key factor of our lifestyle. A lifestyle that we would like to maintain. Many people are taking responsibility for their own health. They pay attention to themselves, listen to the needs of their own body, and actively control their own health. MEDISANA supports them with new products for the recording and evaluation of the personal vital signs. This provides more safety for your own lifestyle.



Simon, 45
Regularly checks his blood pressure.

Regular checks provide a higher degree of security.

Blood pressure is vital to the for the entire body. Regular checks helps to detect possible diseases at an early stage.

Especially for the elderly, the regular control of the blood pressure is a daily necessity, but also a simple and effective health supplement. The upper arm blood pressure monitor BU 535 makes it even easier. A large, easy to read display and above all the additional voice output of the blood pressure values make the measurement a very pleasant exercise. The precise indication of the values, the storage of the data and the calculation of the average value of the last 3 measured values provide a reassuring security.

Dr. med. Johannes Wimmer describes the benefits of the VitDock App, which can be used to better understand the coherency's between the measured values and thus promotes a healthier way of living.



NEW Upper arm blood pressure monitor BU 535 VOICE

- Accurate upper arm blood pressure measurement
- Voice output in 6 languages
- Arrhythmia indicator
- 120 storage locations for 2 users
- Calculation of the average value of the last 3 stored values
- Classification of the WHO* measurement values
- Certified medical product**

Certified Medical Device

WHO

2x120

XL

FLAT DESIGN

3 years Warranty



With VOICE OUTPUT

The blood pressure measurement is simply safe: Excellent readability through extra large digits.



High accuracy leads to reliable measurements.

Measuring the blood pressure, and interpreting it correctly is a routine examination for the doctor. The MTP Pro upper arm blood pressure monitor also makes this possible for the measurement at home.

According to the current definition of the guidelines, hypertension is indicated when either a systolic value of 140 or a diastolic value of 90 is reached or exceeded. With the regular blood pressure measurement, you can control yourself and inform the doctor about deviations at an early stage. Over 1 million people rely on the accurate measurement of the MTP. The application is intuitive and secure, the evaluation of the data is meaningful and reliable according to WHO* criteria.



NEW Upper arm blood pressure monitor MTP Pro

- Accurate blood pressure measurement
- NEW: WHO* indicator
- NEW: Arrhythmia indicator
- Easy one-touch operation
- Optional with large upper arm cuff of 22 - 32 cm
- Calculation of the average value of the all stored values
- Certified medical product**

ETM TESTMAGAZIN · URTEIL	
MEDISANA Blutdruckmessgerät MTP	
GUT	89,3 %
Im Test: 9 Oberarm-Blutdruckmessgeräte Heft 02/2015	
Testurteile: 4x sehr gut, 5x gut	



Simply and safe.

Simon, 45
Everything is totally fine.



NEW Wrist blood pressure monitor BW 335

- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- Display indicator: Systole, diastole, pulse, date, time
- Classification of the WHO measurement values*
- Certified medical product**




Blood pressure measurement on the wrist is a convenient alternative for the measurement on the upper arm.

Everyone should be aware his blood pressure. Due to the large number of measurement values, the regular measurement is an important supplement to the examination by a physician. Thereby, the measurement on the wrist with the BW 335 is the safe and comfortable alternative to the measurement on the upper arm. The positive: The measured values correspond to the blood pressure in daily life and thus provides the doctor with important information.

The simplicity of use makes measurement even more convenient.

* World Health Organisation ** Certified according to MDD (Medical Device Directive)



Maria, 35
Feels comfortable with her figure.

NEW

Glass personal scale PS 425 with voice output

- Voice output in German, English, French and Spanish
- Tread surface made of 6 mm safety glass with rounded corners
- Easy to read LCD-display
- Convenient "step-on" operation



Just a word with you – your perfect weight

The talking personal scale PS 425 is a great help for the elderly, the visually impaired and the blind. It will provide information about the body weight loud and clear.

With
**VOICE
OUTPUT**



The daily weighing process will be greatly eased for elderly people. Impairments regarding the balance or vision make it difficult to weigh. With a tread surface made of 6 mm safety glass and the rounded corners the personal scale PS 425 provides a stable base to stand on. The easy-to-read LCD-display, as well as the voice output can absolutely be relied on when weighing.

The day will turn out good. The weight is right.



INFO Nutrition

Healthy diet –
the best way to control weight.

There are many types of diets, and not all are useful. A balanced and healthy diet is the best means of weight control. Here are some tips: Drink up to 2 liters daily. Eat lots of vegetables and fruit per day. Base your food intake on whole grain and unsaturated fatty acids, and very important: eat conscious and with a broad variety. The scale will show you the success.

The daily weight control is a good contribution to a health-conscious life.

Good mood - the best motivation for stepping onto the scale.

This motivates you on your way to achieving the desired weight. 3 different inspiring photographs are used to change the Motif glass scale.

For many people it is a necessity to regularly control their own weight. This also includes the morning walk to the scale. With the Motif glass scale PS 405 you will always have an overview of your weight. Regardless of whether you want to gain or lose weight. The 4 high-precision DMS sensors provide precise information.

NEW Motif glass scale PS 405

- With optically beautiful exchangeable photo motifs
- Tread surface made of 6 mm safety glass with rounded corners
- Easy to read LCD-display
- Convenient "step-on" operation
- Extra flat design



With exchangeable **PHOTO MOTIFES**



3 modern photographs are available for selection and can be easily exchanged.

- max 180 kg
- LCD
- FLAT DESIGN
- 3 years Warranty



Weighing can be so easy.

Aesthetics and design have long been part of the bathroom. The Design glass scale PS 435 sets a strong accent in the modern wellness-minded bathroom.

Precision and design are not a contradiction, but rather form a beautiful combination with the PS 435 Design glass scale. The 4 DMS sensors ensure the unrivaled pre-

cision of the measurement results, the extra-flat shape, and the special design with the invisible LED digits make the scale an absolute highlight in the bathroom.

NEW Design glass scale PS 435

- Invisible LED-display - only displayed when weighing
- Convenient "step-on" operation
- 4 high precision strain gauge sensors for accurate results

- max 180 kg
- LED
- FLAT DESIGN
- 3 years Warranty



Invisible LED-Display - only displayed when weighing.

I place a great value on style - even when it comes to the accessories in the bathroom.





Simply a safe measurement with infrared.

Benjamin, 6 weeks
Well protected.

Multifunctional thermometer TM 750 connect

- Accurate and very quick measurement on the forehead and ear
- Data transmission via Bluetooth® Smart (4.0) to the VitaDock+ App
- Optical fever alarm due through a colour change from green to red
- Extra function: Measurement of ambient temperature, liquids, surfaces
- Automatic storage of the last 30 measured values
- Certified medical product**





Optical
FEVER-
ALARM

With babies and infants, measuring fever is not always easy. The multifunctional thermometer TM 750 connect allows the uncomplicated measurement on the forehead and in the ear.

Especially with small children, the fever measurement must be fast and absolutely reliable. The Bluetooth multifunctional thermometer TM 750 connect therefore shows, in addition to the body temperature, an elevated temperature through an optical display. As an extra function the device offers the measurement of the ambient temperature, liquids and surfaces - and sends all data via Bluetooth for evaluation to the VitaDock+ App.

Download the VitaDock+ App here:



Contact-free and precise measurement.

Everything fully under control. The contact-free measurement of the body temperature with infrared is very accurate.

Fever and elevated temperatures develop quite differently throughout the course of a day. A regular measurement helps to maintain an overview, and to assess the right time when a doctor should be consulted. The TM 770 stores the last 32 measured values. This makes it possible to precisely follow a development.

Infrared Thermometer TM 770

- Contact-free and very quick measurement on the forehead
- Optical fever alarm: green is ok, red means an elevated temperature
- Automatic storage of the last 32 measured values
- Certified medical product**



Contact-free
THERMOMETER



Finn, 8
Now healthy again.



Margot, 64
Wants precise
measurements.

Blood glucose meter MediTouch 2 connect

- Precise measuring accuracy according to ISO 15197:2013, thanks to a new test strip technology
- Data transmission via Bluetooth® Smart (4.0) to the VitaDock+ App and VitaDock online
- Acoustic tester reminder
- Memory for 480 measurement results
- Certified medical product**



A proper blood glucose control is important for the treatment of Diabetes, in order to avoid subsequent health impairments.




Life is good.
The values are in check.



The blood glucose level should be measured regularly and in a structured fashion. With the acoustic reminder function of the MediTouch 2, the control task will be very easy.

The all-in package for an accurate and simple measurement of the blood glucose value.

The better the awareness of the body's health condition throughout the course of the day, the better the metabolism can be influenced positively. Structured measurements with precise evaluations help to assess the impact of the factors diet, exercise or stress. The MediTouch 2 uses a new test strip technology (GDH-FAD) for the measurement. Additionally, it is possible to evaluate all measurement data comprehensively in VitaDock online, and with the VitaDock+ App.



INFO Diabetes

Why measure blood glucose?

Blood glucose is the amount of glucose in the blood. It provides information regarding the metabolic situation. Deviations of the blood glucose concentration upwards or downwards can provide information on how therapeutic measures are taken effect. With the modern blood glucose meters from MEDISANA, it is possible to conduct a reliable and professional measurement with almost no effort.

Living well despite diabetes.

Upper arm blood pressure monitor Upper arm blood pressure monitor Upper arm blood pressure monitor

BU 550 CONNECT



BU 530 CONNECT



MTP PRO



Our Premium Connect device

- Accurate blood pressure measurement on the upper arm
- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
- 500 memory spaces for each of 2 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

Our Connect device with top price-quality ratio

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

Our million-fold bestseller – now with a WHO scale and arrhythmia indicator

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- Higher accuracy of measurements for severely fluctuating blood pressure thanks to patented 3MAM technology (weighted average of 3 consecutive measurements)

- Arrhythmia indicator
- Very large digits for easy reading
- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Measurement for guests without storage of reading possible
- Battery operation (power adapter optional available)
- Ultra slim design
- Very large cuff (arm circumferences 22-42 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AAA)
- Certified medical device***

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Measurement for guests without storage of reading possible
- Battery operation (power adapter optional available)
- Very large cuff (arm circumferences 22-36 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AAA)
- Certified medical device***

- 99 memory spaces for each of 2 users
- Extra large, easy-to-read display
- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of all stored data
- 1-button operation
- Battery operation (power adapter optional available)
- Large cuff (arm circumferences 22-32 cm)
- Incl. storage bag
- Incl. blood pressure pass
- Incl. batteries (4 x 1.5V AA)
- Certified medical device***

Size: approx. 9 x 14 x 3.3 cm
Weight: approx. 0.35 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40-230 mmHg
 Pulse: 40-199 beats/min.
Batteries: 4 x 1.5V (AAA)
Article No.: 51290
EAN: 40 15588 51290 2
Warranty: 3 Years

Size: approx. 15 x 10 x 6 cm
Weight: approx. 0.4 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 0-299 mmHg
 Pulse: 40-199 beats/min.
Batteries: 4 x 1.5V (AAA)
Article No.: 51174
EAN: 40 15588 51174 5
Warranty: 3 Years

Size: approx. 12 x 16,3 x 6,6 cm
Weight: approx. 0.5 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 30-280 mmHg
 Pulse: 40-200 beats/min.
Batteries: 4 x 1.5 V (AA)
Article No.: 51090
EAN: 40 15588 51090 8
Warranty: 3 Years

Upper arm blood pressure monitor Upper arm blood pressure monitor Upper arm blood pressure monitor



BU 535 VOICE



BU 535



BU 510



Our new ones - now with voice output

- Accurate blood pressure measurement on the upper arm
- Voice output in selectable languages: DE, GB, NL, ES, FR, IT
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading

Our new ones - now with XL display

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading

Our classic

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 90 memory spaces for each of 2 users
- Very large digits for easy reading (up to 22 mm height)

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Very large cuff (arm circumferences 22-36 cm)
- Including storage bag
- Incl. batteries (4 x 1.5V AAA)
- Certified medical device***

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Very large cuff (arm circumferences 22-36 cm)
- Including storage bag
- Incl. batteries (4 x 1.5V AAA)
- Certified medical device***

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of all stored data
- Battery operation (power adapter optional available)
- Very large cuff (arm circumferences 22-36 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AA)
- Certified medical device***

Size: approx. 15 x 9,5 x 4,1 cm
Weight: approx. 0,23 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40-260 mmHg
 Pulse: 40-180 beats/min.
Batteries: 4 x 1,5 V (AAA)
Article No.: 51178
EAN: 40 15588 51178 3
Warranty: 3 Years

Size: approx. ca. 15 x 9,5 x 4,1 cm
Weight: approx. 0,23 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40-260 mmHg
 Pulse: 40-180 beats/min.
Batteries: 4 x 1,5 V (AAA)
Article No.: 51176
EAN: 40 15588 51176 9
Warranty: 3 Years

Size: approx. 13 x 11 x 6 cm
Weight: approx. 0.3 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40-260 mmHg
 Pulse: 40-180 beats/min.
Batteries: 4 x 1.5 V (AA)
Article No.: 51160
EAN: 40 15588 51160 8
Warranty: 3 Years

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3 **World Health Organisation *** Certification according to MDD (Medical Device Directive)

Upper arm blood pressure monitor

MTS



- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 60 memory spaces for each of 2 users
- Very large digits for easy reading
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Automatic switch-off after 1 min.
- Battery operation (power adapter optional available)
- Large cuff (arm circumferences 22-30 cm)
- Incl. storage bag
- Incl. batteries (4 x 1.5V AA)
- Certified medical device***

Size: approx. 19.5×22.5×26.5 cm
Weight: approx. 0.3 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 0-300 mmHg
 pulse: 30-180 beats/min.
Batteries: 4×1.5 V (AA)
Article No.: 51152
EAN: 40 15588 51152 3
Warranty: 3 Years

Wrist blood pressure monitor

BW 300 CONNECT



NEW

Our Premium Connect device for the wrist

- Accurate blood pressure measurement on the wrist
- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
- 180 memory spaces for each of 2 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- Arrhythmia indicator
- Very large digits for easy reading
- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Measurement for guests without storage of reading possible
- Ultra slim design
- Cuff for wrist circumferences 13.5-21.5 cm
- With rechargeable lithium-ion battery and USB charging cable
- Certified medical device***

Size: approx. 6.5×8×1.3 cm
Weight: approx. 0.1 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40-230 mmHg
 pulse: 40-199 beats/min.
Batteries: Lithium ion accumulator
Article No.: 51294
EAN: 40 15588 51294 0
Warranty: 3 Years

Wrist blood pressure monitor

BW 335



Flat design

- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading
- Display items: systolic, diastolic, pulse, date, time
- Ultra slim design

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Cuff for wrist circumferences 14-19.5 cm
- Incl. storage box
- Incl. batteries (2×1.5V AAA)
- Certified medical device***

Size: approx. 8.9×6×3.1 cm
Weight: approx. 0.07 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 0-299 mmHg
 pulse: 40-180 beats/min.
Batteries: 2×1.5 V (AAA)
Article No.: 51077
EAN: 40 15588 51077 9
Warranty: 3 Years

Wrist blood pressure monitor

BW 333



Ultra-flat design with XL cuff

- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Extra large, easy-to-read display
- Display items: systolic, diastolic, pulse, date, time
- Ultra slim design – only 1,5 cm

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Automatic switch-off
- Battery level indicator
- Cuff for wrist circumferences 13.5-23 cm
- Incl. battery (2×1.5V AAA)
- Certified medical device***

Size: approx. 8.4×7.7×1.5 cm
Weight: approx. 0.12 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 30-250 mmHg
 pulse: 40-180 beats/min.
Batteries: 2×1.5 V (AAA)
Article No.: 51075
EAN: 40 15588 51075 5
Warranty: 3 Years

Wrist blood pressure monitor

BW 320



NEW

Elegant design for the top price-quality ratio

- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Extra large, easy-to-read display
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Battery level indicator
- Cuff for wrist circumferences 12.5 - 21.5 cm
- Incl. battery (2×1.5V AAA)
- Certified medical device***

Size: approx. 7.5×6.4×2.9 cm
Weight: approx. 0.17 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 30-280 mmHg
 Puls: 40-200 beats/min.
Batteries: 2×1.5 V (AAA)
Article No.: 51074
EAN: 40 15588 51074 8
Warranty: 3 Years

Wrist blood pressure monitor

BW 315



NEW

Our bestseller


- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Extra large, easy-to-read display
- Display items: systolic, diastolic, pulse, date, time

- Measurements are classified according to the WHO** evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Battery level indicator
- Cuff for wrist circumferences 12.5 - 21.5 cm
- Incl. battery (2×1.5V AAA)
- Certified medical device***

Size: approx. 6.3×6.3×2.9 cm
Weight: approx. 0.11 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 30-280 mmHg
 Puls: 40-200 beats/min.
Batteries: 2×1.5 V (AAA)
Article No.: 51072
EAN: 40 15588 51072 4
Warranty: 3 Years

Body analysis scale

BS 440 CONNECT



Connect device with ITO electrodes


- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- High-quality ITO electrodes
- Convenient „step-on“ switch
- Large easy-to-read LCD-display with backlight
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)
- Maximal load capacity: 180 kg
- Incl. batteries (4×1.5 V AAA)

Size: approx. 32×32×2.7 cm
Weight: approx. 2 kg
Display size: 7×4.9 cm
Weight display: 5–180 kg
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 4×1.5 V (AAA)
Article No.: 40423
EAN: 40 15588 40423 8
Warranty: 3 Years

Body analysis scale

BS 430 CONNECT




- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- High-quality stainless steel electrodes
- Convenient „step-on“ switch
- Large easy-to-read LCD-display with backlight
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)
- Maximal load capacity: 180 kg
- Incl. batteries (4×1.5 V AAA)

Size: approx. 32.5×32.5×2.5 cm
Weight: approx. 2 kg
Display size: 7×4.9 cm
Weight display: 5–180 kg
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 4×1.5 V (AAA)
Article No.: 40422
EAN: 40 15588 40422 1
Warranty: 3 Years

Body analysis scale

BS 444 CONNECT



- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- High-quality stainless steel electrodes
- Convenient „step-on“ switch
- Large easy-to-read LCD-display
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)
- Maximal load capacity: 180 kg
- Incl. batteries (3×1.5 V AAA)

Size: approx. 30×30×2.3 cm
Weight: approx. 1.6 kg
Display size: 7.3×4.8 cm
Weight display: 5–180 kg
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5 V (AAA)
Article No.: 40444
EAN: 40 15588 40444 3
Warranty: 3 Years

Body analysis scale

TARGETSCALE 3



Q3/2017

Now new under VitaDock+

- Precise determination of weight (up to 180 kg), Body Mass Index (BMI), kcal per day, calculates percentage of body fat and water as well as muscle and bone mass
- Multiuser: Automatic detection of up to 8 users
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- Lighted rings depict the difference between current and your target weight
- Large easy-to-read LCD-display with white luminous digits
- Athlete mode for sportspersons
- Stores up to 240 measurements in the scale memory
- Maximal load capacity: 180 kg
- Incl. batteries (4×1.5 V AA)

Size: approx. ø 36×4.75 cm
Weight: approx. 2 kg
Display size: 7×5 cm
Weight display: 5–180 kg
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 4×1.5 V (AA)
Article No.: 40413
EAN: 40 15588 40413 9
Warranty: 3 Years

Design glass scale

PS 435



NEW



Q4/2017

Noble design

- Invisible LED-display (only appears during the weighing)
- Large, easy-to-read LED luminous digits (28 mm)
- Convenient „step-on“ switch
- Energy-saving automatic cut-off after 10 seconds

- Treads of the standing area made of 6 mm extra clear safety glass with rounded corners
- Four high-precision strain gauge sensors for accurate readings
- Extra slim design
- Toggles between kg, lb and st
- Overload indication
- Low battery indicator
- Maximal load capacity: 180 kg
- Incl. batteries (3×1.5 V AAA)

Size: approx. 30×30×2.1 cm
Weight: approx. 1.5 kg
Display size: 7×5 cm
Weight display: 5–180 kg
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5V (AAA)
Article No.: 40434
EAN: 40 15588 40434 4
Warranty: 3 Years

Glass personal scale

PS 425 WITH VOICE OUTPUT



NEW



Q4/2017

Our talking scale

- Voice output in four languages: DE, GB, FR, ES
- Easily readable LCD-display unit
- Convenient „step-on“ switch
- Energy-saving automatic cut-off

- Treads of the standing area made of 6 mm safety glass with rounded corners
- Four high-precision strain gauge sensors for accurate readings
- Toggles between kg and lb
- Overload indication
- Low battery indicator
- Maximal load capacity: 180 kg
- Incl. batteries (3×1.5 V AAA)

Size: approx. 31×31×2.4 cm
Weight: approx. 1.8 kg
Weight display: 4–180 kg
Graduation: 100 g, 0.2 lb
Batteries: 3×1.5V (AAA)
Article No.: 40448
EAN: 40 15588 40448 1
Warranty: 3 Years

Personal scale

PS 430



Anti-slip surface

- Anti-slip surface for safe foothold
- Energy saving auto-off in six seconds
- Easily readable LCD-display

- Convenient „step-on“ switch
- Treads of the standing area made of 6 mm safety glass with rounded corners
- Extra slim design
- Toggles between kg, lb and st
- Overload indication
- Low battery display
- Four high-precision strain gauge sensors for accurate results
- Maximal load capacity: 180 kg
- Incl. battery (1×CR2032)

Size: approx. 30×30×2.1 cm
Weight: approx. 1.5 kg
Display size: 3.3×7.5 cm
Weight display: bis 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Battery: 1×CR2032
Article No.: 40458
EAN: 40 15588 40458 0
Warranty: 3 Years

Motif glass scale

PS 405



NEW



Foto-design glass scale with 2 interchangeable designs

- Energy-saving automatic cut-off
- Easily readable LCD-display
- Tread area of 6 mm extra clear safety glass

- 4 high-precision strain gauge sensors for accurate results
- Convenient „step-on“ switch
- Toggles between kg, lb and st
- Overload indication
- Low battery display
- Extra slim design
- Maximal load capacity: 180 kg
- Incl. batteries (3×1.5 V AAA)



incl. 2 design films for changing

Size: approx. 30×30×2.1 cm
Weight: approx. 1.5 kg
Display size: 2.9×7.4 cm
Weight display: bis 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5 V (AAA)
Article No.: 40405
EAN: 40 15588 40405 4
Warranty: 3 Years

XL personal scale

PS 460



- Particularly large tread area
- Elegant design of high quality stainless steel
- Large easy-to-read LCD-display with red luminous digits
- Maximal load capacity: 200 kg

- Automatic switch-off after 8 seconds
- Convenient „step-on“ switch
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- Incl. batteries (6×1.5 V AAA)

Size: approx. 35×34×2 cm
Weight: approx. 2.6 kg
Display size: 8.8×4.3 cm
Weight display: bis 200 kg, 440 lb, 31 st 6.9 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 6×1.5 V (AAA)
Article No.: 40433
EAN: 40 15588 40433 7
Warranty: 3 Years

Body analysis scale

ISA



- Measurement of weight, body fat, body water and muscle mass
- Determination of the BMI value
- Large easy-to-read LCD-display with backlight

- 10 memory slots
- Automatic switch-off after 8 seconds
- Convenient „step-on“ switch
- Tread area of 8 mm safety glass
- Modern blue LED edge lighting
- High-quality stainless steel electrodes
- Non-slip tread area
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- Maximal load capacity: 180 kg
- Incl. batteries (3×1.5 V AA)

Size: approx. 31.9×31.9×2.15 cm
Weight: approx. 2.6 kg
Display size: 10.5×4.4 cm
Weight display: bis 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5 V (AA)
Article No.: 40480
EAN: 40 15588 40480 1
Warranty: 3 Years

Body analysis scale

PSM



- Measurement of weight, body fat, body water and muscle mass
- Integrated calorie analysis
- Large easy-to-read LCD-display with white luminous digits

- 10 memory slots
- Automatic switch-off after 8 seconds
- Convenient „step-on“ switch
- Tread area of 8 mm safety glass
- High-quality stainless steel electrodes
- Toggles between kg, lb and st
- Four high-precision strain gauge sensors for accurate readings
- Maximal load capacity: 180 kg
- Incl. batteries (3×1.5 V AAA)

Size: approx. 32×32×2.6 cm
Weight: approx. 2.5 kg
Display size: 10.5×4.4 cm
Weight display: bis 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5 V (AAA)
Article No.: 40446
EAN: 40 15588 40446 7
Warranty: 3 Years

Personal scale

PSD



- Trendy retro design
- Battery-free, mechanical operation
- Large full-view analogue dial

- Non-slip tread area
- Solid metal casing

Size: approx. 28.5×42.5×8.5 cm
Weight: approx. 3.6 kg
Weight display: bis 150 kg
Graduation: 500 g
Article No.: 40461
EAN: 40 15588 40461 0
Warranty: 3 Years

Sleep monitor

SLEEPSPACE



- Monitors depth and duration of sleep, heart rate, respiration rate and movements when asleep
- Accurate and comprehensive analysis of sleeping habits

- Ultra-thin piezoelectric band (1.5 mm) for a comfortable fit between the mattress and the sheets
- Patented activation system
- Rechargeable lithium-ion battery
- Incl. USB charging cable

Sleepspace App functions

- Automatic data transfer via Bluetooth Smart to the Sleepspace App*
- Saves and displays all measurement results clearly over time: day, week or month views

Size sensor: approx. 8,8×6,5×1,5 cm
Size band: approx. 79,5×6,5×0,15 cm
Weight: approx. 160 g
Battery: Rechargeable 3,7V lithium-ion battery
Article No.: 79490
EAN: 40 15588 79490 2
Warranty: 3 Years



INFO Sleep monitor Sleepspace

Sleep laboratory for use at home

In comparison to measurements in a medical sleep laboratory, Sleepspace provides just as reliable data for a sleep analysis at home. The only difference: you are in a familiar environment and you can carry out the sleep analysis conveniently over a prolonged period of time.

Infrared thermometer

TM 770



NEW



Q1/2018

- Precise measurement of the body temperature on the forehead in seconds
- Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100° C)
- Visual fever alarm with colour change – green: temperature normal, red: raised temperature

- Acoustic signal at the beginning and end of temperature reading
- Automatic storage of the last 32 measurement results
- Automatic power off after 30 seconds
- Unit selectable: °C or °F
- Incl. battery (2×1,5 AA)
- Certified medical device**

Size: approx. 9,3×15,3×4,1 cm
Weight: approx. 160 g
Measuring range: Body: 32-42.9°C
 Object: 0-100°C
Accuracy: Forehead: ± 0,3°C (32-42.9°C)
 Object: ± 0,1°C (0-100°C)
Battery: 2×1.5V (AA)
Article No.: 76150
EAN: 40 15588 76150 8
Warranty: 3 Years

Multifunctional thermometer

TM 750 CONNECT



NEW



Q1/2018

- Connect device**
- Precise measurement of the body temperature in the ear and on the forehead in seconds
 - Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100° C)
 - Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- Automatic storage of the last 30 measurement results
- Visual fever alarm with colour change – green: temperature normal, red: raised temperature
- Acoustic signal at the beginning and end of temperature reading
- Display of date and time
- Incl. storage pouch
- Incl. batteries (2×1.5V AAA)
- Certified medical device**

Size: approx. 13,4×4,2×4,7 cm
Weight: approx. 65 g
Measuring range: Body: 34-43°C
 Object: 0-100°C
Accuracy: Forehead: ± 0,3°C (34-43°C)
 Ear: ± 0,3°C (35,5-42°C)
 Object: ≤30 ±1.5°C / >30 ±5%
Batteries: 2×1.5V (AAA)
Article No.: 76145
EAN: 40 15588 76145 4
Warranty: 3 Years

Infrared multi-functional thermometer

TM 750



- Precise measurement of the body temperature in the ear and on the forehead in seconds
- Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100° C)
- Automatic storage of the last 30 measurement results

- Visual fever alarm with colour change – green: temperature normal, red: raised temperature
- Acoustic signal at the beginning and end of temperature reading
- Display of date and time
- Incl. storage pouch
- Incl. batteries 2×1.5V (AAA)
- Certified medical device**

Size: ca. 13,4×4,2×4,7 cm
Weight: ca. 65 g
Measuring range: Body: 34-43°C
 Object: 0-100°C
Accuracy: Forehead: ± 0,3°C (34-43°C)
 Ear: ± 0,3°C (35,5-42°C)
 Object: ≤30 ±1.5°C / >30 ±5%
Batteries: 2×1.5V (AAA)
Article No.: 76140
EAN: 40 15588 76140 9
Warranty: 3 Years

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3 ** Certification according to MDD (Medical Device Directive)

Infrared thermometer

FTN



Contactless fever thermometer

- Precise non-contact measurement of the body temperature
- Additional function: measurement of ambient temperature, liquids and surfaces
- Measurement in a few seconds

- Automatic storage of the last 30 measurement results
- Visual fever alarm with colour change in the display
- Incl. storage box
- Incl. battery (2×1.5V AAA)
- Certified medical device**

Size: approx. 14.7×3.8×2.1 cm
Weight: approx. 48 g
Measuring range: Körper: 34-42.2°C
 Objekt: 0-100°C
Accuracy: 0.1°C
Battery: 2×1.5V (AAA)
Article No.: 76120
EAN: 40 15588 76120 1
Warranty: 3 Years

Bluetooth thermometer

TM 740 CONNECT



NEW



Q2/2018

Ideal for babies and toddlers

- Continuous real-time monitoring of body temperature
- Coloured LEDs on the thermometer indicate the temperature range – green: normal temperature, Yellow: increased temperature., Red: Fever (from 38 °C), Blue: Bluetooth connection to the smart-phone has been set up
- Bluetooth® Smart data transfer to the VitaTherm app for iOS and Android*

- Very flat thermometer for pleasant wearing comfort
 - Where to use: on the skin, in the vicinity of the armpit
 - Universal user group: for small children to adults
 - Incl. battery (1×CR2032)
 - Certified medical device**
- Functions of the app**
- Neatly arranged display of measurement values in chronological order
 - Alarm is triggered when exceeding the set temperature zone
 - Reminding function for intake of medication
 - Alarm when Bluetooth® connection is disconnected (radius max. 10 m)
 - Unit of measurement can be selected (°C or °F)



Size: approx. 4.2×4.2×0.9 cm
Weight: approx. 16 g
Measuring range: 32-42°C
Accuracy: ± 0.1°C (35-42°C)
Battery: 1×CR2032
Article No.: 76163
EAN: 40 15588 76163 8
Warranty: 3 Years

Bluetooth thermometer

TM 735 CONNECT



- Continuous real-time monitoring of body temperature
- Universal user group: for small children to adults
- Bluetooth® Smart data transfer to the ThermoApp for iOS and Android*

- Where to use: on the skin, in the armpit
 - If using on children below 2 years of age, consult a doctor
 - Incl. battery (1x CR2025)
 - Certified medical device**
- Functions of the app**
- Neatly arranged display of measurement values in chronological order
 - Alarm is triggered when exceeding the set temperature zone
 - Reminding function for intake of medication
 - Alarm when Bluetooth® connection is disconnected (radius max. 10 m)
 - Unit of measurement can be selected (°C or °F)



Size: approx. 6×3×0.6 cm
Weight: approx. 16 g
Measuring range: 25-45°C
Accuracy: ± 0.05°C (35-38.5°C) und ± 0.1°C (< 35°C und > 38.5°C)
Battery: 1×CR2025
Article No.: 76165
EAN: 40 15588 76165 2
Warranty: 3 Years

Thermometer

TM 705



NEW



Q2/2018

Our classic. Now with traffic-light style scale and fever indicator

- Precise measurement of the body temperature: oral, axillary, rectal
- Flexible tip
- Measurement in 10 seconds
- Water resistant
- With fever indicator: green: 35.8–37.1 °C; yellow: 37.2–37.8 °C; red: 37.9–42.9 °C

- Acoustic fever alarm
- Acoustic signal at the beginning and end of temperature reading
- Automatic storage of the last measurement result
- Automatic power off after 10 minutes
- Digital LCD-display
- Conversion from °C to °F
- Easy battery replacement
- Incl. storage box
- Incl. battery (1x 1.5V AG3/LR41)
- Certified medical device**

Size: approx. 13.9×2.2×1.3 cm
Weight: approx. 11 g
Measuring range: 32-42.9°C
Accuracy: 0.1°C
Battery: 1x 1.5V (AG3/LR41)
Article No.: 77042
EAN: 40 15588 77042 5
Warranty: 3 Years

Thermometer

TM 700



Our gentle classic

- Precise measurement of the body temperature: oral, axillary, rectal
- Flexible tip
- Measurement in 10 seconds
- Water resistant

- Acoustic fever alarm
- Acoustic signal at the beginning and end of temperature reading
- Automatic storage of the last measurement result
- Automatic power off after 10 minutes
- Digital LCD-display
- Conversion from °C to °F
- Easy battery replacement
- Incl. storage box
- Incl. battery (1x 1.5V AG3/LR41)
- Certified medical device**

Size: approx. 13.9×2.2×1.3 cm
Weight: approx. 11 g
Measuring range: 32-42.9°C
Accuracy: 0.1°C
Battery: 1x 1.5V (AG3/LR41)
Article No.: 77040
EAN: 40 15588 77040 1
Warranty: 3 Years

Thermometer

FTC



Top price-quality ratio

- Precise measurement of the body temperature: oral, axillary, rectal
- Acoustic signal at the beginning and end of temperature reading
- Water resistant

- Acoustic fever alarm
- Automatic storage of the last measurement result
- Automatic power off after 8 minutes
- Digital LCD-display
- Easy battery replacement
- Incl. storage box
- Incl. battery (1x 1.5V AG3/LR41)
- Certified medical device**

Size: approx. 12.8×1×1.9 cm
Weight: approx. 10 g
Measuring range: 32-42°C
Accuracy: 0.1°C
Battery: 1x 1.5V (AG3/LR41)
Article No.: 77030
EAN: 40 15588 77030 2
Warranty: 3 Years

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3 ** Certification according to MDD (Medical Device Directive)

Blood glucose meter

MEDITOUCH 2 CONNECT




Blood glucose meter

MEDITOUCH 2 CONNECT



Blood glucose meter

MEDITOUCH 2



Connect device in a carbon design

- Precise measurement accuracy according to the new ISO 15197: 2013 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- Acoustic test reminders
- 480 memory slots
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

Connect device

- Precise measurement accuracy according to the new ISO 15197: 2013 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- Acoustic test reminders
- 480 memory slots
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

Our classic

- Precise measurement accuracy according to the new ISO 15197: 2013 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- Acoustic test reminders
- 480 memory slots
- USB port for data transfer to VitaDock Online and to the VitaDock+ app for iOS and Android

- Measuring unit individually adjustable: mg/dL or mmol/L
- USB port for data transfer to VitaDock Online
- Pre and post meal markers
- Automatic ketone warning
- Measuring time within 5 seconds
- 0,6 µl blood
- Extra large, easy-to-read display
- Incl. batteries (2×3V CR2032)
- Certified medical device**

- USB port for data transfer to VitaDock Online
- Pre and post meal markers
- Automatic ketone warning
- Measuring time within 5 seconds
- 0,6 µl blood
- Extra large, easy-to-read display
- Incl. diabetic diary
- Incl. batteries (2×3V CR2032)
- Certified medical device**

- Calculates average values for 7/14/30/90 days
- Pre and post meal markers
- Automatic ketone warning
- Measuring time within 5 seconds
- 0,6 µl blood
- Extra large, easy-to-read display
- Incl. batteries (2×3V CR2032)
- Certified medical device**



Starter set incl. 10 test strips, lancing device, AST cap, 10 lancets, control solution and bag

Starter set incl. 10 test strips, lancing device, AST cap, 10 lancets, control solution and bag

Starter set incl. 10 test strips, lancing device, AST cap, 10 lancets, control solution and bag

Size: approx. 5×9.8×1.1 cm
Weight: approx. 47 g
Batteries: 2×3V (CR 2032)
Article No.: 79062
EAN: 40 15588 79062 1
Warranty: 3 Years

Size: approx. 5×9.8×1.1 cm
Weight: approx. 47 g
Batteries: 2×3V (CR 2032)
mg/dL Article No.: 79044
EAN: 40 15588 79044 7
mmol/L Article No.: 79046
EAN: 40 15588 79046 1
dual Article No.: 79048
EAN: 40 15588 79048 5
Warranty: 3 Years

Size: approx. 5×9.8×1.1 cm
Weight: approx. 30 g
Batteries: 2×3V (CR 2032)
mg/dL Article No.: 79030
EAN: 40 15588 79030 0
mmol/L Article No.: 79034
EAN: 40 15588 79034 8
Warranty: 3 Years

Pulse oximeter

PM 150 CONNECT



Pulse oximeter

PM 100



- Pulse oximeter to measure the saturation of oxygen in blood (SpO2 %) and the heart rate (pulse)
- View measurements on a LED display
- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*

- Pulse oximeter to measure the saturation of oxygen in blood (SpO2 %) and the heart rate (pulse)
- View measurements on a OLED display
- Display brightness adjustable over 10 levels
- Simple „one touch“ operation – also ideal for travelling

- Simple „one touch“ operation – also ideal for travelling
- Automatic switch-off after 8 seconds
- Incl. battery (1x 1.5V AAA)
- Certified medical device**

- Oxygen saturation (SpO2) and pulse rate can be displayed as bars or waves
- Viewable from 6 different angles
- Automatic switch-off after 8 seconds
- Battery level indicator
- Carrying strap included
- Incl. batteries (2×1.5V AAA)
- Certified medical device**

Size: approx. 7.6×5.6×3.1 cm
Weight: approx. 55 g
Measuring range: SpO2: 70% ~ 99%
 pulse: 30-235 beats/min.
Battery: 1×1.5V (AAA)
Article No.: 79457
EAN: 40 15588 79457 5
Warranty: 3 Years

Size: approx. 5.8×3.4×3.5 cm
Weight: approx. 53 g
Measuring range: SpO2: 70% ~ 99%
 pulse: 30-235 beats/min.
Batteries: 2×1.5V (AAA)
Article No.: 79455
EAN: 40 15588 79455 1
Warranty: 3 Years



The journey is the destination.
And progress its companion.

Sporty assistants – Made for Life.

Alina, 31
Is focused on her goals.

Tina, 27
Is monitoring her calorie consumption.

Sport and fitness have long been accepted as important mainstream factors of a healthy society. For many people, sport is more than just a leisure activity. It is an integral part of someone's life style, and is just a part of the daily routine as eating, drinking and working. In addition, the targeted optimization of the own performance and fitness is currently a very powerful trend, and demonstrates the new self-confidence when it comes to the subject of health. MEDISANA utilizes innovative fitness and medical devices to promote this ambition.

Relaxation and recovery.

The fascia and trigger points can now be self-treated even more precisely.

Muscle tensions are not just annoying and painful after the workout, but they can also harden if they are not treated properly. With the PowerRoll-Ball it is possible to effectively release muscle and fascia tensions. With four gentle intensity levels, you can also regenerate the deeper muscle sections.



70 minutes running time for an effective muscle training



INFO

Fascia training

Smooth fascia are the basis for fitness and a good posture. Because the fascia envelope the muscles. They transmit the muscle power and communicate with the nervous system. They are the sensory organs for the movement, and also protect the muscles. Training your fascia is also a strengthening of the muscles.

The PowerRoll-Ball helps to accurately relieve hardened points in deeper muscle areas

MEDISANA Sport-Sponsorship



Svenja Thoes,
Triathlete, 26 years

3rd place: Ironman 70.3 Calgary
3rd place: Challenge Heilbronn
1st place AG, and 15th in the women's overall Ironman 140.6 Kona, Hawaii



Andreas Dreiz,
Triathlete, 28 years

2nd place: 70.3 Mallorca
4th place: The Championship 2017
2nd place: Ironman 70.3. Oceanside 2016
1st place: Ironman 70.3. Mallorca 2016

Svenja Thoes, 26 Triathlete
Knows what effectively helps her against tensions.



NEW PowerRoll-Ball

- With vibration - also reaches the deeper muscles
- 4 intensity levels
- LED-display
- Targeted fascia training



Provides a long-term relaxation and stimulation.

What can help against tensions or blockages in the muscles?
The targeted training with the PowerRoll will yield a lasting effect.

NEW PowerRoll (31 cm)
PowerRoll-XT (45 cm)

- With vibration – also reaches the deeper muscle areas
- Includes training DVD & training plan
- 4 intensity levels
- 2 integrated batteries

Hard 4
3 years Warranty

New functional surface



Anyone who demands a lot from his body during the workout, should also provide it with optimal support. The PowerRoll activates and regenerates the muscles, but it's also suitable for strength and stability exercises. Through this, the entire muscular system is optimally supported will beyond the training phase. The PowerRoll is available in lengths of 31 and 45 cm.



A trained body feels good.



Andreas Dreitz, 28 Triathlete
Knows the effective secret:
no tensions, fast muscle regeneration.



The PowerRoll optimally prepares the body for the training.

Soft vibrations
and fascia training in one.



This will make training fun again.

The holistic fascia training with the PowerRoll-Slim acts specifically against adhesions and local hardening of the muscles.

The Massage roller has different surfaces and degrees of hardness, with which the feeling of massaging fingers and palms is imitated realistically. Supported by the vibration effect, you can perform specific exercises for the relaxation and activation of the fascia and muscles.



With its extra pouch you can even reach inaccessible muscle region.

NEW PowerRoll-Slim

- Targeted Fascia training
- 4 intensity levels
- With 3-stage capacity display (battery indicator)



With the MEDISANA PowerRoll-Slim it is possible to train the joints, back, spine, legs, neck and buttocks in a targeted manner.



Get rid of anything, that causes tension.

Celine, 19
Notices successes.

A single unit for many applications. The VarioBall for the activation and regeneration.

Every body is different. Luckily, the VarioRoll can be adjusted through a clearance range of 20 to 24 cm between the balls, for a perfect adaption to the different muscle areas. Tensions, trigger points and adhesions are resolved faster and easier.


VarioRoll

- For the highly accurate self-massage, for example, back muscles, neck, etc.
- For the specific application of fasciae, trigger points, muscles
- Especially gentle to the spine thanks to a pull-out mechanism
- Includes training plan



NEU EMS Body Trainer

- To build up the muscles
- Modern, sporty design
- Auto programs: Massage, acupuncture, tapping, combination
- Set: 3 devices with a pad for the belly and 2 pads for upper arm/leg



The EMS Body Trainer also works on individual muscle areas in a targeted manner.

Quickly achieve more.

Muscles support and strengthen the body. They are the basic pillars of the bodies mobility at any age. With the EMS stimulation method, it is possible to gently build up muscles.

In physiotherapy, the electrical stimulation has long been used for the muscle build up and regeneration. With the EMS Body Trainer you can now also do this at home. The muscles on the abdomen, legs or arms are stimulated with 15 intensity levels. Through this, the regular training can be optimally supplemented.



**Foam roller
POWERROLL-SLIM**



Q4/2017

- PowerRoll-Slim with depth vibration – also reaches the deeper muscles
- For activation and regeneration as well as for strength and stability exercises and fasciae training
- 4 intensity levels
- Simple One-touch operation
- Hardness grade: medium

- LED indication of the selected intensity level
- Surface structure supports the massage action
- Surface is resistant to moisture
- Integrated rechargeable Lithium-Ion battery (charger included)
- Max. operating time with fully charged battery: ca. 180 min. (depending on the selected intensity level)
- Auto-Stop after 15 minutes of operation
- Capacity information



Includes pouch and training plan

Size: approx. 31 cm, Ø 8,5 cm
Weight: approx. 0,8 kg
Battery: Lithium-Ion battery
Article No.: 79534
EAN: 40 15588 79534 3
Warranty: 3 Years

**Foam roller
POWERROLL**



Q1/2018

New design

- PowerRoll with depth vibration – also reaches the deeper muscles
- For activation and regeneration as well as for strength and stability exercises and fasciae training
- 4 intensity levels
- Simple One-touch operation
- Hardness grade: hard

- LED indication of the selected intensity level
- Surface structure supports the massage action
- Surface is resistant to moisture
- Integrated rechargeable Lithium-Ion battery (charger included)
- Max. operating time with fully charged battery: ca. 180 min. (depending on the selected intensity level)
- Auto-Stop after 10 minutes of operation
- Capacity information



Includes training DVD and training plan

Size: approx. 31 cm, Ø 15 cm
Weight: approx. 0,75 kg
Battery: Lithium-Ion battery
Article No.: 79530
EAN: 40 15588 79530 5
Warranty: 3 Years

**Foam roller
POWERROLL-XT**



Q1/2018

New design & extra long

- PowerRoll-XT with depth vibration – also reaches the deeper muscles
- For activation and regeneration as well as for strength and stability exercises and fasciae training
- 4 intensity levels
- Simple One-touch operation
- Hardness grade: hard

- LED indication of the selected intensity level
- Surface structure supports the massage action
- Surface is resistant to moisture
- Integrated rechargeable Lithium-Ion battery (charger included)
- Max. operating time with fully charged battery: ca. 180 min. (depending on the selected intensity level)
- Auto-Stop after 10 minutes of operation
- Capacity information



Includes training DVD and training plan

Size: approx. 45 cm, Ø 15 cm
Weight: approx. 0,9 kg
Battery: Lithium-Ion battery
Article No.: 79532
EAN: 40 15588 79532 9
Warranty: 3 Years



**Foam roller
VARIOROLL**



Ultra flexible for any body position

- Particular suitable for punctual self-massage – for example, of the back, neck etc.
- For the targeted application of fascia, trigger points and muscles
- Very gentle on the spinal column thanks to 3-stage size setting – choosable distance between the spheres (20 – 24 cm)
- For extended applications such as the arms, legs etc.

- For activation and regeneration as well as for strength and stability exercises and fasciae training
- Resistant to moisture
- Designed in Germany



Includes training plan

Size: approx. 20 x 10 x 10 cm, Ø 15 cm
Weight: approx. 0,13 kg
Article No.: 79516
EAN: 40 15588 79516 9
Warranty: 3 Years



















Fitness mat



- Fitness mat for MEDISANA massage roller
- Impact-cushioning (shock-absorbent), soft surface
- Slightly corrugated surface

- Including a strap for easy transport
- Surface resistant to moisture – hygienic to wipe

Size: approx. 60 x 190 x 1,5 cm
Weight: approx. 0,85 kg
black Article No.: 79477
EAN: 40 15588 79477 3
Warranty: 3 Years

Foam ball POWERROLL-BALL	EMS pad BODY TRAINER	Activity tracker VIFIT CONNECT MX3	Activity tracker VIFIT TOUCH	Activity tracker VIFIT TOUCH	Activity tracker VIFIT TOUCH
					
					
<p>Q4/2017</p>	<p>Q4/2017</p>				
<p>Modern design</p> <ul style="list-style-type: none"> Particular suitable for punctual self-massage – for example, of the back, neck etc. For the targeted application of fascia, trigger points and muscles For extended applications such as the arms, legs etc. 4 different intensity levels Capacity information 	<ul style="list-style-type: none"> To build the muscles Modern, sporty design Auto programs: Massage, acupuncture, tapping, combination Set including 3 Pads (1 Abdomen pad, 2 Arm pads) 	<ul style="list-style-type: none"> Activity and Sleep Tracker in one High-quality OLED-display Easy 1-button operation Stores the results of the last 7 days Flight mode (Bluetooth deactivated) Splash proof Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android* 	<ul style="list-style-type: none"> Activity and Sleep Tracker in one High-quality OLED display Displays incoming calls or messages Stores the results of the last 15 days Individual daily goals and movement reminder Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android* 	<ul style="list-style-type: none"> Activity and Sleep Tracker in one High-quality OLED display Displays incoming calls or messages Stores the results of the last 15 days Individual daily goals and movement reminder Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android* 	<ul style="list-style-type: none"> Activity and Sleep Tracker in one High-quality OLED display Displays incoming calls or messages Stores the results of the last 15 days Individual daily goals and movement reminder Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*
<ul style="list-style-type: none"> LED indication of the selected intensity level For activation and regeneration as well as for strength and stability exercises Resistant to moisture Designed in Germany Timer: 10 min. Volume: <60dBA Runtime: 120 minutes at level 4 	<ul style="list-style-type: none"> Intensity adjustable in 15 steps Automatic shut-off after 20 min. PU leather pads Including 3 Controllers and 3 batteries (CR2032)  <p>Includes 2 Pads for arms and legs</p>	<ul style="list-style-type: none"> Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in % Sleep mode: Tracks movement activity during sleep and duration of sleep Motion sensor for activity measurement: counts actual steps Vibration alarm clock Scratch and abrasion resistant housing Battery life: up to 7 days Rechargeable Lithium-Ion battery Display of time, date and battery indicator  <p>Includes 3 wristband in 2 lengths (black, aqua, coral)</p>	<ul style="list-style-type: none"> Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in % Sleep mode: Tracks movement activity during sleep and duration of sleep Motion sensor for activity measurement: counts actual steps Splash proof Vibration alarm clock Battery life: 5 – 7 days Rechargeable Lithium-Ion battery Display of time, date and battery indicator Including USB cable 	<ul style="list-style-type: none"> Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in % Sleep mode: Tracks movement activity during sleep and duration of sleep Motion sensor for activity measurement: counts actual steps Splash proof Vibration alarm clock Battery life: 5 – 7 days Rechargeable Lithium-Ion battery Display of time, date and battery indicator Including USB cable 	<ul style="list-style-type: none"> Activity mode: Tracks steps, distance, burned calories, activity time in hours and daily target in % Sleep mode: Tracks movement activity during sleep and duration of sleep Motion sensor for activity measurement: counts actual steps Splash proof Vibration alarm clock Battery life: 5 – 7 days Rechargeable Lithium-Ion battery Display of time, date and battery indicator Including USB cable
<p>Size: Ø 13 cm Weight: approx. 0,82 kg Battery: Lithium-Ion battery Article No.: 79510 EAN: 40 15588 79510 2 Warranty: 3 Years</p>	<p>Size: approx. 17 x 18,5 cm (Abdomen pad) approx. 5 x 20 cm (Arm pad) Weight: approx. 30 g (Abdomen pad) approx. 13 g (Arm pad) Battery: 1 x CR2032 Article No.: 79524 EAN: 40 15588 79524 4 Warranty: 3 Years</p>	<p>Size: approx. 2 x 1 x 21,5 cm (short wristband) approx. 2 x 1 x 25,2 cm (long wristband) Weight: approx. 22 g (incl. wristband) Battery: Lithium-Ion battery Article No.: 79790 EAN: 40 15588 79790 3 Warranty: 3 Years</p>	<p>Size: approx. 23 x 2 x 0,9 cm Weight: approx. 22 g Battery: Lithium-Ion battery Article No.: 79486 EAN: 40 15588 79486 5 Warranty: 3 Years</p>	<p>Size: approx. 23 x 2 x 0,9 cm Weight: approx. 22 g Battery: Lithium-Ion battery Article No.: 79487 EAN: 40 15588 79487 2 Warranty: 3 Years</p>	<p>Size: approx. 23 x 2 x 0,9 cm Weight: approx. 22 g Battery: Lithium-Ion battery Article No.: 79488 EAN: 40 15588 79488 9 Warranty: 3 Years</p>

* Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3

Enjoy life, totally relaxed.

Sophisticated wellness products – Made for Life.



Sarah, 49
Likes it warm
and cosy.



George, 52
Likes to feel invincible.

The last few years have shown that wellness has become a component of everyday life for many people. The aim is a well equalized life-balance, in which the body and spirit exist in harmony. Enjoying your life, taking advantage of the perks, while doing something positive for your own health at the same time. The demand for wellness products is correspondingly high. As a market leader in the field of massage products, MEDISANA supports the pursuit of contentment and well-being through versatile innovations.



Susanne, 32
Relaxed like never before.

NEW Shiatsu massage seat cover MC 826

- 3 massage zones selectable for the shiatsu massage: entire, upper or lower back
- Acupressure massage with 3 intensity levels in the waist and thigh area
- Separate, height-adjustable neck massage
- Switchable red light and heating function

ELECTRICALLY HEIGHT ADJUSTABLE
NECK MASSAGE



- HD 3D
- washable
- Red Light + Heat
- 3 years Warranty

Vitalizes the good body feeling.

The MC 826 Shiatsu massage seat cover is able to convince with its three highly effective massage techniques: Shiatsu, tapping and acupressure massage.

After a stressful day it's very soothing to just lean back and get a massage from the shiatsu massage seat cover. This helps to relieve tensions and to promote the feeling of well-being throughout the entire body. The new tapping massage works in the back area, the 3-stage vibration massage acts on the buttocks and thighs, and the blood circulation of the waist is stimulated with a gentle acupressure technique.

Follows the shape of your body.

During the massage, the MC 824 supports the back in the truest sense of the word. Thanks to the innovative swivel technology, the massage cover adapts to the spine during movements.

15 minutes should absolutely be invested into an effective Shiatsu spot massage. Thanks to the innovative swivel mechanism, the Shiatsu massage seat cover MC 824 follows the shape of your body. The massage effect is thereby intensified.




NEW Shiatsu massage seat cover MC 824

- Pleasant massage along the spine with new swivel function
- 3 massage zones selectable for the shiatsu massage: total, upper or lower back
- Intensive shiatsu, roller or spot massage selectable in the back area
- Switchable red light and heating function
- Automatic shut-off after 15 minutes

INFO
Swivel function

The innovative swivel mechanism follows during the massage the shape of your body.



- HD 3D
- washable
- Red Light + Heat
- 3 years Warranty

Gentle to the muscles. Relaxing for the back.

The alternation of the massage movements increases the effect. The highly precise tapping massage is supplemented by the circular Shiatsu massage of the MC 818.

Christine, 23
Is always relaxed.

Through its highly precise pressure massage action, the tapping massage promotes the blood circulation, and also supports the general cell regeneration, as well as the relaxation of the back muscles.

NEW Comfort massage seat cover MC 818

- Tapping massage for the entire back with 2 intensity levels
- 3 massage zones selectable for the shiatsu massage: entire, upper or lower back
- Selective spot massage selectable
- Switchable red light and heating function



When the back is in need of relaxation.

MEDISANA has a superior expertise in Shiatsu massage seat covers. With a strong demand on the market, the massage seat cushions are able to impress with a targeted and intensive depth effect.

NEW

Shiatsu massage seat cover MC 815

- 3 massage zones can be selected: entire, upper or lower back
- Separate vibration massage in the seat surface
- 4 rotating massage heads
- Switchable red light and heating function
- Easy care material




Back pain is a national disease. Most of the pain is not acute, but could ultimately lead to tensions or unhealthy postures. The MC 815 provides relaxation for the back.

Martin, 42
Deeply relaxed like never before.



Our top seller with acupressure function.

Consumer organizations have tested the Shiatsu acupressure massage seat cover MC 825 with a grade of "good". Many of the MEDISANA customers can enjoy the relaxing effect on the back and entire body.



Shiatsu acupressure massage seat cover MC 825

- 3 massage zones selectable for the shiatsu massage: entire, upper or lower back
- Acupressure massage with 3 intensity levels in the waist
- Separate, height-adjustable neck massage
- Switchable red light and heating function



The Shiatsu massage is known for its deep relaxing effect. Combined with the new acupressure massage for the waist, as well as a height-adjustable neck massage, the MC 825 effectively counteracts stubborn muscle hardening.

Removes tensions from the back.

Tensions can lead to protective postures, and thus to incorrect postures of the back in the long run. The Shiatsu massage seat cover MCN is there to counteract this.

Movement reliably prevents tensions. The MCN is highly supportive here. Just 15 minutes provide a noticeable relaxation in the entire back.

Our **BESTSELLER**

Shiatsu massage seat cover MCN

- 3 massage zones selectable for the shiatsu massage: entire, upper or lower back
- Separate, neck massage
- Switchable red light and heating function
- Removable washable neck cover



Relaxed feet all day long.

A vitalizing timeout for your feet. The effective foot massage activates the circulatory system, and acts highly pleasant on the muscle areas to be treated.

The feet have nerve endings, which act in very different places throughout the body. Through this, the feet are in a permanent interaction with the entire body. The FM 885 foot massager stimulates these important nerve points, and will therefore provide you and your body with relaxing and soothing many moments.

NEU Shiatsu foot massager FM 885

- Intensive foot reflexology massage
- Shiatsu massage of the foot soles and heels
- 3-stage compression massage of the foot top sides
- Switchable heat in the Shiatsu massage heads
- Acts blood circulation-promoting and invigorating



Lena, 43
Is gaining new energy.



Recreation for the feet.

Someone who sits, stands or walks all day, will feel his feet in the evening. They'll be tired and "stressed". An invigorating foot reflexology massage can perform "miracles".

The foot reflexology massage is a traditional form of therapy that relieves pain, improves blood flow and stimulates well-being throughout the body. The FM 888 foot massager combines different reflexology massages, such as the shiatsu and compression massage, thus providing a holistic feeling that radiates from the feet into the entire body.



Susanne, 38
Appreciates the relaxing effect.

NEW Foot massager HM 888

- Intensive foot reflexology massage
- 2 different massage types selectable
- Shiatsu massage of the foot soles and heels
- Compression massage of the foot top sides
- Suitable up to shoe size 46



62
 Relaxation for
 inbetween-and on the go

The neck and shoulders quickly become tensioned after a long day in the office. The Heat cape HP 626 quickly releases tensions thanks to the Turbo-Heat heating function, and offers a lasting relaxing effect.



A blessing: To freely move the neck and shoulders.

Mobile
 HEAT

NEW

Heat cape HP 626

- Heat cape with rechargeable batteries for indoor and outdoor use Application
- 4 temperature levels
- Battery life approx. 5 hrs.
- Ideal for home or on the road use



Heat where it is needed. The HP 626 distributes the heat evenly across the neck and shoulder thanks to the matrix sensor system. The Turbo-Heat heating function can be switched on with powerful 100 watts in order to release the tensions as quickly as possible. The relaxation effect can be felt immediately.

Heat contributes
 to quality-of-life.

For many people, warmth is a prerequisite for wellbeing. The heating blanket with a cosy factor distributes the heat quickly and evenly – thus ensuring a maximum relaxation effect.



Maren, 39
 Loves cold winter days.

Warmth relaxes, and releases tension as well as blockages. It makes you feel good about yourself. The matrix sensor system and the flexible filament provide a uniform heat distribution. Thanks to 120

watts, the heating blanket heats up very quickly, and can be adjusted precisely to the comfort temperature thanks to the 4 temperature settings.

NEW

Cosy heating blanket
 XXL HB 675

XXL
 SIZE
 200 x 150 cm

- Flexible heating wire and matrix sensor system for an even heat distribution
- Fast warm-up with 120 watts
- Safety through an overheating protection
- Reversible cover in a cuddly 2-color design
- Öko-Tex® Standard 100



**Virtual reality-
Shiatsu massager chair**

SVR 450



NEW

VRX



Q2/2018

The new massage experience for the home

- Novel virtual reality massage experience. Integrated Bluetooth function to connect virtual reality goggles.
- Control via cable remote control or via Bluetooth
- Diverse massage: Shiatsu, tapping, rolling massage and vibration massage in the seat pad
- Targeted spot massage in the back area (Shiatsu and tapping massage)

- High quality design, easy to clean material
- Loadable up to 150 kg
- High-definition massage: targeted, intensive massage for tangible sense of well-being
- 3 selectable massage zones: total, upper or lower back
- Additional heat function
- High backrest
- Handy remote control
- Pleasant seating comfort
- An eye-catcher in every apartment-even without using the massage

Size: 88 x 67 x 107 cm
Weight: 15,5 kg
Power supply: 220-240V~ 50/60Hz
 12V DC, 3.000mA
Power output: 50 Watts
Artikel-No.: 88720
EAN: 40 15588 88720 8
Warranty: 3 Years



INFO **Virtual reality-
Shiatsu massager
chair SVR 420**

Holistic relaxation with tapping and acupressure massage

A highlight in every apartment. The Virtual reality-Shiatsu massage chair SVR 420 combines 4 different massage techniques, and a virtual reality glasses into a holistic massage experience. The relaxation effect is intensified by the high-definition massage.

**Virtual reality-
Massage seat cover**

VR 850



NEW

VRX



Q1/2018

The new massage experience for the home

- Novel virtual reality massage experience. Integrated Bluetooth function to connect virtual reality goggles.
- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Acupressure massage in the waist area with 3 selectable intensity levels
- 2 massage programs for the back selectable (shiatsu and tapping massage)

- Neck and shoulder massage (Height-adjustable)
- Separate vibration massage in the seat with 3 selectable intensity levels
- Targeted spot massage only for the tapping massage
- For tapping massage 3 selectable programs as well as in 3 intensity levels selectable
- Automatic height-adjustable neck massage
- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Relaxed massage for the entire upper or lower back
- Removable washable neck cover
- Control via cable remote control or via Bluetooth
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 120 x 26 x 46 cm
Weight: approx. 14 kg
Power supply: 220-240V~ 50/60 Hz
Power output: approx. 60 Watts
Artikel-No.: 88710
EAN: 40 15588 88710 9
Warranty: 3 Years

**Shiatsu massage
seat Cover**

MC 826



NEW



Q4/2017

Holistic relaxation with tapping and acupressure massage

- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Acupressure massage in the waist area with 3 selectable intensity levels
- 2 massage programs for the back selectable (shiatsu and tapping massage)

- Neck and shoulder massage (Height-adjustable)
- Separate vibration massage in the seat with 3 selectable intensity levels
- Targeted spot massage only for the tapping massage
- For tapping massage 3 selectable programs as well as 3 intensity levels selectable
- Automatic height-adjustable neck massage
- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Relaxed massage for the entire upper or lower back
- Removable washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 120 x 26 x 46 cm
Weight: approx. 14 kg
Power supply: 220-240V~ 50/60 Hz
Power output: approx. 60 Watts
Artikel-No.: 88926
EAN: 40 15588 88926 4
Warranty: 3 Years

**Shiatsu acupressure
massage seat cover**

MC 825



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Novel acupressure massage in the waist and thigh area
- Separate, height-adjustable neck massage

- 3 selectable massage zones: total, upper or lower back
- 3 intensity levels of acupressure massage selectable
- Intensive shiatsu massage in the back area
- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 110 x 42,5 x 18 cm
Weight: approx. 5,9 kg
Power supply: 220-240V~ 50/60 Hz
Power output: 40 Watts
Artikel-No.: 88939
EAN: 40 15588 88939 4
Warranty: 3 Years

**Shiatsu
massage seat cover**

MC 824



NEW



Q4/2017

The new feel-good massage - follows the individual back contour

- Pleasant massage along the spine with new rocking technique
- 3 massage types for the back area: Shiatsu, rolling and spot massage
- Shiatsu and rolling massage in 2 speeds

- 3 massage zones for shiatsu massage selectable: total, upper back or shoulder area
- Separate vibration massage in the seat
- Additional heat function
- Handy cable remote control
- Auto-Stop after 15 minutes of operation
- Removable, washable protective cover and pillowcases
- Easy to clean material
- Convenient velcro for fastening at the back of the chair
- 2 Demo-programs

Size: approx. 135 x 48 x 21 cm
Weight: approx. 7,8 kg
Power supply: 220-240V~ 50/60 Hz
Power output: approx. 60 Watts
Artikel-No.: 88921
EAN: 40 15588 88921 9
Warranty: 3 Years

**Shiatsu
massage seat cover**

MC 822



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- 3 types of massage in the back area: Shiatsu, rolling and spot massage
- Separate, height-adjustable neck massage

- 3 selectable massage zones: total, upper or lower back
- Separate vibration massage in the seat with 3 selectable intensity levels
- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 118 x 40 x 9 cm
Weight: approx. 6,5 kg
Power supply: 220-240V~ 50/60Hz
Power output: 40 Watts
Artikel-No.: 88922
EAN: 40 15588 88922 6
Warranty: 3 Years

**Comfort
massage seat cover**

MC 818



NEW



Q4/2017

Intensive relaxation with tapping massage

- Novel tapping massage in the entire back area with two intensity levels
- 3 selectable massage zones: entire, upper or lower back or point-specific spot massage

- Timer function: 5/10/15 minutes
- Additional red light and heat function
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 102 x 43 x 15 cm
Weight: approx. 3,4 kg
Power supply: 220-240V~ 50/60 Hz
Power output: 30 Watts
Artikel-No.: 88918
EAN: 40 15588 88918 9
Warranty: 3 Years

**Shiatsu
massage seat cover**

MC 815



NEW



Q4/2017

- High-definition massage: targeted, intensive massage for tangible sense of well-being
- 3 selectable massage zones: total, upper or lower back

- 4 rotating massage heads
- Additional red light and heat function
- Timer function: 5/10/15 minutes
- Separate vibration massage in the seat with 3 selectable intensity levels
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair
- Demo function

Size: approx. 110 x 47,2 x 11 cm
Weight: approx. 4,5 kg
Power supply: 220-240V~ 50/60Hz
Power output: 24 Watts
Artikel-No.: 88916
EAN: 40 15588 88916 5
Warranty: 3 Years

**Shiatsu massage
seat cover with gel cushions**

MC 830



Soft back massage

- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Massage heads made of high-quality gel for a softer and healthier massage
- Gentle to hips and spine and shoulder blades
- 3 selectable massage zones: total, upper or lower back

- Separate vibration massage in the seat with 3 selectable intensity levels
- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 105 x 43 x 8,5 cm
Weight: approx. 3,4 kg
Power supply: 100-240V~ 50/60Hz
Power output: 18 Watts
Artikel-No.: 88944
EAN: 40 15588 88944 8
Warranty: 3 Years



**Shiatsu
massage seat cover**

MCN PRO



Our bestseller now with height-adjustable neck massage

- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Separate adjustable neck massage
- Separate vibration massage in the seat with 3 selectable intensity levels

- 3 selectable massage zones: total, upper or lower back
- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 118 x 40 x 9 cm
Weight: approx. 6,5 kg
Power supply: 220-240V~ 50/60Hz
Power output: 48 Watts

Artikel-No.: 88970
EAN: 40 15588 88970 7

Warranty: 3 Years

**Shiatsu
massage seat cover**

MCN



Our bestseller

- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Separate neck massage
- Separate vibration massage in the seat with 3 selectable intensity levels

- High-definition massage: targeted, intensive massage for tangible sense of well-being
- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Removable washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 118 x 40 x 9 cm
Weight: approx. 6,5 kg
Power supply: 220-240V/50-60 Hz
Power output: 48 Watts

Artikel-No.: 88930
EAN: 40 15588 88930 1

Warranty: 3 Years

**Shiatsu
massage seat cover**

RBI



- High-definition massage: targeted, intensive massage for tangible sense of well-being
- 3 selectable massage zones: total, upper or lower back
- Separate vibration massage in the seat with 3 selectable intensity levels

- Additional red light and heat function
- Auto-Stop after 15 minutes of operation
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

Size: approx. 105 x 45 x 9 cm
Weight: approx. 4,5 kg
Power supply: 100-240V~ 50/60Hz
Power output: 18 Watts

Artikel-No.: 88914
EAN: 40 15588 88914 1

Warranty: 3 Years

**Vibration
massage seat cover**

MC 810



- 4 massage zones in the back and seat area
- 3 selectable massage intensities
- 9 selectable massage programmes
- Intensive heat function with or without massage possible

- 6 vibration motors
- Auto-Stop after 15 minutes of operation
- Handy cable remote control
- Easy to clean material
- Very handy and easy to store
- Also ideal for office use or in the car
- Incl. car power adapter

Size: approx. 107 x 46 x 4 cm
Weight: approx. 1,4 kg
Power supply: 100-240V~ 50/60Hz, 12 V DC 2500mA
Power output: 30 Watts

Artikel-No.: 88937
EAN: 40 15588 88937 0

Warranty: 3 Years

**Vibration
massage seat cover**

MCH



- 4 massage zones in the back and seat area
- 3 selectable massage intensities
- 5 selectable massage programmes
- Additional heat function

- 5 vibration motors
- Auto-Stop after 15 minutes of operation
- Handy cable remote control
- Easy to clean material
- Very handy and easy to store
- Also ideal for office use or in the car
- Incl. car power adapter

Size: approx. 107 x 47 x 3,5 cm
Weight: approx. 1,1 kg
Power supply: 100-240V~ 50/60Hz, 12V DC 800mA
Power output: 9,6 Watts

Artikel-No.: 88935
EAN: 40 15588 88935 6

Warranty: 3 Years

Presentation chair



- Ergonomic swivel chair for presentation of massage seat covers
- Optimum geometry for seat back massage cushion

- With attachment for remote control
- Suitable for the massage seat covers by Medisana
- Easy to assemble

Size: approx. 73 x 66 x 111 cm
Weight: approx. 13,3 kg
Max. load capacity: max. 120 kg

Artikel-No.: 88913
EAN: 40 15588 88913 4

Warranty: 3 Years

Neck massager

Neck massager

Neck massager

Shiatsu neck massager

Neck massager with tapping function

Shiatsu massage cushion

MNV

NM 870

NM 865

NM 860

MNT

MC 840



- 2-levelled vibration massage for the neck and shoulders
- Optional red light and warmth function
- Mains and battery powered – perfect for travelling too

- Neck massage - also ideal for travel**
- Neck Massager made of comfortable molded foam – fits perfectly to the contours of the head and neck
 - Pleasant vibration massage for the neck and shoulders
 - Ideal for travelling and at home

- 2-levelled vibration massage for the neck and shoulders
- Optional red light and warmth function
- Mains and battery powered – perfect for travelling too
- Comfortable use by means of practical Velcro fastening

- For an optimal neck massage**
- Intensive shiatsu massage in the neck area; also applicable on shoulder, back, abdomen, upper and lower leg
 - 4 rotating massage heads
 - Massage intensity individually adjustable through pulling the holding straps

- Vigorous tapping massage
- 5 selectable massage intensities
- 3 selectable massage programs
- Ergonomic shape

- Soothing Shiatsu massage for neck, shoulder, back and legs
- 4 rotating massage heads
- Additional red light and heat function
- Particularly flat design for flexible, convenient operation

- Simple to use and cordless thanks to the integrated operating panel

- With integrated „One Touch“- control panel
- Extra soft and fluffy fleece cover
- With integrated storage bag – Space-saving storage or on the go
- Including 2 batteries (1.5 V AA)

- Simple to use and cordless thanks to the integrated operating panel

- 2 levels of rotation speed selectable
- Additional light and heat function in the upper massage heads
- Simple to use and cordless thanks to the integrated operating panel
- Removable, washable cover over the massage heads

- Suitable for neck, shoulders, legs and hips
- 2 massage heads
- Position of the individual tapping elements are easy movable
- Handy remote control

- Handy cable remote control
- Convenient velcro for fastening at the back of the chair
- Easy to clean material

Size: approx. 34 x 24 x 13,5 cm
Weight: approx. 0,8 kg
Power supply: 100-240V~ 50/60Hz, 6V DC 600mA oder mit 4 x AA (1,5V) Batterien
Power output: 3,6 Watts
Artikel-No.: 88941
EAN: 40 15588 88941 7
Warranty: 3 Years

Size: approx. 28 x 26 x 9 cm
Weight: approx. 0,4 kg
Batteries: 2 x 1,5 V (AA) Batteries
Artikel-No.: 88946
EAN: 40 15588 88946 2
Warranty: 3 Years

Size: approx. 35,5 x 24 x 12 cm
Weight: approx. 0,4 kg
Power supply: 100-240V~ 50/60Hz, 6V DC 600mA or with 4 x AA (1,5V) Batteries
Power output: 3,6 Watts
Artikel-No.: 88945
EAN: 40 15588 88945 5
Warranty: 3 Years

Size: approx. 58 x 34 x 20 cm
Weight: approx. 1,55 kg
Power supply: 100-240V ~ 50/60Hz, 12V DC 1,500 mA
Power output: 18 Watts
Artikel-No.: 88942
EAN: 40 15588 88942 4
Warranty: 3 Years

Size: approx. 61 x 43 x 6 cm
Weight: approx. 2,6 kg
Power supply: 230 V, 50 Hz
Power output: 75 Watts
Artikel-No.: 88940
EAN: 40 15588 88940 0
Warranty: 3 Years

Size: approx. 31 x 23 x 11 cm
Weight: approx. 1,3 kg
Power supply: 100-240V ~ 50/60 Hz 12V DC 1000 mA
Power output: 12 Watts
Artikel-No.: 88949
EAN: 40 15588 88949 3
Warranty: 3 Years



Massage mat

MM 825



- 4 selectable massage zones: upper and lower back, upper and lower leg
- 2 intensity levels for each massage zone
- Comfortable vibration massage for the whole body
- Smooth heat function in the back area

- 5 selectable massage programmes
- Auto-Stop after 15 minutes of operation
- Handy cable remote control
- Extra soft and cosy fleece cover
- Cushion in the head area of the mat for lying comfortably

Size: approx. 180 x 62 x 10 cm
Weight: approx. 1,8 kg
Power supply: 100-240V~ 50/60Hz, 12V DC 1200mA
Power output: 14,4 Watts
Artikel-No.: 88955
EAN: 40 15588 88955 4
Warranty: 3 Years

Hand massager

HM 886



- Red light and heat function
- 2 selectable intensity levels
- Different massage angles adjustable
- Including 5 massage attachments

- Handy design for a convenient self-massage



Including 5 massage attachments

Size: approx. 41.1 x 10.2 x 12.2 cm
Weight: approx. 0,3 kg
Power supply: 220-240V
Power output: 6,5 Watts
Artikel-No.: 88268
EAN: 40 15588 88268 5
Warranty: 3 Years

Hand massager

HM 858



- Vigorous vibration massage
- Continuously adjustment of massage intensity
- Additional red light and heat function
- Various massage attachments: acupressure and spot massage

- Easy to use thanks to its ergonomic shape
- Non-slip grip
- Extra-long cable



Including 4 massage attachments

Size: approx. 32,5 x 17 x 14,3 cm
Weight: approx. 1,8 kg
Power supply: 220-240V~ 50/60Hz
Power output: 28 Watts
Artikel-No.: 88292
EAN: 40 15588 88292 0
Warranty: 3 Years

Mini-hand massager

HM 840



- Comfortable vibration massage
- With removable and washable massage pad made of natural loofah
- Waterproof for use in the shower or when bathing

- Intensive peeling effect when using with loofah pad
- Light and handy design for convenient selfmassage
- Including 3 batteries (1.5 V AAA)



Including loofah pad

Size: approx. 8 x 8 x 7 cm
Weight: approx. 0,2 kg
Batteries: 3 x 1,5 V (AAA)
Artikel-No.: 88521
EAN: 40 15588 88521 1
Warranty: 3 Years



INFO Hand massager HM 886

Absolutely precise release of tensions. The hand massager HM 886 will support this.

What is the best way to relieve tensions? The hand massage device HM 886 offers various massage attachments for the effective acupressure or point massage. You can increase the massage effect by switching on red light and warmth. The massage intensity is infinitely adjustable.

INFO Shiatsu foot massager FM 883

A blessing for the feet. A few minutes on the Shiatsu foot massage device FM 883, and the tension in the feet will be released.

The Shiatsu massage acts upon the foot reflex zones and releases tension in the foot muscles with 18 gentle rotating massage heads. Thereby, the invigorating effect travels all the way up to the back. Red light and heat can be switched on for an additional support.

Foot massager

FM 888



NEW



- Intensive foot reflexology massage
- 2 different selectable massage programmes
- Shiatsu massage of the foot soles and heels
- Upper foot compression massage
- Comfortable Size (suitable up to shoe size 46)

- Additional red light and heat function
- 2 intensity levels selectable for the compression massage
- Promotes the blood circulation and invigorates tired feet
- Removable, washable inner lining
- Easy to clean material

Size: approx. 41,7 x 35,7 x 22,4 cm
Weight: approx. 4,1 kg
Power supply: 220-240V~ 50/60Hz
Power output: 50 Watts

Artikel-No.: 88398
EAN: 40 15588 88398 9

Warranty: 3 Years

Shiatsu foot massager

FM 885



NEW



- Shiatsu massage of the foot soles and heels
- Compression massage of the foot topsides in 3 intensity levels
- Switchable heat in the Shiatsu massage heads

- Space-saving softcover
- Acts blood circulation promoting and invigorating
- Comfortable to use with the foot

Size: approx. 40,2 x 37,5 x 17,7 cm
Weight: approx. 3,5 kg
Power supply: 220-240V~ 50/60Hz
Power output: 28 Watts

Artikel-No.: 88397
EAN: 40 15588 88397 2

Warranty: 3 Years

Shiatsu foot massager

FM 883



NEW



- Soothing Shiatsu massage
- Suitable for feet and back
- 2 selectable massage intensities
- Additional red light and heat function

- 6 x 3 rotating massage heads
- Two height adjustable feet to adjust the massage position
- Handy cable rewind under the device
- Simple one-button operation

Size: approx. 40,5 x 36,5 x 12 cm
Weight: approx. 3,3 kg
Power supply: 230V~ 50/60Hz
Power output: 30 Watts

Artikel-No.: 88396
EAN: 40 15588 88396 5

Warranty: 3 Years

Foot spa

FS 885 COMFORT



- Relaxing spa and vibration massage
- 5 in 1: Bubble- and vibrating massage, red light and heat function and massage rollers in the middle of the foot bubble bath
- Water can be heated to the preferred temperature (20°C to 48°C)
- Timer function for 20/30/40/50/60 min.

- Massage rollers in the middle of the foot bubble bath to massage the inner side of the foot
- Digital display for timer setting and water temperature
- Red light and heat function
- Ergonomically arranged air bubble nozzles
- Removable massage rollers in the area of the soles
- Good heat insulation through double-walled casing
- Intensive reflexology
- Comfortable Size (suitable up to shoe size 46)

Size: approx. 38,6 x 24,3 x 45,9 cm
Weight: approx. 2,9 kg
Power supply: 220-240V~ 50Hz
Power output: 390 Watts
Artikel-No.: 88378
EAN: 40 15588 88378 1
Warranty: 3 Years

Foot spa

FS 883



- Relaxing spa and vibration massage
- Pleasant red light area
- 3 in 1: water heating, bubble- and vibrating massage
- Function for keeping the water warm
- 3 exchangeable pedicure attachments (Cornea remover, brush, massage attachment)

- Massage rollers for foot reflexology massage
- Massage-aiding footbed
- Dry massage function
- Non-slip rubber feet
- Removable splash protection
- Comfortable Size (suitable up to shoe size 46)

Size: approx. 41 x 35 x 18 cm
Weight: approx. 1,5 kg
Power supply: 220-240V~ 50/60Hz
Power output: 80 Watts
Artikel-No.: 88387
EAN: 40 15588 88387 3
Warranty: 3 Yearse

Foot spa

FS 881



- Relaxing spa and vibration massage
- Soothing foot spa
- 3 in 1: water heating, bubble- and vibrating massage
- Intensive reflexology

- I Includes massage attachment
- Non-slip rubber feet
- Automatic heat function
- Non-slip rubber feet
- Comfortable Size (suitable up to shoe size 44)

Size: approx. 39 x 33 x 13 cm
Weight: approx. 1,22 kg
Power supply: 220-240V~ 50Hz
Power output: 60 Watts
Artikel-No.: 88380
EAN: 40 15588 88380 4
Warranty: 3 Years



Bath spa

BBS



- Relaxing air bubble bath to stimulate blood circulation and to relieve muscles
- Massage/bubble function with 3 intensity levels and combination program
- Intensive heat function selectable
- Extra quiet

- Automatic heat function
- Intensity levels adjustable on the device or via remote control
- Timer function: 10/20/30 minutes
- Non-slip soft mat for lying comfort
- Suitable for any bath tub
- Integrated fragrance dispenser for scented oils
- Double safety valve for air hose and device
- Incl. 3 aroma dispensers
- Incl. remote control

Size: Soft mat: approx. 120 x 36 cm
 Air hose: approx. 2,4 m
 Motor: approx. 31 x 21 x 31 cm
 approx. 4,3 kg (total)
Weight:
Power supply: 220-240V~ 50Hz
 Remote control: 3 V CR2032
Power output: 570 Watts
Artikel-No.: 88386
EAN: 40 15588 88386 6
Warranty: 3 Years

Bath spa

MBH



- Relaxing air bubble bath to stimulate blood circulation and to relieve muscles
- Massage/bubble function with 2 intensity levels and combination program
- Auto-Stop after 20 minutes of operation

- Automatic heat function
- Intensity levels adjustable on the device or via remote control
- Non-slip soft mat for lying comfort
- Suitable for any bath tub
- Integrated fragrance dispenser for scented oils
- Double safety valve for air hose and device
- Incl. 3 aroma dispensers
- Incl. remote control

Size: Soft mat: approx. 120 x 36 cm
 Air hose: approx. 2,4 m
 Motor: approx. 25 x 16 x 17 cm
 approx. 2,7 kg (total)
Weight:
Power supply: 220-240V~ 50Hz
 Remote control: 3 V CR2032
Power output: 450 Watts
Artikel-No.: 88375
EAN: 40 15588 88375 0
Warranty: 3 Years

Heat pad

HP 605



- Ultra-soft breathable cushion
- Safety system with overheating protection
- Quick heating with 100 Watts
- 4 temperature settings

- Automatic switch-off after 90 minutes
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Removable cover of cotton
- Cover machine-washable at 30°C
- Cushion washable thanks to detachable controller

Size: approx. 40 x 33 cm
Weight: approx. 0,5 kg
Power supply: 220-240V~ 50Hz
Power output: 100 Watts

Artikel-No.: 61147
EAN: 40 15588 61147 6
Warranty: 3 Years

Heat pad for neck and back

HP 610



- Ultra-soft breathable cushion
- Safety system with overheating protection
- Quick heating with 100 Watts
- 4 temperature settings

- Automatic switch-off after 90 minutes
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Removable cover of cotton
- Cover machine-washable at 30°C
- Cushion washable thanks to detachable controller
- Perfect fit for the back and neck
- Convenient velcro strap for individual fastening of the cushion

Size: approx. 40 x 60 cm
Weight: approx. 0,5 kg
Power supply: 220-240V~ 50Hz
Power output: 100 Watts

Artikel-No.: 61167
EAN: 40 15588 61167 4
Warranty: 3 Years

Heat pad for neck and shoulder

HP 622



- Flexible heating wires and a matrix sensor system for even heat distribution
- 6 temperature settings
- Ergonomic shape for pleasant warmth at shoulder and neck

- Turbo heat with powerful 100 watts
- Automatic shut-off after 90 minutes
- Handy controller
- Machine washable at 40°C
- Oeko-Tex® standard 100

Size: approx. 56 x 52 cm
Weight: approx. 0,5 kg
Power supply: 220-240V~ 50Hz
Power output: 100 Watts

Artikel-No.: 61155
EAN: 40 15588 61155 1
Warranty: 3 Years

Heat cape with rechargeable batteries

HP 626



NEW



- Mobile warmth**
- Heat cape with rechargeable batteries for indoor and outdoor use
 - 4 temperature settings
 - 2500 mAh Li-Ion battery is enough for 300 minutes of operation at the highest level
 - Individually customizable fit

- Fast direct warming
- Charging power supply for a fast recharging (approx. 2-3 hours)
- Automatic switch-off after 90 minutes
- Overheating protection

Size: approx. 49 x 42 cm
Weight: approx. 0,2 kg
Capacity: 7,2V DC 2500mAh

Artikel-No.: 61165
EAN: 40 15588 61165 0
Warranty: 3 Years

Heat belt with rechargeable batteries

HS 680



- Mobile warmth**
- Heat belt with rechargeable batteries for indoor and outdoor use
 - Soothing heat for the back and abdomen
 - 4 temperature settings
 - Maximum freedom of movement – can also be worn under clothing

- Easy to attach thanks to Velcro system
- Additional elastic belt for stabilization
- Overheating protection
- Automatic switch-off after 90 minutes
- Ideal at home and when travelling
- Long-lasting heat thanks to the powerful lithium-ion battery – min. 3 hours at max. temperature
- Machine washable at 30°C

Size: approx. 125 x 20,5 cm
Weight: approx. 0,5 kg
Capacity: 7,2V DC 2500mAh

Artikel-No.: 61160
EAN: 40 15588 61160 5
Warranty: 3 Years



Heat pad

HP 625 COMFORT



Fastest comfortable warmth

- Ultra-soft breathable cushion
- Innovative 4D heat flow technology – for an effective heat development
- Flexible heating wire and matrix sensor system for an even heat distribution
- 4 temperature levels

- “APS” technology – active protection system: visual and acoustic temperature control
- Turbo heat with powerful 100 Watts
- Automatic switch-off after 90 minutes
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Cover machine-washable at 40°C
- Removable cover of cosy and soft upper material
- Oeko-Tex® standard 100

Size: approx. 46 x 35 cm
Weight: approx. 0,5 kg
Power supply: 230V~ 50Hz
Power output: 100 Watts

Artikel-No.: 61140
EAN: 40 15588 61140 7

Warranty: 3 Years

Heat pad for shoulder and back

HP 630



Fastest comfortable warmth

- Ultra-soft breathable cushion
- Innovative 4D heat flow technology – for an effective heat development
- Flexible heating wire and matrix sensor system for an even heat distribution
- 4 temperature levels

- “APS” technology – active protection system: visual and acoustic temperature control
- Turbo heat with powerful 100 Watts
- Automatic switch-off after 90 minutes
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Machine-washable at 40°C
- Ergonomic shape for pleasant warmth at shoulder and neck
- Cosy and soft upper material
- Oeko-Tex® standard 100

Size: approx. 55 x 65 cm
Weight: approx. 0,76 kg
Power supply: 230V~ 50Hz
Power output: 100 Watts

Artikel-No.: 61157
EAN: 40 15588 61157 5

Warranty: 3 Years

Cosy heating blanket XXL

HB 675



XXL Version

- 4 temperature settings
- Quick heating with 120 Watts
- Extra big: 200 x 150 cm
- Automatic switch-off after 180 minutes

- Safety system with overheating protection
- Handy controller
- Reversible blanket with 2-colour design
- Pleasantly soft upper material
- Washable thanks to detachable controller
- Oeko-Tex® standard 100

Size: approx. 200 x 150 cm
Weight: approx. 2,2 kg
Power supply: 220-240V~ 50Hz
Power output: 120 Watts

Artikel-No.: 60230
EAN: 40 15588 60230 6

Warranty: 3 Years

INFO Cosy heating blanket
 XXL HB 675

Heat XXL. The elegant and stylish cosy heating blanket is suitable for many occasions.

The cosy heating blanket XXL HB 675 offers ample space for the relaxing cuddling. Whether it's on a cold winter evening on the sofa, or on a cool summer night in the garden - it always delivers the right degree of warmth with four temperature levels.

Cosy heating blanket

HDW



- Quick heating with 120 Watts
- 4 temperature settings
- Automatic switch-off after 180 minutes

- Safety system with overheating protection
- Handy controller
- Reversible blanket with 2-colour design
- Pleasantly soft upper material
- Washable thanks to detachable controller
- Oeko-Tex® standard 100

Size: approx. 180 x 130 cm
Weight: approx. 1,7 kg
Power supply: 220-240V~ 50Hz
Power output: 120 Watts

Artikel-No.: 60227
EAN: 40 15588 60227 6

Warranty: 3 Years



Foot warmer

FWS



- Quick heating with 100 Watts
- 3 temperature settings
- Removable fluffy soft lining
- Large foot space (suitable up to shoe size 46)

- Safety system with overheating protection
- Extra-soft breathable core
- Automatic switch-off after 90 minutes
- Ergonomically designed controller with LED control display
- Lining machine washable at 30°C
- Oeko-Tex® standard 100

Size: approx. 30 x 30 x 23 cm
Weight: approx. 0,6 kg
Power supply: 220-240 V, 50/60 Hz
Power output: 100 Watts

Artikel-No.: 60257
EAN: 40 15588 60257 3
Warranty: 3 Years

INFO Foot warmer FWS

When temperatures start to fall, your feet are the first to get cold. The FWS foot warmers heat the feet up quickly and directly, thereby also boosting the circulation.

Especially women often have cold feet. This can have many causes. The focus is on the blood circulation. The FWS foot warmer can help with three adjustable temperature settings. It heats up quickly thanks to the 100 watt performance, and stays warm up to 90 minutes.



Heated underblanket

HU 662



- Extra warm foot zone
- Safety system with overheating protection
- 6 temperature settings
- Quick heating with 100 Watts

- 2 warming zones for body and feet
- Automatic switch-off after 180 minutes
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Washable thanks to detachable controller
- Pleasantly soft upper material
- Suitable for mattress sizes 200 x 90-100 cm
- Oeko-Tex® standard 100

Size: approx. 150 x 80 cm
Weight: approx. 1 kg
Power supply: 230V~ 50Hz
Power output: 100 Watts

Artikel-No.: 61220
EAN: 40 15588 61220 6
Warranty: 3 Years

Heated underblanket

HU 665



- Our bestseller**
- Safety system with overheating protection
- 3 temperature settings
- Washable thanks to detachable controller

- Pleasantly soft upper material
- Automatic switch-off after 180 minutes
- Handy controller
- Suitable for mattress sizes 200 x 90-100 cm
- Oeko-Tex® standard 100

Size: approx. 150 x 80 cm
Weight: approx. 1 kg
Power supply: 230V~ 50Hz

Power output: 60 Watts
Artikel-No.: 60217
EAN: 40 15588 60217 7

Warranty: 3 Years

Heated underblanket

HU 650



- Convenient heated underblanket with stretch function for easy attachment to the mattress
- 4 temperature settings
- Quick heating with 100 Watts
- Upper material of high-quality cotton

- Safety system with overheating protection
- 2 warming zones for body and feet
- Automatic switch-off after 180 minutes
- Ergonomically designed controller with LED control display
- Washable thanks to detachable controller
- Suitable for mattress sizes 200 x 90-100 cm

Size: approx. 200 x 90-100 cm
Weight: approx. 1,4 kg
Power supply: 230V~ 50Hz
Power output: 100 Watts

Artikel-No.: 61210
EAN: 40 15588 61210 7

Warranty: 3 Years

A good feeling, to be able to help yourself

Effective therapy products – Made for Life.

Paula, 6
Has a sensitive lung.

Bernd, 40
No longer feels tired
on dark days.

Many people know that they can do a lot for their health themselves. Now they are taking responsibility, and handling their health in a self-determined manner. Wellness and fitness are an essential part of this, and a lot can also be achieved in the therapeutic area. The aim is not to replace the physician. But one's own well-being can also be taken into the own hands. MEDISANA supports people to prudently monitor their own health, with therapeutic products which provide a gentle and specific assistance.

Light on for a good mood.


People shouldn't be depressed by gloomy weather or long winter days.

Energising daylight applied in a specific manner lifts the mood, and acts against the so-called "winter

blues". The Daylight Lamp LT 460 are able to effectively replenish the positive mood.

Daylight lamp LT 460

- Energising daylight
- Reduction of winter depression
- Uniform light distribution with LED technology
- 2 light intensities can be selected
- Certified medical product**



Small and flat like a tablet.



Lea, 33
Feels well.



Light lifts the mood, and provides energy as well as the joy of life.

Mia, 4
Needs light,
to be cheerful.

INFO Healthy light

What is healthy light?

The quality of light greatly influences the well-being, increases productivity and the quality of life. Healthy light is like sunlight. It has a red portion that is responsible for the melatonin secretion in the body. The hormone controls the day-night rhythm.

* Certified according to MDD (Medical Device Directive)

Peter, 39
Now enjoys getting up in the morning.

Awakened by the sun.

A comfortable way to start the day. The Light Alarm Clock WL 444 will put a smile on your face.

Sofie, 33
Likes it romantic.



A good mood is independent from the weather.

The proper light spectrum will determine whether or not we feel awake rested and powerful. It helps us to feel fit and healthy. The Light alarm clock WL 444 with the simulated sunrise will revive your spirits.




Wellness light with a colour change in 7 colours



Wake up light WL 444


- Relaxed awakening through a gentle start into the day
- Integrated FM Radio
- Durable, energy-saving LED lighting



For a better trimmer tummy.

The abdominal muscles are a "problem zone" for many people. With electrical muscle stimulation (EMS) you train your abdominal muscles in a gentle and lasting fashion.

Christian, 37
In best form.



NEW

Abdominal muscle stimulator AM 880

- Flexible neoprene waist-belt with inserted electrodes
- 6 abdominal muscle training programs, each 30 min long
- 40 intensity levels selectable
- Higher wearing comfort thanks to a flat design



In physiotherapy, the electrical muscle stimulation has long been used for the build-up and regeneration of saggy muscles. With the AM 880 Abdominal muscle stimulator, this effect can also be use at home. The front abdominal muscles are stimulated in a targeted manner, through light muscle contortions which promote the blood circulation.



Improves breathing instantly.

People with sensitive lungs and bronchia react particularly sensitive to the cold and humid season. Regular inhalations alleviate the symptoms.

For children and adults with asthma or bronchitis, the regular inhaling is part of everyday life during the seasons of the year. Thanks to a very fine nebulization technology and a high respirable aerosol content, the IN 500 compact is able to precisely apply its beneficial effects. The additional nose adapter gives children more freedom of movement during the inhalation process.

Inhaler IN 500 compact

- Targeted treatment of diseases of the upper and lower respiratory tract
- Fogging with compressed air technology
- Particularly effective through a high respirable aerosol content
- Comprehensive accessories: Inhaler with mask, mouthpiece, nose adapter, extra long tube
- Certified medical product**



Our BESTSELLER

Viktoria, 23
Can breathe again.



Tough against lice. Gentle to the scalp.

Lice usually originate from kindergarten or school. Unfortunately they do not go away by themselves. The Electrical Lice Comb gently gets rid of lice and nits.

Now with LED LIGHT

NEW

Electric lice comb LC 870


- Integrated LED light for a better discoverability and removal of scalp lice
- No side effects, since no chemical agents are used
- Also suitable for dogs and cats

If children in the kindergarten or school age suffer from itching, skin redness and sleep disorders, the diagnosis is often: Lice. In order to avoid that combing the hair from becomes a painful and unloved procedure, the Electric lice comb LC 870 is equipped with an LED light. It kills scalp lice and nits during the combing through electrical impulses, without the use of chemical agents.

The hair-friendly way, to combat lice.



* World Health Organisation ** Certified according to MDD (Medical Device Directive)



Daylight lamp

LT 480



- Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- UV Filter
- Flicker-free

- Light intensity: up to 10,000 Lux
- Comfortable one-button operation
- Particularly silent
- Folding base (wall mounting possible)
- Two 36 W lamps
- Certified medical device*

Size: approx. 48×34×10 cm
Weight: approx. 2 kg
Light area: 31×45 cm
Power supply: 220-240V~, 50 Hz
Power output: 72 Watts
Cable length: approx. 1.5 m
Article No.: 45224
EAN: 40 15588 45224 6
Warranty: 3 Years

Daylight lamp

LT 470



- Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- UV Filter
- Flicker-free

- Light intensity: up to 10,000 Lux
- Comfortable one-button operation
- Particularly silent
- Folding base (wall mounting possible)
- Energy-saving 45W lamp
- Certified medical device*

Size: approx. 27.3×20.7×7.4 cm
Weight: approx. 0.8 kg
Light area: 20×21 cm
Power supply: 220-240V~, 50 Hz
Power output: 45 Watts
Cable length: approx. 1.3 m
Article No.: 45222
EAN: 40 15588 45222 2
Warranty: 3 Years

Daylight lamp

LT 460



- Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- Uniform light distribution with LED technology
- 2 light intensities can be selected

- Light intensity: up to 10,000 Lux
- Comfortable one-button operation
- Can be used horizontally and vertically
- Ideal for the home and on the road
- LED technology
- Including storage bag
- Certified medical device*

Size: approx. 26×15×1.7 cm
Weight: approx. 0.5 kg
Light area: 20.5×12 cm
Power supply: Input: 100-240V~, 50/60 Hz, 500mA; Output: 12 V- 1 A
Power output: 12 Watts
Cable length: approx. 1.8 m
Article No.: 45220
EAN: 40 15588 45220 8
Warranty: 3 Years



Wake up light

WL 460



- Natural awakening for a soft start in the day
- Natural falling asleep through red light** (Timer can be set: 30, 60 or 90 min)
- Integrated radio
- Long-lasting, energy-saving LED lighting

- Increasing light intensity: 30 minutes before alarm
- 2 alarms can be adjusted
- Snooze function and zizz button (extra timer for alarm after 15 to 90 minutes)
- 8 natural sounds with increasing volume: tweeting birds, croaking frogs, running brook, rushing waves, crackling fireplace, falling rain, chirping crickets, Asian wooden gong
- Colour-changing wellness light with 7 colours
- Large, easy-to-read display with dimmer function (adjustable light + dark)
- MP3 player connection is possible
- Incl. antenna

Size: approx. 19×14.5×27.5 cm
Weight: approx. 0.9 kg
Power supply: 100-240 V/ 50-60 Hz, 12V DC 1A
Power output: 12 Watts
Cable length: approx. 2 m
Article No.: 45115
EAN: 40 15588 45115 7
Warranty: 3 Years

Wake up light

WL 444



- Natural awakening for a soft start in the day
- Alarm function through an integrated FM radio or 8 natural sounds
- Long-lasting, energy-saving LED lighting

- Increasing light intensity: 30 minutes before alarm
- 2 alarms can be adjusted
- Touch button: snooze function
- Backup battery
- 8 radio-memory slots
- Light intensity can be adjusted
- Colour-changing wellness light with 7 colours
- Blue LCD display with backlight
- Incl. antenna

Size: approx. 17.5×17.5×14.5 cm
Weight: approx. 0.37 kg
Power supply: Netzadapter Input: 100-240V~, 50/60Hz, 350mA Output: 12V = 1000mA
Article No.: 45110
EAN: 40 15588 45110 2
Warranty: 3 Years

* Certification according to MDD (Medical Device Directive)

** Red light stimulates the discharge of sleep hormones and makes it easier to fall asleep

<p>Inhaler IN 550 PRO</p>	<p>Inhaler IN 500 COMPACT</p>	<p>Ultrasonic inhaler IN 525</p>	<p>Ultrasonic Inhaler USC</p>	<p>Electric Lice Comb LC 870</p>	<p>Electric Lice Comb LC 860</p>
		 NEW		 NEW	 NEW
					
	<p style="text-align: right;">Q4/2017</p>	<p style="text-align: right;">Q1/2018</p>	<p style="text-align: right;">Q1/2018</p>	<p style="text-align: right;">Q1/2018</p>	<p style="text-align: right;">Q1/2018</p>
<p>3 Fogging attachments for a more targeted application</p> <ul style="list-style-type: none"> • Nebulisation with compressed air technology • Including 3 attachments to produce different particle sizes <ul style="list-style-type: none"> - 3,0 – 4,0 µm for the deeper areas of the lung - 3,5 – 4,5 µm for respiratory tract and lungs front - 4,5 – 6,0 µm for mouth and throat • For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold) 	<p>Our bestseller</p> <ul style="list-style-type: none"> • Nebulisation with compressed air technology • For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold) 	<p>Small inhaler – ideal for when you're on the go</p> <ul style="list-style-type: none"> • High efficiency through innovative micro-membrane nebulisation (mesh technology) with ultrafine droplets • For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold) • Comprehensive set of accessories • Use for a short time in the supine position is also possible 	<ul style="list-style-type: none"> • High efficiency through innovative micro-membrane nebulisation (mesh technology) with ultrafine droplets • For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold) • Comprehensive set of accessories 	<p>Now with LED light</p> <ul style="list-style-type: none"> • Simple combing frees hair from head lice and nits • No side effects, and no use of chemical agents • Integrated LED light for an easier detection and removal of lice 	<ul style="list-style-type: none"> • Simple combing frees hair from head lice and nits • No side effects, and no use of chemical agents
<ul style="list-style-type: none"> • Comprehensive set of accessories: Inhalation with mask, mouthpiece or nose adapter – ideal for children • Highly efficient due to very precise nebulisation • Particularly effective thanks to high proportion of respirable aerosols • Extra long hose – 2 m • Integrated accessory tray • Certified medical device* 	<ul style="list-style-type: none"> • Comprehensive set of accessories: Inhalation with mask, mouthpiece or nose adapter – ideal for children • Highly efficient due to very precise nebulisation • Particularly effective thanks to high proportion of respirable aerosols • Extra long hose – 2 m • Certified medical device* 	<ul style="list-style-type: none"> • Particularly effective thanks to high proportion of respirable aerosols • Inhalation with mask or mouthpiece – also eminently suitable for children • Low-noise nebulisation • Low power consumption (no more than 0.9 Watts) • Automatic switch-off with an empty solution tank • Incl. storage case • Incl. 2 batteries (1.5 V AA) • Certified medical device* 	<ul style="list-style-type: none"> • Particularly effective thanks to high proportion of respirable aerosols • Inhalation with mask or mouthpiece – also eminently suitable for children • Low-noise nebulisation • Incl. storage case • Incl. 4 batteries (1.5 V AA) • Certified medical device* 	<ul style="list-style-type: none"> • The lice are destroyed by an electrical current upon contact with the teeth of the comb • Also suitable for dogs and cats • Incl. 3 batteries (1.5 V AAA) 	<ul style="list-style-type: none"> • The lice are destroyed by an electrical current upon contact with the teeth of the comb • Also suitable for dogs and cats • Incl. batterie (1.5 V AA)
 <p>Incl. mouthpiece, masks and nose adapters for children and adults, connecting piece, 3 attachments to produce different particle sizes, 5 replacement filters, 2 replacement fuses</p>	 <p>Incl. mouthpiece, masks and nose adapters for children and adults, connecting piece, 5 replacement filters</p>	 <p>Incl. mouthpiece and masks for children and adults</p>	 <p>Incl. mouthpiece, masks for children and adults and connecting piece</p>		
<p>Size: approx. 33×20×13.8 cm Weight: approx. 2 kg Power supply: 230V~/50 Hz Power output: 60 Watts Filling volume: max. 10 ml Nebulizing output: ≤ 0,2 ml/min Particle size: (3.0-4.0/3.5-4.5/4.5-6.0 µm) Hose length: approx. 2 m Article No.: 54530 EAN: 40 15588 54530 6 Warranty: 3 Years</p>	<p>Size: approx. 19.5×13.5×19.2 cm Weight: approx. 1.3 kg Power supply: 230V~/50 Hz Power output: 60 Watts Filling volume: max. 10 ml Nebulizing output: ≤ 0,2 ml/min Particle size: >60% ≤ 5 µm Hose length: approx. 2 m Article No.: 54520 EAN: 40 15588 54520 7 Warranty: 3 Years</p>	<p>Size: approx. 7.2 × 4.0 × 10.8 cm Weight: approx. 0.1 kg (without batt.) Batteries: 2×1.5 V (AA) Power output: 0,9 Watts Filling volume: max. 8 ml Nebulizing output: ≤ 0,2 ml/min Particle size: Average <5 µm Article No.: 54115 EAN: 40 15588 54115 5 Warranty: 3 Years</p>	<p>Size: approx. 14.2×7.5×7.5 cm Weight: approx. 0.3 kg (without batt.) Batteries: 4×1.5 V (AA) Filling volume: max. 8 ml Nebulizing output: ≤ 0,4 ml/min Particle size: Average <4,8 µm Article No.: 54105 EAN: 40 15588 54105 6 Warranty: 3 Years</p>	<p>Size: approx. 20.9×6.6×3.1 cm Weight: approx. 0.1 kg Batteries: 3×1.5V (AAA) Article No.: 41019 EAN: 40 15588 41019 2 Warranty: 3 Years</p>	<p>Size: approx. 19.3×6.6×2.9 cm Weight: ca. 0.095 kg Battery: 1×1.5V (AA) Article No.: 41017 EAN: 40 15588 41017 8 Warranty: 3 Years</p>

* Certification according to MDD (Medical Device Directive)

Abdominal muscle stimulator

AM 880



Q4/2017

- Flexible Neoprene waist belt with incorporated electrodes
- Choice of 6 abdominal muscle training programmes of 30 minutes each
- Intensity can be set over 50 levels
- Can be worn under clothing thanks to the flat design

- Automatic switch-off at the end of the programme
- Simple operation with a large LCD display
- Device is integrated into the belt
- 1 channel with 4 electrodes
- Waist belt with Velcro closer
- Incl. extension belt for waist circumference up to 140 cm
- Water contact electrodes, no need for contact gel
- Incl. 3 batteries (1.5 V AAA)

Size: approx. 110×16×3 cm
Weight: approx. 0.37 kg
Batteries: 3×1.5 V (AAA)
Channels: 1
Curve form: Bi-phase square wave pulse
Intensity: 0–70mA
Frequency: 30–90 Hz
Pulse width: 200–300 µs
Article No.: 88322
EAN: 40 15588 88322 4
Warranty: 3 Years

Abdominal muscle stimulator

BOB



- Flexible abdominal belt with inbuilt electrodes
- 10 pre-installed programmes
- Intensity can be set over 25 levels
- Can be worn under clothing thanks to the flat design

- Automatic switch-off at the end of the programme
- Simple operation with a large LCD display
- Device is integrated into the belt
- 1 channel with 2 electrodes
- Waist belt with Velcro closer
- Incl. 3 batteries (1.5 V AAA)
- Certified medical device*

Size: approx. 110×14.4×2.2 cm
Weight: approx. 0.23 kg
Batteries: 3×1.5 V (AAA)
Channels: 1
Curve form: Bi-phase square wave pulse
Intensity: 0–90mA
Frequency: 3–100 Hz
Pulse width: 150–200 µs
Article No.: 88320
EAN: 40 15588 88320 0
Warranty: 3 Years

Electronic body toner

BT 850



- Stimulates blood flow in the area of use
- Relaxes tense muscles
- Helps to build up muscles
- Compatible with: Bluetooth Smart 4.0 from iPhone 4s and Android 4.3 onwards

- 5 different application programmes for the waist, shoulder and neck, back, limbs
- 7 intensity levels to choose from
- Different user accounts can be saved with individual settings (name, age, height, weight, gender)
- With timer function
- Including 4 self-adhesive electrodes, storage box and replacement battery
- Only works in conjunction with the app
- Incl. storage box
- Incl. 2 batteries (CR2032)



Incl. 4 self-adhesive electrodes

Size: Box: 12×6.5×2.5 cm
 Electrodes: approx. ø 40mm
Weight: approx. 25 g (without storage box)
Batteries: 1×3 V (CR2032)
Cable length: approx. 40 cm
Article No.: 88330
EAN: 40 15588 88330 9
Warranty: 3 Years

Infrared lamp

IR 885



- Intense infrared light to relax muscles or treat colds
- Stimulates blood circulation
- Relieves pain
- Durable lamp (halogen)

- Adjustable tilt angle (35°)
- On/Off button
- Timer (15 to 1 minute)
- Fan caster for cooling
- Cable rewind
- Carrying handle
- Large-surface application
- Incl. lamp (halogen)
- Certified medical device*

Size: approx. 19.8×18.6×29.5 cm
Weight: approx. 1.5 kg
Power supply: 230 V/50 Hz
Power output: 300 Watts
Article No.: 88257
EAN: 40 15588 88257 9
Warranty: 3 Years

Infrared lamp

IRL



- Our bestseller**
- Intense infrared light to relax muscles or treat colds
 - Stimulates blood circulation
 - Relieves pain
 - Durable lamp (Philips)

- Effective heat therapy through pressed glass bulb
- Adjustable tilt angle
- Made in EU
- Certified medical device*

Size: approx. 19×14×22 cm
Weight: approx. 0.65 kg
Power supply: 230 V/50 Hz
Power output: 150 Watts
Article No.: 88254
EAN: 40 15588 88254 8
Warranty: 3 Years

Infrared lamp

IRH



- Intense infrared light to relax muscles or treat colds
- Stimulates blood circulation
- Relieves pain
- Durable lamp (Philips)

- Adjustable tilt angle
- Made in EU
- Certified medical device*

Size: approx. 16.5×13.5×20.5 cm
Weight: approx. 0.36 kg
Power supply: 230-240 V/50-60 Hz
Power output: 100 Watts
Article No.: 88258
EAN: 40 15588 88258 6
Warranty: 3 Years

* Certification according to MDD (Medical Device Directive)

Enjoying a healthy life at home.

Soothing ambient air – Made for Life.

Sabine, 38
Enjoys the freshness.

Mia, 9
Enjoys the day asthma-free.

Elias, 41
Feels fit for everyday life.

Linus, 7
Always sleeps with the light on.

In an ever-increasing hectic world, the home is becoming an important place of peace and wellbeing. This is where new life energy is replenished. Therefore, many people place a great emphasis topic of health. They design their home according to health-conscious criteria, and recognize that this can also be achieved with “small” sustainable measures. MEDISANA creates an atmosphere of easy breathing and relaxing, through effective and soothing indoor air innovations.

NEW Aroma diffuser AD 625

- Decorative aroma diffuser for scented oils made of bamboo
- Generates a fine mist using ultrasound technology, through which the scent is distributed through the air throughout the room
- Wellness light with a colour change in 6 colours



Made of real bamboo



Good air, good mood

Many studies have shown: Fragrances influence our well-being. The nose will clearly sense it when relaxation and peace permeate the air.

Nathalie, 24
This simply feels good.



Depending on the mood and situation, the aroma diffuser AD 625 made of bamboo can be used to distribute stimulating or relaxing fragrances throughout the room. The ultrasonic technology distributes the aromatized air very quietly in the room in the form of a fine mist. Aroma oils as well as the wellness light in six changing colours create a pleasant "wellness" atmosphere.



Good room air.
Good conversation climate.



Healthy room air, healthy living comfort.



Annette, 38
Sometimes likes to stay at home.

SANITIZING
PREHEAT
FUNCTION

NEW Air humidifier AH 661

- Air humidifier with ultrasonic technology to improve the room air
- Hygienically clean thanks to an integrated pre-heating function to reduce germs and bacteria
- Timer for programming the operating time: 2/4/6 hours



Good, healthy air has a humidity of 40 - 60%.

People inhale about 20,000 litres of air every day. It is good if this air is as clean as possible. This is ensured by the integrated pre-heating function in the AH 661 air humidifier. It significantly reduces the amount of germs and bacteria in the discharged air. For this purpose, the air humidifier scents the air with relaxing aromas – and if desired up to 14 hours.

There's something in the air.

The fragrance makes the difference.

Stimulating, calming or activating scents. There are many fragrance oils with very specific scents. The The air humidifier AH 662 conveys the moisture into the room unobtrusively through a micro fine atomisation.

Air humidifier AH 662

- Air humidifier with ultrasound technology to improve the breathable air
- Highly effective through a micro fine atomisation
- Aroma compartment for fragrance oils



AROME
DISPENSER
FOR FRAGRANCE OILS



Tim, 29
Has a cleared head again.

Air humidifier

AH 661



NEW



Q4/2017

Germ-reducing pre-heating

- Air humidifier with ultrasound technology to improve the breathing air quality
- Hygienically clean through an integrated pre-heating function to reduce germs and bacteria
- Pleasant fragrance with aroma-compartment – suitable for aromatic oils
- Timer for programming the operating time: 2/4/6 hours

- 3 Fogging intensities selectable
- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Water tank 3,500 ml
- Touch Control Buttons
- Low-noise and energy-saving

Size: approx. 18.5×18.5×28.1 cm
Weight: approx. 1.1 kg
Power supply: 200-240V/50Hz
Power output: 75 Watt
Nebulizing output: cold 250 ml/h, warm 350 ml/h
Operating time: approx. 10-14 hrs
Article No.: 60052
EAN: 40 15588 60052 4
Warranty: 3 Years

Air humidifier

AH 662



- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Integrated fragrance dispenser for fragrant oils

- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Water level indicator (full=blue, empty=red)
- Water tank 900 ml
- Incl. cleaning brush

Size: approx. 15×15×22 cm
Weight: approx. 0.45 kg
Power supply: 100-240V~, 50/60Hz, 1.6A
Power output: 30 Watts
Nebulizing output: max. 60 ml/h
Operating time: approx. 15 hrs
Article No.: 60077
EAN: 40 15588 60077 7
Warranty: 3 Years

Air humidifier

AH 660



- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Steplessly adjustable misting levels

- Removable tank
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving: 30 watts
- Incl. water filter
- Water tank 4,500 ml
- For rooms up to approx. 30 m²



Incl. water filter

Size: approx. 24×17×36 cm
Weight: approx. 1.3 kg
Power supply: 220-240 V/50-60 Hz
Power output: 30 Watts
Nebulizing output: max. 300 ml/h
Operating time: approx. 15 hrs
Article No.: 60067
EAN: 40 15588 60067 8
Warranty: 3 Years

Air humidifier

UHW



Our bestseller

- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Steplessly adjustable misting levels

- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Water tank 4,200 ml
- For rooms up to approx. 30 m²

Size: approx. 26,5×33×20 cm
Weight: approx. 1.46 kg
Power supply: 220-240 V/50-60 Hz
Power output: 30 Watts
Nebulizing output: max. 300 ml/h
Operating time: variable, depending on setting
Article No.: 60065
EAN: 40 15588 60065 4
Warranty: 3 Years

Air humidifier with timer function

ULTRABREEZE



- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Timer for programming the misting start and duration
- Pre-heating to kill germs and bacteria

- Switchable lights
- 2 intensity levels to choose from
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Integrated clock with alarm function
- Water tank 1,300 ml

Size: approx. 19.6×13.5×25.6 cm
Weight: approx. 1.1 kg
Power supply: 100 – 240V/50 Hz/1.5 A, 24 V DC 1.5A
Power output: 36 Watts
Nebulizing output: max. 65 ml/h
Operating time: Level 1: approx. 20 hrs
 Level 2: approx. 30 hrs
Article No.: 60050
EAN: 40 15588 60050 0
Warranty: 3 Years

Air humidifier

MEDIBREEZE



- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- 2 intensity levels to choose from

- Switchable lights
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Water tank 680 ml

Size: approx. 16×11×20 cm
Weight: approx. 1 kg
Power supply: 100-240 V/50-60 Hz/0.5 A, 13.5V DC 1A
Power output: 14 Watts
Nebulizing output: max. 80 ml/h
Operating time: Level 1: approx. 9 hrs
 Level 2: approx. 6 hrs
Article No.: 60003
EAN: 40 15588 60003 6
Warranty: 3 Years

Hygrometer

HG 100



- Displays humidity and room temperature for a healthy climate indoors
- Saves the highest and lowest values
- For hanging on the wall (magnetic mount) or standing upright

- Temperature measurement range: -10°C ~ +50°C (+14° ~ +122 F°)
- Humidity measurement range: 20% ~ 90%
- Comfort zone indicator (smiley display)
- Time display
- Alarm function
- Incl. battery (1.5 V AAA)

Size: approx. 8.2×7×2 cm
Weight: approx. 0.07 kg
Battery: 1×1.5 V (AAA)
Article No.: 60079
EAN: 40 15588 60079 1
Warranty: 3 Years



Room air purifier

MEDISANA AIR®



- Cleans room air using an antibacterial particulate filter (F9>PM2.5) from 3M® and additional activated carbon filtering; allergy-friendly
- Frees the air of particulates, dust, pollens, odours and bacteria
- Particularly effective and quiet with low power consumption

- Simple and safe to use thanks to the LED filter change indicator
- Switches automatically between the day and night modes (even quieter)
- Operated using the included power supply unit and USB cable
- For rooms up to approx. 60 m²
- Incl. USB cable
- Digital motor with an extra long operating life – 5 year warranty



Incl. 3M® filter

Size: approx. 14.8×14.8×30.7 cm
Weight: approx. 0.6 kg
Power supply: DC 5V/300mA
Power output: 1,5 Watts
Purifying output: 60m²/24h
Cable length: approx. 1 m
Article No.: 60300
EAN: 40 15588 60300 6
Warranty: 5 Years

Aroma diffusor

AD 640



- Aroma diffusor for fragrant oils
- Using ultrasound technology, creates a fine mist of fragrance that spreads throughout your entire room
- Wellness light with a colour change in 3 colours
- Timer for programming operation: 30/60/120/180 minutes

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Touch control buttons
- LED time display
- Incl. measure cup

Size: approx. 11.6×11.2×16.2 cm
Weight: approx. 0.43 kg
Power supply: 100-240V/50-60Hz, DC 24V/12W
Water tank: 100 ml
Nebulizing output: max. 15 ml/h
Operating time: approx. 7 hrs
Article No.: 60086
EAN: 40 15588 60086 9
Warranty: 3 Years

Aroma diffusor

AD 625



NEW



Q4/2017

- Our stylish new addition – made of real bamboo**
- Decorative bamboo aroma diffusor for scented oils
 - Using ultrasound technology, creates a fine mist of fragrance that spreads throughout your entire room
 - Wellness light with a colour change in 6 colours

- Stylish living accessory without visible control elements
- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving

Size: approx. ø16,6×7,2 cm
Weight: approx. 0,68 kg
Power supply: Mains adaptor
Water tank: 100 ml
Nebulizing output: max. 16 ml/h
Operating time: approx. 6 hrs
Article No.: 60083
EAN: 40 15588 60083 8
Warranty: 3 Years

Aroma diffusor

AD 620



- Our bestseller**
- Aroma diffusor for fragrant oils
 - Using ultrasound technology, creates a fine mist of fragrance that spreads throughout your entire room
 - Wellness light with a colour change in 5 colours
 - Timer for programming operation: 2 and 4 hours

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Incl. measure cup

Size: approx. 13.2×13.2×14.5 cm
Weight: approx. 0.25 kg
Power supply: 100-240V/50-60Hz, DC 24V/500mA; 12W
Water tank: 300 ml
Nebulizing output: max. 30 ml/h
Operating time: approx. 10 hrs
Article No.: 60082
EAN: 40 15588 60082 1
Warranty: 3 Years

Enjoy the feeling of doing something good for yourself.

Smart care products – Made for Life.

Anna, 32
Feels one with herself

Marie, 9
Has beautiful skin

The satisfaction with life greatly depends on someone's own body self-awareness. If I feel good in my skin, I feel even more beautiful. Every second woman consciously takes care of her body and grooms herself intensely. However, the desire to do something good for the own body and to take care of yourself, has also been an important issue for men for quite some time. Body care is an integral part of a natural lifestyle. MEDISANA offers care products that meet all requirements for care and well-being.



Quartz cartridge with 200,000 light pulses

NEW Hair removal system IPL 850

- Durable quartz cartridge with 200,000 light pulses
- Effective for a permanent hair reduction
- Kills the hair root deep under the skin
- Automatic skin sensor in 6 light intensity levels
- Integrated UV filter to protect the skin



So simple, and safe to smooth skin.



The IPL 850 provides a very gentle way to ensure that women have a beautiful and well-groomed skin throughout the year.

Hair removal just like in the cosmetic studio. The IPL 850 enables the professional hair removal at home. Dermatologically compatible, it gently and quickly removes the hair even on small body surfaces. With its skin type sensor, the system measures skin pigmentation before every light pulse, and thus reduces the hair in a very gentle manner.

The pleasant feeling of velvety smooth skin.

This is how gentle and lasting a permanent hair removal can be. The IPL 840 noticeably reduces hair growth visibly after just a few treatments, and above all without pain.

Celine, 29
Always an eye-catcher.



NEW Hair removal system IPL 840

- Effective for a permanent hair reduction
- Disables the hair root deep under the skin
- Automatic skin tone sensor in 6 light intensity levels
- Durable glass cartridge with 100,000 light pulses
- Integrated UV filter to protect the skin



Finest alignment for smooth skin: the automatic skin tone sensor and the 6 different light intensity levels.

Velvety smooth skin - what every woman wants. Shaving, growing or epilating don't all equally agree with each other. The IPL 840 hair removal system works dermatologically compatible with the IPL technology (Intense Pulse Light). The light pulses deactivate the root of the hair deep under the skin. IPL is also used by skin doctors for the hair removal because of its good compatibility. After only a few weeks the hair growth is reduced by 50%.

Well-groomed nails for right-and lefthanders.

Well-groomed hands and feet enhance the personal well-being and prevent diseases.

Attractive hands and feet are beautiful at any age, and in any season. But care requires time and the proper tools. The MP 815 gently and very thoroughly takes care of the nails, cuticle and cornea. With seven high-quality attachments for the different applications, two speed

levels, and the left and right rotation, the MP 815 can be adapted to the individual care needs very easily. This prevent a proper treatment of a discoloration and splintering of the nails.



NEW

Manicure and pedicure device MP 815

- Care and treatment of nails, cuticles and the cornea
- 7 high-quality accessory parts for different applications
- Easy operation, and an easy exchange of the attachments
- Left/right operation, and 2 speed levels



Easy to change, quick to use:
7 attachments for every application purpose.



Lara, 29
Loves her hands.

NEW

UV nail hardener ND 870

- UV nail hardener for a professional nagel-design
- Simple operation with a 4-stage timer (90 sec / 120 sec / 30 min / infinity)
- Removable base plate for exchanging the florescent tubes



More individual design options for the nails.

Individual nail design is a major trend. With the nail hardener ND 870 you can design nails in your own home.

Healthy and well-groomed nails are naturally beautiful. They can be even more individualized with their own naildesign. The ND 870 is ideally suited for the artistic nail design with UV gel and liquid acrylic.



Cellulite massager AC 850

- Effective massage for the stimulation of connective tissue
- Promotes the blood circulation
- 3 x 2 rotating massage rollers (exchangeable)
- 2 massage intensities selectable
- Particularly versatile thanks to the mains operation



Intensive
**ROLLER
MASSAGE**

Sophie, 25
Likes to show
her legs.



Noticeably better looking skin.

The regular massage with the AC 850 will ensure a noticeably tighter skin.

Cellulite is caused by the uneven distribution of body fat in the subcutaneous tissue, condensed connective tissue and a lack of blood circulation. The Cellulite Massager AC 850 promotes the blood circulation and completely stimulate the skin on legs and buttocks with three times two massage rollers.

Gentle on the skin. Tough on cellulite.

A massage with the AC 850 is a thoroughly stimulating experience.

The Cellulite Massager doesn't just stimulate the skin one-sided, but acts vitalizing through suction, vibration, and supporting red light. This makes the massage much more effective. The blood circulation is stimulated and the connective tissue is tightened.



Cellulite massager AC 855


- Effective vacuum & roll massage for the stimulation of connective tissue
- Red light the blood circulation for a noticeably tighter skin
- Infinitely adjustable suction performance



With
**VACUUM-
MASSAGE**

Suction, vibration, red light - all for the better circulation of the connective tissue.

Hannah, 27
Everything is perfect.



Deep
PURITY



Face cleansing brush FB 885

- Face cleaning brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products
- Rechargeable lithium-ion battery

Including 4 attachments
(Sensitive Skin, Normal Skin,
Peeling, Make-Up)-




Facial cleansing is the best foundation
for a proper care of the face.

Clean from the very start.

Our face is exposed to many stimuli and environmental challenges throughout the day. It is therefore beneficial if the facial skin is cleansed regularly.

Sun, wind, rain, dust. Our facial skin is exposed to many forms of stress every day. With the facial cleaning brush FB 885, the face can be thoroughly cleaned within 60 seconds. The gentle vibration removes dead skin cells, stimulates the circulation and prepares the skin perfectly for care products.



Relaxation for the face.

The wellness program for your face.
The facial sauna FSS will ensure healthy beauty.



Facial sauna FSS

- For cosmetic use
- Diffuser for an even steam distribution
- With aroma attachment for soothing applications
- Additional attachment for the nose



Healthy facial skin is naturally beautiful.

A beautiful, healthy skin complexion appears well-maintained and exudes attractiveness and self-confidence. The Facial sauna FSS relaxes the skin, opens and cleans the pores. The warm moisture thereby acts like a gentle spa. The skin feels healthy and refreshed.

Hair removal system

IPL 850



High-quality quartz cartridge with 200,000 flashes

- Effective technology for permanent hair reduction
- Suitable for the face, arms, legs, underarms, bikini zone, back, chest
- Automatic skin tone sensor in 6 light intensity levels (4 manually selectable levels)
- LCD display shows the adjusted energy intensity as well as the number of used light pulses

- Long-life quartz cartridge with 200,000 light pulses
- Disables the hair root deep under the skin and prevents further hair growth
- Easy and safe to use at home
- Safe and effective for the treatment of skin types I – VII and black to red hair
- Integrated UV filter
- Flash-glide function for comfortable use
- Sensor starts treatment when skin contact is made
- Automatic shut-off after 3 min.
- Spare cartridge:
Article No.: 88589 / EAN-No.: 40 15588 88589 1

Size: approx. 13,8 x 8 x 4,5 cm
Weight: approx. 280 g
Lighting area: 2,5 cm²
Intensity: 5 J/cm²
Wavelength: 470 – 1.200 nm
Power supply: 100–240V, 50/60Hz, 1,0A
Article No.: 88588
EAN: 40 15588 88588 4
Warranty: 3 Years

Hair removal system

IPL 840



- Effective technology for permanent hair reduction
- Suitable for the face, arms, legs, underarms, bikini zone, back, chest
- Automatic skin tone sensor in 6 light intensity levels (4 manually selectable levels)
- Durable glass cartridge with 100,000 light pulses

- Disables the hair root deep under the skin and prevents further hair growth
- Easy and safe to use at home
- Safe and effective for the treatment of skin types I – VII and black to red hair
- Integrated UV filter
- Flash-glide function for comfortable use
- Sensor starts treatment when skin contact is made
- Automatic shut-off after 3 min.
- Spare cartridge:
Article No.: 88586 / EAN-No.: 40 15588 88586 0

Size: approx. 13,8 x 8 x 4,5 cm
Weight: approx. 280 g
Lighting area: 2,7 cm²
Intensity: 5 J/cm²
Wavelength: 470 – 1.200 nm
Power supply: 100–240V, 50/60Hz, 1,0A
Article No.: 88585
EAN: 40 15588 88585 3
Warranty: 3 Years



Hair removal system

SILHOUETTE IPL 800



- Effective technology for permanent hair reduction
- Suitable for the face, arms, legs, underarms, bikini zone, back, chest
- Incl. 2 cartridges with 50,000 light pulses each

- Disables the hair root deep under the skin and prevents further hair growth
- Easy and safe to use at home
- 5 light intensity levels
- Safe and effective for the treatment of skin types I – IV and black to light brown hair
- Integrated UV filter
- Flash-glide function for comfortable use
- Sensor starts treatment when skin contact is made
- Automatic shut-off after 10 min.
- Mains adapter included



Replacement cartridge with 50.000 light pulses incl.

Size: approx. 14 x 8,1 x 4,8 cm
Weight: approx. 0,2 kg
Lighting area: 2,7 cm²
Intensity: 4 J/cm²
Wavelength: 475-1.200 nm
Power supply: 100-240V, 50/60Hz, 0,8A-DC 12V, 3A
Article No.: 88580
EAN: 40 15588 88580 8
Warranty: 3 Years

Cellulite massager

AC 855



Vacuum massage

- Effective vacuum and rolling massage for stimulating the tissue for tangibly tauter skin
- Infrared heating function stimulates the circulation
- Additional vibration enhances the massage effect
- Stimulates the lymph circulation and strengthens the connective tissue

- Supports the treatment of cellulite
- The massaging intensity is continuously variable
- Light and handy design for a comfortable self-massage
- Single-hand operation
- Especially powerful thanks to mains operation
- Removable handle

Size: approx. 12,6 x 10,6 x 8,2 cm
Weight: approx. 0,26 kg
Power supply: 110-240 V, 50/60 Hz; 9 V DC 600 mA
Power output: 5,4 Watts
Cable length: 1,80 m
Article No.: 88545
EAN: 40 15588 88545 7
Warranty: 3 Years

Cellulite massager

AC 850



Intensive rolling massage

- Effective massage to stimulate the connective tissue for noticeably firmer skin
- Can help improve blood circulation
- 2 selectable massage intensities
- Light and handy design for convenient self-massage

- 3 x 2 rotating massage rollers (removable)
- Individually adjustable hand strap
- Particularly powerful thanks to mains operation
- Easy cleaning of the removable massage rollers

Size: approx. 13 x 11 x 8,5 cm
Weight: approx. 0,35 kg
Power supply: 110-240 V, 50/60 Hz; 12 V DC 1200 mA
Power output: 14,4 Watts
Cable length: 2 m
Article No.: 88540
EAN: 40 15588 88540 2
Warranty: 3 Years

Manicure and pedicure device

MP 840



- Care and treatment of nails, cuticles and slight calluses
- With LED light for accurate application
- 3 speed levels for the individual requirements of care-giving

- Clockwise and counter-clockwise rotation
- High quality and durable attachments
- Stops automatically when the pressure is too high – also suitable for diabetics
- Additional nail polishing attachment with 2 smoothing and polishing rollers
- Natural, durable nail gloss (without additional care products)
- Simple to use, easy replacement of the attachments
- With storage bag



Including 7 attachments: Sapphire cone (coarse), Sapphire disc (fine), Cylindrical burr, Flame burr, Felt cone, Needle burr as well as additional attachment for nail polishing

Size: approx. 3,8 x 16,1 cm
Weight: approx. 0,27 kg
Power supply: 230 V, 50 Hz, DC18V/500mA
Rotation speed: 3 levels: 5.000/6.000/7.000
Power output: 9 Watts
Cable length: 3 m
Article No.: 85155
EAN: 40 15588 85155 1
Warranty: 3 Years

Manicure and pedicure device

MPS



- Care and treatment of nails, cuticles and slight calluses
- Continuous speed adjustment

- Clockwise and counter-clockwise rotation
- High quality and durable sapphire attachments
- Stops automatically when the pressure is too high – also suitable for diabetics
- With storage bag



Including 4 attachments: Sapphire cone (coarse), Sapphire disc (fine), Flame burr, Felt cone

Size: approx. 15 x 4 x 4 cm
Weight: approx. 0,13 kg
Power supply: 100-240 V, 50/60 Hz, 18V DC 500 mA
Rotation speed: 3800-5500 U/min
Article No.: 85126
EAN: 40 15588 85126 1
Warranty: 3 Years

Manicure and pedicure device

MP 815



NEW



- **Top price-quality ratio**
- Care and treatment of nails, cuticles and slight calluses
- 2 speed levels for individual needs of nail care
- High quality and durable attachments

- Clockwise and counter-clockwise rotation
- Simple to use, easy replacement of the attachments
- With storage bag




Including 7 attachments: Sapphire disc (fine), Sapphire cone (coarse), Sapphire callus remover, Sapphire cone, Felt cone, Cylindrical burr and Flame burr

Size: approx. 15 x 3,8 x 3,7 cm
Weight: approx. 101 g
Power supply: 230 V, 50 Hz, DC18V/ 500mA
Rotation speed: 2 levels: 3.800/4.600
Cable length: 1,5 m
Article No.: 85153
EAN: 40 15588 85153 7
Warranty: 3 Years

Electric callus remover

CR 860



- Gently and effectively remove excess calluses
- For smoother hands, feet and elbows
- Apply to dry skin

- 2 speed levels
- Ergonomic, non-slip handle
- Including cleaning brush, protective cover, storage bag and batteries



Including 2 attachments (fine and coarse)

Size: approx. 15,7 x 6,7 x 4,6 cm
Weight: approx. 0,175 kg (including batteries)
Batteries: 2 x 1,5 V (AA)
Rotation speed: 2 levels: 1.600/1.800 U/min
Article No.: 88570
EAN: 40 15588 88570 9
Warranty: 3 Years

Electric nail polisher

NP 860



- **Long lasting shine without chemicals**
- For glossy, natural-looking nails
- Professional nail care at home
- Durable nail gloss (without additional care products)

- Simple to use on fingernails and toenails
- Small and handy – ideal even when travelling
- 2 speed levels for individual needs of nail care
- Ergonomic design
- Smoothing and polishing – in two steps to nails with the “Wow” effect
- Including 2 batteries (1.5 V AA)



Including 2 smoothing and polishing rollers

Size: approx. 15,5 x 3,6 x 3,3 cm
Weight: approx. 0,15 kg (including batteries)
Batteries: 2 x 1,5 V (AA)
Rotation speed: 2 levels: 2.000/2.500 U/min
Article No.: 85160
EAN: 40 15588 85160 5
Warranty: 3 Years

UV nail hardener

ND 870



NEW



- Professional nail design for home use
- Ideal for the curing of gel and liquid acrylic
- 360 degree reflective coating for a even curing
- Easy operation thanks to 4-step timer (90sec / 120sec / 30min / infinite)

- Detachable floor section for an easy lamp replacement
- 4 non-slip rubber feet for a firm standing

Size: approx. 253 x 221 x 113 mm
Weight: approx. 1,036 kg
Power supply: 230V~, 50 Hz; 36W
Bulb wave length: λ=365 nm
Article No.: 85470
EAN: 40 15588 85470 5
Warranty: 3 Years

Q3/2017

Q1/2018

3in1 Cosmetic mirror

CM 845



3in1 Function

- Illuminated cosmetic mirror including wall bracket and stand
- 2 swivelling mirrors: Normal and 5x magnification
- Uniform ring illumination with 17 white LED-lights
- Diameter of mirror approx. 13 cm

- High-quality chrome finish
- Including wall bracket, stand, batteries and screws for mounting onto the wall
- Including 4 batteries (1.5 V AAA)



Includes wall bracket and stand

Size: approx. 17,5 x 24/47 x 36,5 cm
Weight: approx. 1,55 kg
Diameter of mirror: approx. ø 13 cm
Batteries: 4 x 1,5 V (AAA)
Article No.: 88552
EAN: 40 15588 88552 5
Warranty: 3 Years

3in1 Cosmetic mirror

CM 840



Our bestseller

- Illuminated cosmetic mirror
- 2 swivelling mirrors: Normal and 5x magnification
- Uniform ring illumination with 18 LED-lights
- Diameter of mirror approx. 13 cm

- High-quality chrome finish
- Stable base
- Including 4 batteries (1.5 V AAA)

Size: approx. 20,5 x 13,5 x 31,5 cm
Weight: approx. 0,95 kg
Diameter of mirror: approx. ø 13 cm
Batteries: 4 x 1,5 V (AAA)
Article No.: 88550
EAN: 40 15588 88550 1
Warranty: 3 Years

Cosmetic mirror

CM 835



- Illuminated cosmetic mirror
- 2 swivelling mirrors: Normal and 5x magnification
- Uniform ring illumination with 12 white LED-lights
- Diameter of mirror approx. 12 cm

- High-quality chrome finish
- Stable base
- Including 3 batteries (1.5 V AAA)

Size: approx. 16,9 x 8 x 19 cm
Weight: approx. 0,375 kg
Diameter of mirror: approx. ø 12 cm
Batteries: 3 x 1,5 V (AAA)
Article No.: 88554
EAN: 40 15588 88554 9
Warranty: 3 Years

Facial cleansing brush

FB 885



- Facial cleansing brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products
- Rechargeable lithium-ion battery

- Cleanses the skin in 60 seconds (Timer)
- 4 speed settings
- Automatic shut-off after 1 min.
- Including 4 cleansing attachments (Sensitive skin, regular skin, peeling, cosmetics sponge)
- Splash-proof
- Including protective cover, storage bag and charge base



Including 4 cleansing attachments (Sensitive skin, regular skin, peeling, cosmetics sponge)

Size: approx. 18,6 x 5,0 x 2,7 cm
Weight: approx. 0,27 kg
Batteries: Lithium-ion battery
Power output: 500 mAh
Cable length: 1,30 m
Article No.: 88565
EAN: 40 15588 88565 5
Warranty: 3 Years

Facial cleansing brush

FB 880



- Facial cleansing brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products

- Cleanses the skin in 60 seconds (Timer)
- 2 speed settings: normal or intensive cleansing
- Automatic shut-off after 1 min.
- Including 3 cleansing attachments (Sensitive skin, regular skin, peeling)
- Batteries last up to 70 min.
- Splash-proof
- Including protective cover, storage bag and batteries (1,5 V AA)



Including 3 cleansing attachments (Sensitive skin, regular skin, peeling)

Size: approx. 15,7 x 5,7 x 7,9 cm
Weight: approx. 0,22 kg
Batteries: 2 x 1,5 V (AA)
Article No.: 88560
EAN: 40 15588 88560 0
Warranty: 3 Years

Facial sauna

FSS



- Opens the pores, cleans the skin and relaxes facial lines with humid heat
- Additional attachment for soothing aroma therapies
- Easy cleaning and handling

- For cosmetic use
- Diffuser for even steam distribution
- Additional nose piece
- High-quality heating element
- Indicator light
- Measuring cup for easy filling



Including nose piece and measuring cup

Size: approx. 30 x 20 x 19,5 cm
Weight: approx. 0,56 kg
Power supply: 100-240 V / 50-60 Hz, max. 100 Watts
Article No.: 88245
EAN: 40 15588 88245 6
Warranty: 3 Years

Esther, 31
Likes trying out new ingredients.

Lukas, 34
Is very accurate when cooking.

Lily, 3
Is not very accurate.

Staying healthy, cooking healthier,
and enjoying consciously.

Accurate kitchen helpers – Made for Life.

Also in the future everything concerning food and the kitchen will focus on indulgence. A health-conscious diet without additives is in high demand. Many seasonal and regional products are highly popular. The common meal is an important experience during the day, which is deliberately “celebrated”. Whether it is with the family, friends or alone. People who want to know about the ingredients in the food, they use good recipes and utilize high-quality kitchen appliances. MEDISANA supports this experience with kitchen appliances, which bring even more joy into the kitchen.

Digital Kitchen Scale with App KS 250

- Determination and display of calories, protein, fat, cholesterol, carbohydrates, fibres and sodium in an app
- It is possible to create own foods
- Sensor touch buttons
- Large LCD-display
- High quality stainless steel weighing surface



**SMARTER
WEIGHTING**
WITH
NUTRITIONAL VALUE APP



The high-quality kitchen scale is not just a visual highlight. In conjunction with the associated app, you will be provided with reliable information about the most important nutritional values such as calories, protein, fat, cholesterol, carbohydrates, dietary fibre and sodium. How? Select a food, and then weigh it. At the same time, it will maintain a nutritional diary of all past values.



Healthy food, create a good mood. The KS 250 provides the certainty that the diet will be balanced.

How balanced the diet is, will be output the nutritional value App.



Strictly speaking,
for the best results.

Enjoy your meal - thanks to good nutritional values. The kitchen scale with App helps to specify what and how much.

Lisa, 6
Loves cake.

Emma, 7
Is already a skilled cook.



Digital kitchen scale

KS 200



- Energy-saving automatic shut-off function
- Overload indicator
- High-quality steel weighting surface
- Anti-slip feet

- Practical add-weight function (tare)
- On/off switch
- Low battery display
- 3 kg max load capacity
- Including battery (CR2032 3V)

Size: approx. 16 x 11,7 x 2,7 cm
Weight: approx. 0,24 kg
Weighting surface: approx. 11 x 11 cm
LCD Display: approx. 1,6 x 4,6 cm
Weight Display: g, oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40469
EAN: 40 15588 40469 6
Warranty: 3 Years



Digital glass kitchen scale

KS 210



- Touch sensor keys
- Energy-saving automatic shut-off function
- Flat construction
- Glass weighting surface – easy to clean

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- 5 kg / 5,000 ml max load capacity
- Including battery (CR2032 3V)

Size: approx. 22,4 x 18,4 x 1,7 cm
Weight: approx. 0,6 kg
Weighting surface: approx. 22,4 x 18,4 cm
LCD Display: approx. 5,5 x 2,5 cm
Weight Display: g, kg, ml, fl.oz and lb.oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40463
EAN: 40 15588 40463 4
Warranty: 3 Years

Digital glass kitchen scale

KS 210



- Touch sensor keys
- Energy-saving automatic shut-off function
- Flat construction
- Glass weighting surface – easy to clean
- With photo motif

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- 5 kg / 5,000 ml max load capacity
- Including battery (CR2032 3V)

Size: approx. 22,4 x 18,4 x 1,7 cm
Weight: approx. 0,6 kg
Weighting surface: approx. 22,4 x 18,4 cm
LCD Display: approx. 5,5 x 2,5 cm
Weight Display: g, kg, ml, fl.oz and lb.oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40471
EAN: 40 15588 40471 9
Warranty: 3 Years

Digital glass kitchen scale

KS 210



- Touch sensor keys
- Energy-saving automatic shut-off function
- Flat construction
- Glass weighting surface – easy to clean
- With photo motif

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- 5 kg / 5,000 ml max load capacity
- Including battery (CR2032 3V)

Size: approx. 22,4 x 18,4 x 1,7 cm
Weight: approx. 0,6 kg
Weighting surface: approx. 22,4 x 18,4 cm
LCD Display: approx. 5,5 x 2,5 cm
Weight Display: g, kg, ml, fl.oz and lb.oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40472
EAN: 40 15588 40472 6
Warranty: 3 Years

Digital glass kitchen scale

KS 210



- Touch sensor keys
- Energy-saving automatic shut-off function
- Flat construction
- Glass weighting surface – easy to clean
- With photo motif

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- 5 kg / 5,000 ml max load capacity
- Including battery (CR2032 3V)

Size: approx. 22,4 x 18,4 x 1,7 cm
Weight: approx. 0,6 kg
Weighting surface: approx. 22,4 x 18,4 cm
LCD Display: approx. 5,5 x 2,5 cm
Weight Display: g, kg, ml, fl.oz and lb.oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40473
EAN: 40 15588 40473 3
Warranty: 3 Years

Digital glass kitchen scale

KS 220



- With time and temperature
- Energy-saving automatic shut-off function
- Timer Function
- Glass weighting surface – easy to clean

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- 5 kg max load capacity
- Including battery (CR2032 3V)

Size: approx. 15 x 22,7 x 3,6 cm
Weight: approx. 0,54 kg
Weighting surface: approx. Ø 15 cm
LCD Display: approx. 2,4 x 5,3 cm
Weight Display: g, kg, oz and lb.oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40467
EAN: 40 15588 40467 2
Warranty: 3 Years

Digital kitchen scale

KS 230



- Touch sensor keys
- Energy-saving automatic shut-off function
- Large, easy-to-read LCD display with white luminous digits
- High-quality stainless steel weighting surface
- Hanging fixture

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- 5 kg / 5,000 ml max load capacity
- Flat construction
- Including battery (CR2032 3V)

Size: approx. 20,2 x 21,8 x 1,6 cm
Weight: approx. 0,34 kg
Weighting surface: approx. 20,2 x 18,0 cm
LCD Display: approx. 2,3 x 5,3 cm
Weight Display: g, ml, oz, fl.oz and lb.oz
Graduation: 1 g
Batteries: 2 x CR2032 3V
Article No.: 40465
EAN: 40 15588 40465 8
Warranty: 3 Years

Digital glass kitchen scale XL

KS 240



- **Multi-functional utilization up to 20kg**
- Touch sensor keys
- HOLD / Memory function
- Energy-saving automatic shut-off function
- Glass weighting surface – easy to clean
- Hanging fixture

- Practical add-weight function (tare)
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- 20 kg / 20,000 ml max load capacity
- Flat construction
- Including battery (CR2032 3V)

Size: approx. 21 x 29,7 x 1,7 cm
Weight: approx. 0,82 kg
Weighting surface: approx. 21 x 29,7 cm
LCD Display: approx. 6,5 x 3,5 cm
Weight Display: g, kg, ml, fl.oz and lb.oz
Graduation: 1 g
Battery: 1 x CR2032 3V
Article No.: 40468
EAN: 40 15588 40468 9
Warranty: 3 Years

INFO Digital glass kitchen scale XL KS 240

Every kitchen should have them. Nicely shaped and able to take up to 20 kg.

Cooking and baking become even better. The KS 240 sensor-touch touch buttons quickly makes itself indispensable in the kitchen. Up to 20 kg can be weighed precisely. The overload display reacts at once to an overload. The glass surface is easy to clean. The weight displays at the top and on the side also help when weighing letters or parcels.

Digital kitchen scale with App

KS 250



- **Kitchen scale with nutritional value app**
- Determination and display of calories, protein, fat, cholesterol, carbohydrates, fibre and sodium in an app
- Sensor touch keypads
- Large, easy-to-read LCD-display with white luminous digits
- High-quality stainless steel weighting surface

- Practical add-weight function (tare)
- Energy-saving automatic shut-off function
- Overload indicator
- Anti-slip feet
- On/off switch
- Low battery display
- 5 kg / 5,000 ml max load capacity
- Flat construction
- Including 2 batteries (CR2032 3V)

- **App functions**
- Foods are simply selected in the app and then weighed
- Nutrition diary: keeping a record of all the measured data by date and time
- Multifunctional nutritional value determination
- More than 500 foods saved in the app
- Creating your own foods possible

Size: approx. 20,2 x 18,3 x 1,7 cm
Weight: approx. 0,34 kg
Weighting surface: approx. 20 x 18 cm
LCD Display: approx. 2,3 x 5,3 cm
Weight Display: g, oz, lb.oz, ml, fl.oz
Graduation: 1 g
Batteries: 2 x CR2032 3V
Article No.: 40474
EAN: 40 15588 40474 0
Warranty: 3 Years

Upper arm cuffs

for blood pressure monitors



MTM | MTP

Large cuff (32–42 cm)
Article No.: 51037
EAN: 40 15588 51037 3

Normal cuff (22–32 cm)
Article No.: 51038 (included in package)
EAN: 40 15588 51038 0

Small cuff (17–22 cm)
Article No.: 51034
EAN: 40 15588 51034 2

BU 535 | BU 535 Voice | MTV | MTC | MTD | MTS | MTR | CardioCompact

Large cuff (30–42 cm)
Article No.: 51126
EAN: 40 15588 51126 4

Normal cuff (22–30 cm)
Article No.: 51135 (included in package)
EAN: 40 15588 51135 6

MTX

Large „Easy form“ cuff (22–42 cm)
Article No.: 51087
EAN: 40 15588 51087 8

Normal „Easy form“ cuff (22–32 cm)
Article No.: 51088
EAN: 40 15588 51088 5

BU 550 Connect | BU 575 Connect

Normal cuff (22–42 cm)
Article No.: 51299
EAN: 40 15588 51299 5

BU 510 | BU 530 Connect

Large cuff (32–42 cm)
Article No.: 51169
EAN: 40 15588 51169 1

Normal cuff (22–36 cm)
Article No.: 51168 (included in package)
EAN: 40 15588 51168 4

CARDIODOCK | CARDIODOCK 2

Large „Easy form“ cuff (30–42 cm)
Article No.: 51284
EAN: 40 15588 51284 1

Normal „Easy form“ cuff (22–32 cm)
Article No.: 51283
EAN: 40 15588 51283 4

Mains adapters

for blood pressure monitors



Netzadapter für BU 550 Connect

Article No.: 51013
EAN: 40 15588 51013 7

Netzadapter für BU 575 Connect | BU 300 Connect

Article No.: 51013
EAN: 40 15588 51013 7

Netzadapter für MTM | MTP | MTX

Article No.: 51036
EAN: 40 15588 51036 6

Netzadapter für CardioCompact

Article No.: 51095
EAN: 40 15588 51095 3

Netzadapter für MTV | MTC | MTD | MTS | BU 510 Connect | BU 530 Connect

Article No.: 51125
EAN: 40 15588 51125 7

Protective sleeves

for thermometer



FTF

50 pieces
Article No.: 76111
EAN: 40 15588 76111 9

Accessories

for blood pressure monitors



Control solution for MediTouch 2 | MediTouch 2 Connect | MediTouch 2 Carbon

Article No.: 79039
EAN: 40 15588 79039 3

Lancets for MediTouch 2 | MediTouch 2 Connect | MediTouch 2 Carbon

Article No.: 79028
EAN: 40 15588 79028 7

Test strips for MediTouch 2 | MediTouch 2 Connect | MediTouch 2 Carbon

Article No.: 79038
EAN: 40 15588 79038 6

Lancing Device

Article No.: 79002
EAN: 40 15588 79002 7

USB cable for Meditouch 2

Article No.: 79049
EAN: 40 15588 79049 2

Wristbands

for activity tracker



ViFit Connect

Black
Article No.: 79417
EAN: 40 15588 79417 9

Orange
Article No.: 79418
EAN: 40 15588 79418 6

Blue
Article No.: 79470
EAN: 40 15588 79470 4

Green
Article No.: 79471
EAN: 40 15588 79471 1

Red
Article No.: 79472
EAN: 40 15588 79472 8

Pink
Article No.: 79473
EAN: 40 15588 79473 5

Spare part sets

for EMS Pad



Body Trainer

Spare part sets (2 pads for the abdomen)
Artikel-Nr.: 79525
EAN: 40 15588 79525 1

Spare part sets (4 pads for arms/legs)
Artikel-Nr.: 79526
EAN: 40 15588 79417 8

Accessories

for bath spas



BBS | MBH

Air hose
Article No.: 30666
EAN: 40 15588 30666 2

Soft mat
Article No.: 88383
EAN: 40 15588 88383 5

Fragrance dispenser with 3 pads
Article No.: 30669
EAN: 40 15588 30669 3

BBS

Remote control
Article No.: 88384
EAN: 40 15588 88384 2

MBH

Remote control
Article No.: 88376
EAN: 40 15588 88376 7

Loofah massage pad

for mini handmassager



HM 840

Article No.: 88526
EAN: 40 15588 88526 6



Electrodes
for pain therapy devices



TDD | TDP | TED
4 Multiple electrodes (approx. 3 x 5 cm)
Article No.: 87080
EAN: 40 15588 87080 4

2 Multiple electrodes (approx. 5 x 9 cm)
Article No.: 87100
EAN: 40 15588 87100 9

Note: Not suitable for hair removal system EPD (Article No.: 85550)

TDP
Cable set
Article No.: 88315
EAN: 40 15588 88315 6

BT 850
Adhesive pads (2 pcs.)
Article No.: 88334
EAN: 40 15588 88334 7

Replacement lamp
for infrared lamps



IRH
100 Watts
Article No.: 30101
EAN: 40 15588 30101 8

IRL
150 Watts
Article No.: 88255
EAN: 40 15588 88255 5

Accessories
for Inhalators



IN 500 | IN 550
Connecting piece
Article No.: 54524
EAN: 40 15588 54524 5

Air filter and nebulizer kit
Article No.: 54525
EAN: 40 15588 54525 2

Mouthpiece and nose adapter
Article No.: 54526
EAN: 40 15588 54526 9

Air hose
Article No.: 54527
EAN: 40 15588 54527 6

Mask for children
Article No.: 54528
EAN: 40 15588 54528 3

Mask for adults
Article No.: 54529
EAN: 40 15588 54529 0

Nebulizer attachments (3 pcs.)
Article No.: 54531
EAN: 40 15588 54531 3

USC
Mask set (2 masks + adapter)
Article No.: 54101
EAN: 40 15588 54101 8

Mouthpiece
Article No.: 54102
EAN: 40 15588 54102 5

Covering cap
Article No.: 54103
EAN: 40 15588 54103 2

Nebulizer Kit
Article No.: 54104
EAN: 40 15588 54104 9

Aroma essences
for Aroma diffusers
and air humidifiers



Fragrance Lemon 10 ml
Article No.: 60030
EAN: 40 15588 60030 2

Fragrance Eucalyptus 10 ml
Article No.: 60031
EAN: 40 15588 60031 9

Fragrance Lavender 10 ml
Article No.: 60032
EAN: 40 15588 60032 6

Fragrance Citronella 10 ml
Article No.: 60033
EAN: 40 15588 60033 3

Fragrance Ylang-Ylang 10 ml
Article No.: 60035
EAN: 40 15588 60035 7

Fragrance Peppermint 10 ml
Article No.: 60036
EAN: 40 15588 60036 4

Fragrance Orange 10 ml
Article No.: 60037
EAN: 40 15588 60037 1

Fragrance Pinewood 10 ml
Article No.: 60038
EAN: 40 15588 60038 8

Fragrance Wildberries 10 ml
Article No.: 60039
EAN: 40 15588 60039 5

Fragrance Oceanbreeze 10 ml
Article No.: 60040
EAN: 40 15588 60040 1

Replacement filter
for air purifier



MEDISANA AIR®
1 pc.
Article No.: 60390
EAN: 40 15588 60390 7

2 pcs.
Article No.: 60391
EAN: 40 15588 60391 4

Replacement water filter
for air humidifier



AH 660
2 pcs.
Article No.: 60069
EAN: 40 15588 60069 2

Replacement brushes
for facial cleansing brushes



FB 880
3 pcs.
Article No.: 88561
EAN: 40 15588 88561 7

FB 885
4 pcs.
Article No.: 88566
EAN: 40 15588 88566 2

Replacement cartridge
for hair removal system



IPL 800
Article No.: 88581
EAN: 40 15588 88581 5

IPL 805
Article No.: 88584
EAN: 40 15588 88584 6

IPL 836 | IPL 840
Article No.: 88586
EAN: 40 15588 88586 0

IPL 850
Article No.: 88589
EAN: 40 15588 88589 1

Replacement rolls
for electric callus remover



CR 860
2 pcs.
Article No.: 88572
EAN: 40 15588 88572 3

Replacement rolls
for manicure and pedicure devices



NP 860 | MP 840
4 pcs.
Article No.: 85164
EAN: 40 15588 85164 3

Attachments
for manicure and pedicure devices



Manilux I | Manilux s | MP 820
Set of attachments (10 pieces)*
Article No.: 85450
EAN: 40 15588 85450 7

*contains: sapphire cone (coarse), sapphire disc (coarse), flame burr, sapphire ball burr, felt cone, sapphire cone (long), sapphire disc (fine), needle burr (long), cylindrical burr, drainage grinder

MEDISTYLE S | MEDISTYLE L | MPS | MP 840

Felt cone
Article No.: 85001
EAN: 40 15588 85001 1

Sapphire cone (coarse)
Article No.: 85003
EAN: 40 15588 85003 5

Sapphire disc (coarse)
Article No.: 85004
EAN: 40 15588 85004 2

Sapphire ball burr
Article No.: 85005
EAN: 40 15588 85005 9

Needle burr (long)
Article No.: 85006
EAN: 40 15588 85006 6

Sapphire cone (long)
Article No.: 85007
EAN: 40 15588 85007 3

Cuticle pusher
Article No.: 85008
EAN: 40 15588 85008 0

Cylindrical burr
Article No.: 85009
EAN: 40 15588 85009 7

Sapphire disc (coarse)
Article No.: 85010
EAN: 40 15588 85010 3

Flame burr
Article No.: 85011
EAN: 40 15588 85011 0

Attachment adapter
Article No.: 85015
EAN: 40 15588 85015 8

Grinding set for artificial nails
Article No.: 85460
EAN: 40 15588 85460 6

Sanding Bands (30 pieces)
Article No.: 85461
EAN: 40 15588 85461 3



MEDISANA®

Jagenbergstraße 19

41468 Neuss

Germany

Phone: +49 2131 / 36 68-0

Fax: +49 2131 / 36 68-50 95

E-Mail: info@medisana.de

Internet: www.medisana.de



September 2017